

Quick Reference Guides are not a substitute for the User Manuals. Please refer to the User Manuals until you are an experienced user.

**1 Prepare Reagents**

The following reagents are supplied by NanoString:

Kit	Reagents	Lot Number
GeoMx DSP Protein Slide Prep Kit for FFPE	Buffer S (4°C/RT) Buffer W (4°C)	
GeoMx Morphology Kit for Protein	Morphology markers (4°C) Nuclear Stain—SYTO 13 (-20°C)	
GeoMx Core & optional Module Kits for NGS or nCounter readout	Antibody (Ab) Mix (-80°C)	

**For nCounter assays:** Not all nCounter protein modules are compatible with one another. Each module is assigned a Probe R space. **Do not combine two modules with a common Probe R space** (i.e. same Probe R\_number) in the same experimental run, or the data cannot be interpreted. Instead, use Substitute Probe R (available from NanoString). See nCounter protein modules' Probe R designations in the Appendix of the **GeoMx DSP Manual Slide Preparation User Manual** and use the planning grid at right, if needed.

See the **GeoMx DSP Automated Slide Preparation User Manual (MAN-10151)** for more information.

The following reagents are not supplied by NanoString:

Reagents	
BOND Dewax Solution, 1 L	DEPC-treated water
BOND Epitope Retrieval 1, 1 L	1X TBS
BOND Wash Solution 10X Concentrate, 1 L (Dilute to 1X with DEPC-treated water)	1X TBS-T
4% PFA	RNase AWAY*


R #	R_1	R_2	R_3	R_4	R_5	R_6	R_7	R_8	R_9
Core/Module to be used									

**2 Prepare Samples**

1. Use Apex BOND Adhesive slides. Ensure your tissue sample fits inside the **Scan Area** of the slide (green, in diagram at right) and within the white dots of the BOND slide. Tissue may not overlap gasket (blue) or tip calibration area (red). Important: If the tissue sample is outside of this area, the GeoMx will not process the slide.

**Important:** If the tissue sample is outside of this area, the GeoMx will not process the slide

2. **Bake slides** at 60°C for 30 minutes to 3 hours.



**3 Set up the BOND RX/RX<sup>m</sup>**

To set up the BOND RX/RX<sup>m</sup> Semi-Automated Protein Protocol for the first time, refer to the GeoMx DSP Automated Slide Preparation User Manual (MAN-10151) for instructions.


## Day 1 (continued)

### 4 Start a Run on the BOND RX/RX<sup>m</sup>

1. Fill the 30 mL Home Buffer container with 1X Bond Wash Buffer and place it in slot 1 of the Bond Research Detection System.
2. Fill the Titration Container insert with Buffer W and place in slot 2 of the Bond Research Detection System (150  $\mu$ L per slide + 350  $\mu$ L of dead volume).



3. Load the Bond Research Detection System into the BOND RX/RX<sup>m</sup>.



4. From Slide Setup, select Add Study. Enter study info; dispense volume 150  $\mu$ L; Preparation protocol \*Bake and Dewax. Select OK.
5. Select Add Slide. Enter slide information as shown below, then click Add slide.



A screenshot of the 'Add slide' software interface. It includes fields for 'Tissue TMA', 'Tissue type' (Test tissue, Negative tissue, Positive tissue), 'Dispense volume' (100  $\mu$ L, 150  $\mu$ L), 'Staining mode' (Single, Routine), 'Process' (IHC, ISH), 'Marker' (\*Negative), 'Protocols', 'Staining' (IHC DSP Protocol), 'Preparation' (\*Bake and Dewax), 'HIER' (\*HIER 20 min with ER1), and 'Enzyme' (\*...).

6. Print and attach labels to the slides.
7. Load slides with covertiles into the slide tray. Ensure positioning is correct.



### 5 Load and Run your Slides

1. Load and run your slides. The run will take approximately 3 hours.
2. When run is complete, remove covertiles from one slide at a time and place slide in 1X PBS. Remove excess buffer and draw a closed barrier with a hydrophobic pen. DO NOT let slides dry out.



## Day 1 (continued)

### 6 Primary Antibody Incubation

**Thaw Ab mix**, keeping the aliquot on ice and protected from light. Do not vortex. **Make working antibody solution** (see table). Handling one slide at a time, **blot excess buffer** and place slide in humidity chamber. **Cover tissue with working Ab solution**. Make sure the entire tissue is covered (~ 200  $\mu$ L per slide). **Incubate overnight** in closed humidity chamber at 4°C, protected from light.



**Humidity Chamber**  
 Temp: 4°C  
 Time: Overnight

Working Ab solution for NGS or nCounter assays (n = number of slides):

Core Ab	Module Ab 1	Module Ab 2	Other Modules*	Morph Marker 1	Morph Marker 2	Other Markers	Buffer W	Total Volume
(8 $\mu$ L x n)	(8 $\mu$ L x n)	(8 $\mu$ L x n)	...	(5 $\mu$ L x n)	(5 $\mu$ L x n)	...	(up to 200 $\mu$ L)**	(200 $\mu$ L x n)
_____	+ _____	+ _____		+ _____	+ _____		+ _____	= _____

\*If adding custom barcoded Ab, see Appendix V of the GeoMx DSP Manual Slide Preparation User Manual for instructions.

\*\*If a different number of module Ab or morphology markers are used, Buffer W volume can be adjusted; final volume = 200  $\mu$ L/slide.


## Day 2

### 7 Post-Fix

Staining jars should be thoroughly cleaned with RNase AWAY before use. Handling one slide at a time, wash in 1X TBS-T (see right for time and repetitions). Tap slide on paper towels to remove excess solution. Place slide in humidity chamber. Cover tissue with 4% PFA (~200  $\mu$ L per slide). Incubate 30 minutes at room temp in closed humidity chamber. Wash with 1X TBS-T (see right).



**4% PFA in Humidity Chamber**  
 Room Temp  
 Time: 30 min



### 8 Nuclei Staining

Picofuge thawed SYTO 13 1 min; do not pipette from the bottom of the vial. Dilute 1:10 in 1X TBS then mix by pipetting up and down. Handling one slide at a time, place in humidity chamber. Cover tissue with diluted SYTO 13 (~200  $\mu$ L per slide). Incubate 15 minutes at room temp in closed humidity chamber. Dip slide into fresh 1X TBS-T, then transfer it to fresh 1X TBS-T. Scrape off hydrophobic pen on one slide at a time with a razor blade. DO NOT allow slide to dry out. Store in 1X TBS-T.



**Nuclear stain in Humidity Chamber**  
 Room Temp  
 Time: 15 min



**1X TBS-T**  
 2 washes, then store

## 9 Additional resources in GeoMx DSP Automated Slide Preparation User Manual (MAN-10151)

- Slide storage guidelines
- Fresh frozen and fixed frozen sample preparation
- Alternative morphology marker methods
- Troubleshooting

For more information, please visit [nanosttring.com](https://nanosttring.com)

#### NanoString Technologies, Inc.

530 Fairview Avenue North  
Seattle, Washington 98109

T (888) 358-6266  
F (206) 378-6288

nanosttring.com  
support@nanosttring.com

#### Additional Customer Resources

For more comprehensive information, visit us at [nanosttring.com](https://nanosttring.com) and go to [Support > Product Support](#) to view manuals and other technical product literature.

FOR RESEARCH USE ONLY. Not for use in diagnostic procedures.

© 2022 NanoString Technologies, Inc. All rights reserved. NanoString, NanoString Technologies, the NanoString logo, nCounter and GeoMx are registered trademarks or trademarks of NanoString Technologies, Inc. in the United States and/or other countries. Leica and BOND RX are registered trademarks of Leica Biosystems.

MAY 2022 MK4780

nanosttring®