

## Sanger Tree of Life HMW DNA Fragmentation: Covaris g-TUBE for ULI PacBio

### Authors:

Graeme Oatley, Filipa Sampaio, Lucy Kitchin, Raquel Vionette Do Amaral, Caroline Howard

### Abstract:

This protocol is for centrifugation-mediated fragmentation of HMW DNA from any of the Sanger Tree of Life HMW DNA extraction protocols intended for ultra-low input (ULI) sequencing, into an average fragment size range of 8–10 kb, using the Covaris g-Tube. This process is highly effective for DNA extracts from all of the taxonomic groups covered by the Tree of Life Programme. The output of this protocol is sheared DNA which can be directed towards Fragmented DNA clean up, using either the Manual or Automated SPRI protocols.

### Keywords/Acronyms

HMW DNA fragmentation, g-TUBE, ultra-low input, reference genome, long read sequencing

### Safety Warnings

- The operator must wear a lab coat, powder-free nitrile gloves and safety specs to perform the laboratory procedures in this protocol.
- Waste needs to be collected in a suitable container (e.g. plastic screw-top jar or Biobin) and disposed of in accordance with local regulations.

### Guidelines

- The DNA sheared using this protocol must be intended for ULI PacBio sequencing.
- This protocol is for the shearing of DNA into fragment sizes of 8–10 kb.
- 150 ng of unsheared DNA is ideal for this protocol, but 90–150 ng total input can be used.

### Laboratory Protocol:

1. Label the required number of Covaris g-TUBEs for each DNA sample that will be sheared; ensure to label the tubes both on the lid and on the bottom.
2. Prior to transferring the DNA sample from its original tube, first mix the DNA sample by pipetting carefully with wide-bore pipette tip.
3. Transfer between 100–150  $\mu$ L of the DNA sample to its corresponding labelled g-TUBE; 150  $\mu$ L is the ideal volume.
4. If processing multiple samples, ensure that the volumes are equal by adding EB buffer to normalise the volumes.  
Be aware that the processing time from transfer of the sample into the g-TUBE and the g-TUBE containing the sample into the centrifuge should be no longer than 15 minutes, to avoid the sample migrating into the filter whilst on the bench.

5. Secure the g-TUBEs by tightly sealing the screw-caps, using the g-TUBE prep station provided with the g-TUBEs.
6. Place the g-TUBEs into a bench-top centrifuge. If there is an uneven number of samples, make sure to balance the centrifuge using a spare g-TUBE.
7. Run for 1 minute at 5000 rpm for a shear size of 10 kb; refer to the manufacturer's instructions for conditions to obtain different shear sizes.
8. Repeat the centrifugation in step 7 until all the volume has passed through the filter, up to three times. Check the tube between chambers to ensure that all of your sample has passed through the filter. If the sample has not passed through the filter after the third spin, please refer to the 'Troubleshooting' section in the Guidelines section.
9. Invert the g-TUBEs within the prep station so the fixed base is facing up and the screw-cap is facing down, then repeat steps 6, 7 and 8 so that all of the DNA has passed completely through the g-TUBE filter a total of two times.
10. Once all of the sample has passed through the filter, remove the g-TUBEs from the centrifuge whilst maintaining their current orientation (screw-cap facing down), and place them screw-cap side down into the g-TUBE prep station.
11. Using the g-TUBE prep station to stabilise the tube, carefully unscrew the tube from the lid, leaving the lid in the prep station.
12. Aspirate the sample from the screw-cap using a standard pipette tip and transfer the sheared DNA to a fresh DNA Lo-Bind microcentrifuge tube.
13. Store the sheared DNA at 4°C until further processing.

### **Troubleshooting:**

#### **1) Low concentration samples**

In some cases you may need to centrifuge >150 µL in order to process 150 ng. In this scenario, split your DNA sample into 100 µL or 150 µL subsamples, and centrifuge each subsample through the same g-TUBE. Re-combine the outputs from each centrifugation in preparation for downstream processes.

#### **2) Sample does not spin through the filter**

If, after three spins, the sample has not passed through the filter and still remains in the chamber, try the following:

- Flip the g-TUBE (screw cap down), and spin at 4300 x g for 30 seconds.
- Recover the sample from the cap following the usual procedure (step 11).
- Filter the sample with a 0.45 µm filter.
- Use a new g-TUBE and follow the normal protocol.

**Materials:**

- 1.5 mL DNA Lo-Bind microcentrifuge tubes (Eppendorf Cat. no. 0030 108.051)
- Covaris g-TUBE with Covaris g-TUBE Prep Station (Cat. no. 520079)
- Buffer EB (Qiagen Cat. no. 19086)

**Equipment**

- Pipettes for 0.5 to 1000  $\mu$ L and filtered tips
- Wide bore tips (200  $\mu$ L, filtered if available)
- Eppendorf™ Centrifuge 5425/5425 R (Cat. no. 5405000263)

**References:**

User Manual g-TUBE: [pn\\_010154.pdf \(covaris.com\)](#)