

MANUAL

# CosMx<sup>®</sup> SMI Manual Slide Preparation for RNA Assays

MAN-10184-06

May 2025

**FOR RESEARCH USE ONLY.** Not for use in diagnostic procedures.

Innovation with Integrity

**CosMx<sup>®</sup>**  
Spatial Molecular Imager

Bruker Spatial Biology, Inc.  
500 Fairview Ave N  
Seattle, Washington 98109

[www.brukerspatialbiology.com](http://www.brukerspatialbiology.com)

Tel: 888.358.6266

Fax: 206.378.6288

Technical Support: [support.spatial@bruker.com](mailto:support.spatial@bruker.com)

Customer Service: [customerservice.spatial@bruker.com](mailto:customerservice.spatial@bruker.com)

Sales Contacts

North America: [nasales.bsb@bruker.com](mailto:nasales.bsb@bruker.com)

EMEA: [emeasales.bsb@bruker.com](mailto:emeasales.bsb@bruker.com)

APAC: [apacsales.bsb@bruker.com](mailto:apacsales.bsb@bruker.com)

All other regions: [globalsales.bsb@bruker.com](mailto:globalsales.bsb@bruker.com)



**EMERGO EUROPE**

Westervoortsedijk 60

6827 AT Arnhem

The Netherlands

## Rights, License, and Trademarks

### Use

For Research Use Only. Not for use in diagnostic procedures.

### Intellectual Property Rights

This CosMx<sup>®</sup> Spatial Molecular Imager (SMI) User Manual and its contents are the property of Bruker Spatial Biology, Inc. ("Bruker"), and are intended solely for use by Bruker customers, for the purpose of operating the CosMx SMI System. The CosMx SMI System (including both its software and hardware components), this User Manual, and any other documentation provided to you by Bruker in connection therewith, are subject to patents, copyright, trade secret rights and other intellectual property rights owned by, or licensed to, Bruker. No part of the software or hardware may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into other languages without the prior written consent of Bruker. For a list of patents, see [www.nanostring.com/company/patents](http://www.nanostring.com/company/patents).

### Limited License

Subject to the terms and conditions of the CosMx SMI System contained in the product quotation, Bruker grants you a limited, non-exclusive, non-transferable, non-sublicensable, research use only license to use the proprietary CosMx SMI System only in accordance with the manual and other written instructions provided by Bruker. Except as expressly set forth in the terms and conditions, no right or license, whether express, implied or statutory, is granted by Bruker under any intellectual property right owned by, or licensed to, Bruker by virtue of the supply of the proprietary CosMx SMI System. Without limiting the foregoing, no right or license, whether express, implied or statutory, is granted by Bruker to use the CosMx SMI System with any third party product not supplied or licensed to you by Bruker or recommended for use by Bruker in a manual or other written instruction provided by Bruker.

### Trademarks

Bruker, Bruker Spatial Biology, the NanoString logo, CosMx, AtoMx, GeoMx, and nCounter are trademarks or registered trademarks of Bruker Spatial Biology, Inc., in the United States and/or other countries. All other trademarks and/or service marks not owned by Bruker Spatial Biology, Inc. that appear in this document are the property of their respective owners.

### Open Source Software Licenses

Visit <http://nanostring.com/cosmx-oss> and <http://nanostring.com/atomx-oss> for a list of open source software licenses used in CosMx Spatial Molecular Imaging.

### Copyright

©2025 Bruker Spatial Biology, Inc. All rights reserved.

## Table of Contents

---

<b>Manual Slide Preparation for RNA Assays</b> .....	<b>1</b>
<b>Contacts</b> .....	<b>2</b>
<b>Rights, License, and Trademarks</b> .....	<b>3</b>
<b>Table of Contents</b> .....	<b>4</b>
<b>Changes in this Revision</b> .....	<b>6</b>
<b>Conventions</b> .....	<b>7</b>
<b>Safety</b> .....	<b>7</b>
<b>Introduction to CosMx SMI</b> .....	<b>8</b>
<b>CosMx SMI User Manuals and Resources</b> .....	<b>9</b>
<b>Panel and Cell Segmentation Marker Selection</b> .....	<b>10</b>
<b>RNA FFPE Manual Slide Preparation</b> .....	<b>11</b>
RNA FFPE Manual Slide Preparation Overview .....	12
Equipment, Materials, and Reagents .....	13
Bruker Supplied Reagents .....	16
Tissue Sectioning and Slide Preparation: FFPE Samples .....	20
Prepare RNA FFPE Assay Reagents .....	22
Day 0: Prepare Shelf Stable Reagents and Bake Slides Overnight .....	24
Day 1: Deparaffinize FFPE Tissue Sections (20 minutes) .....	25
Perform Target Retrieval (50 minutes) .....	28
Tissue Permeabilization (40 minutes) .....	30
Fiducial Preparation and Application (20 minutes) .....	33
Post-Fixation (20 minutes) .....	35
NHS-Acetate Preparation and Application (25 minutes) .....	36
In Situ Hybridization (overnight) .....	37
Day 2: Perform Stringent Washes (90 minutes) .....	42
Nuclear and Cell Segmentation Staining (2 hours) .....	43
<b>RNA Fresh Frozen Manual Slide Preparation</b> .....	<b>46</b>
RNA FF Manual Slide Preparation Overview .....	47
Equipment, Materials, and Reagents .....	48
Bruker Supplied Reagents .....	51
Tissue Sectioning and Slide Preparation: Fresh Frozen Samples .....	55

---

Prepare RNA Fresh Frozen Assay Reagents .....	56
Day 1: NBF Fixation and Bake .....	58
Prepare Materials and Tissue for Target Retrieval (1 hour) .....	59
Perform Target Retrieval (50 minutes) .....	61
Tissue Permeabilization (40 minutes) .....	63
Fiducial Preparation and Application (20 minutes) .....	66
Post-Fixation (20 minutes) .....	68
NHS-Acetate Preparation and Application (25 minutes) .....	69
In Situ Hybridization (overnight) .....	70
Day 2: Perform Stringent Washes (90 minutes) .....	75
Nuclear and Cell Segmentation Staining (2 hours) .....	76
<b>Safe Storage Guidelines for RNA Slides .....</b>	<b>79</b>
<b>Flow Cell Assembly .....</b>	<b>80</b>
<b>Appendix I: CosMx SMI Tissue Sectioning Guidelines .....</b>	<b>84</b>
<b>Appendix II: Tissue-Specific Slide Preparation Considerations .....</b>	<b>88</b>
<b>Troubleshooting .....</b>	<b>92</b>
Poor Tissue Adherence to Slide .....	92
Tissue Peeling following Digestion (RNA protocol) .....	92
High Autofluorescence in Tissue (RNA protocol) .....	93
"Fluidic Issue Detected" during Deck Validation (on-instrument) .....	94
"Bottle Disconnected" during Deck Validation (on-instrument) .....	95
Critical Error during Instrument Preparation (on-instrument) .....	96
Blurry Preview Scan (on-instrument) .....	97
<b>Final page .....</b>	<b>98</b>

## Changes in this Revision

Changes in this manual revision include:

- Updated content in preparation for launch of the Whole Transcriptome (WTX) RNA Panel.
- A revised slide preparation method specifically for FFPE human brain tissue, found to improve RNA assay performance for this tissue type on CosMx SMI. The new method for FFPE human brain tissue has the following changes:
  - Slides undergo an additional wash in 70% ethanol for 2 hours (up to overnight).
  - Tissue digestion buffer is prepared in PBS-T (0.5% Tween-20).
  - Fiducials solution is prepared at 0.0005%, half the concentration as in the non-FFPE human brain method.

These changes are written in line into the FFPE protocol beginning [on page 11](#). Look for specific notations and instructions if preparing FFPE human brain tissue.

- Minor edits for clarification.


## Conventions


The following conventions are used throughout this manual and are described for your reference.

**Bold** text is typically used to highlight a specific button, keystroke, or menu option. It may also be used to highlight important text or terms.


Blue underlined text is typically used to highlight links and/or references to other sections of the manual. It may also be used to highlight references to other manuals and/or instructional material.


The gray box indicates general information that may be useful for improving assay performance. The notes may clarify other instructions or provide guidance to improve the efficiency of the assay workflow .

 **WARNING:** This symbol indicates the potential for bodily injury or damage to the instrument if the instructions are not followed correctly. Always carefully read and follow the instructions accompanied by this symbol to avoid potential hazards.

 **IMPORTANT:** This symbol indicates important information that is critical to ensure a successful assay. Following these instructions may help improve the quality of your data.

## Safety

 **WARNING:** Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eye wear, clothing, and gloves. SDSs are available from <https://nanosting.com/resources/safety-data-sheets>.

 **IMPORTANT:** Read all steps before you begin to familiarize yourself with this procedure.

## Introduction to CosMx SMI

The CosMx SMI platform is an integrated system with cyclic in situ hybridization chemistry, a high-resolution imaging readout instrument, and an interactive data analysis and visualization software. The CosMx SMI platform enables rapid quantification and visualization of more than 18,000 RNA and 64 protein analytes. This flexible spatial single-cell solution drives deeper insights into the cell atlas, cell-cell interaction, cellular processes, and biomarker discovery.




Figure 1: CosMx SMI is an integrated system that includes validated reagents and consumables, an instrument for chemistry and readout, and an interactive cloud-based software suite for data analysis.

Sample preparation involves basic in situ hybridization (ISH) processing steps. The protocols are compatible with the recommended glass pathology slides, and do not require complicated tissue expansion or clearing, cDNA synthesis or amplification.

RNA or protein targets in individual cells are identified via hybridization or binding with highly specific probes or antibodies labeled with a unique barcode system. Barcode readout occurs through multiple rounds of reporter probe binding and fluorescence imaging using the CosMx SMI instrument. Each RNA target appears as a distinct bright spot in the sample and is digitally quantified in the image. The data is then migrated to the cloud-based AtoMx<sup>®</sup> Spatial Informatics Platform (SIP) for analysis and visualization. Within AtoMx SIP, users can incorporate custom analysis workflows.

The CosMx Spatial Molecular Imager is the first platform to demonstrate simultaneous single-cell and sub-cellular detection of over 64 proteins on standard, bio-banked, FFPE tissue samples. The CosMx Protein technology uses an antibody-oligonucleotide conjugate to detect each protein's (sub)-cellular localization and quantify its expression level. CosMx oligo-labeled antibodies undergo rigorous QC testing, and site-specific labeling chemistry to select for pure imaging reagents with no unconjugated antibody or free oligonucleotide contamination, which could lead to background noise.

## CosMx SMI User Manuals and Resources

The CosMx SMI workflow is divided into the following user manuals:

Workflow Step 1	<a href="#">CosMx SMI Manual Slide Preparation for RNA Assays</a> MAN-10184  <a href="#">CosMx SMI Manual Slide Preparation for Protein Assays</a> MAN-10185	<a href="#">CosMx SMI Semi-Automated Slide Preparation for RNA Assays</a> MAN-10186  <a href="#">CosMx SMI Semi-Automated Slide Preparation for Protein Assays</a> MAN-10187	<a href="#">CosMx SMI Manual Slide Preparation for Multiomic Assays</a> MAN-10201
Workflow Step 2	<a href="#">CosMx SMI Instrument User Manual</a> MAN-10161		
Workflow Step 3	<a href="#">CosMx SMI Data Analysis User Manual</a> MAN-10162		

User manuals and other documents can be found online in the NanoString University Document Library at <https://university.nanostring.com>.

Instrument and workflow training courses are also available in NanoString University.

For information about the AtoMx Spatial Informatics Platform, please refer to the [AtoMx . SIP . Platform Administration Manual \(MAN-10170\)](#).

Additional data analysis support and resources can be found at <https://github.com/Nanostring-Biostats>.


## Panel and Cell Segmentation Marker Selection

Bruker Spatial Biology currently provides **8 pre-defined panels** for use with CosMx SMI:

- Human Whole Transcriptome Panel (WTX), 19K-plex, RNA
- Human 6K Discovery Panel, 6K-plex, RNA
- Human Universal Cell Characterization Panel, 1000-plex, RNA
- Human Immuno-oncology Panel, 100-plex, RNA
- Mouse Neuroscience Panel, 1000-plex, RNA
- Mouse Universal Cell Characterization Panel, 1000-plex, RNA
- Human Immuno-oncology Panel, 64-plex, Protein
- Mouse Neuroscience Panel, 64-plex, Protein (includes the Mouse Alzheimer's Pathology Module, which must be run with the core Mouse Neuroscience Panel)

In addition, these **custom add-on, swap-in, and stand-alone panels** are available:

- **Add on** up to 200 user-defined genes to the Human Whole Transcriptome Panel or 6K-plex Discovery Panel. NOTE: Custom Add-On probes are only compatible with the panel they are designed for. They are not interchangeable with different panels.
- **Add on** up to 8 user-defined protein targets to the 64-plex Protein Panels.
- **Swap in** 7-10 user-defined genes to the Human Immuno-Oncology 100-plex RNA Panel or 7-50 user-defined genes to any 1000-plex RNA Panel.
- **Stand-alone, de novo RNA panels** can be customized for up to 1000 targets. Collaborate with Bruker's Bioinformatics team to build a made-to-order custom panel. Email Bruker Support at [support.spatial@bruker.com](mailto:support.spatial@bruker.com) or visit <https://nanosttring.com/about-us/contact-us/> for custom panel design.

 **IMPORTANT:** If using any customized panel, email [AtoMx.KitAdmin@Bruker.com](mailto:AtoMx.KitAdmin@Bruker.com) at least one business day prior to the CosMx SMI run to add the custom kit to the Control Center.

### CosMx SMI Cell Segmentation Marker Selection

Bruker Spatial Biology provides panel-specific segmentation marker kits for each analyte. See the section on Bruker Supplied Reagents for available kits. Contact your Applications Scientist to discuss cell segmentation marker selection.

## RNA FFPE Manual Slide Preparation

### Slide Preparation Workflow




Figure 2: CosMx SMI Workflow Overview

**Day 1: Slide Preparation.** Prepare slides and incubate biological targets with RNA specific probes. Prepare manually or using the BOND RX/RX<sup>m</sup> fully automated IHC/ISH stainer from Leica Biosystems (BOND RX/RX<sup>m</sup>).

**Day 2: Process Slides on CosMx SMI.** Remove off-target probes and add cell segmentation markers to each slide. Load assembled flow cells into the CosMx SMI instrument and enter flow cell/study information. Tissue is scanned to capture RNA readout and morphology imaging within user-designated fields of view (FOVs).

**After run completion: Create a Data Analysis study** in the AtoMx Spatial Informatics Platform (SIP) and perform quality-control checks, data analysis, and generate analysis plots.


**Day 0: Prepare Reagents and Bake Tissue Overnight**

- Prepare shelf stable reagents



- Bake slides overnight at 60°C to improve tissue adherence


**Day 1: Slide Preparation**



x2 x2

**Remove Paraffin**


- Xylene and EtOH washes
- Dry Slides 5 min at 60°C



x2


**Target Retrieval**

- 15 mins at 100°C
- H<sub>2</sub>O rinse
- FFPE human brain only: 70% EtOH for 2 hours
- 3 min EtOH wash
- Dry for 30 min-1 hour



**Protease Digestion**

- Apply incubation frame
- Apply digestion buffer
- Incubate at 40°C for 30 mins
- Wash 2X in DEPC treated water



**Apply Fiducials**

- Prepare and apply fiducials
- Incubate for 5 minutes
- Wash with 1X PBS




**Post-Fix Tissue**

- Wash in 10% NBF for 1 min.
- 2 washes, 5 mins each, of NBF Stop Buffer
- Wash 5 min in 1X PBS



**Blocking**

- Prepare and apply NHS-Acetate
- Incubate for 15 mins.
- 2 washes, 5 mins each, in 2X SSC



**Overnight Hybridization**


- Prepare and apply assay specific probes
- Incubate at 37°C overnight



**Add Buffer 4 Enzymes**

- Add P2OX and Catalase to Buffer 4 and leave sealed on bench until instrument loading.

**Day 2: Wash and Stain Slide**



x2

**Stringent Washes**

- 2 stringent washes, 25 mins each.
- 2 washes, 2 mins each, 2X SSC



**Blocking & Nuclear Stain**

- Prepare Nuclear Stain stock and apply to tissue
- Incubate 15 mins at RT
- Wash in 1X PBS for 5 mins




x3

**Segmentation Markers**


- Prepare Segmentation mix and apply to tissue
- Incubate 1 hour at RT
- 3 washes, 5 mins each, 1X PBS

**Prepare Flow Cells and Load Instrument**



**Prepare Flow Cells**


- Use the flow cell assembly tool to assemble the flow cells.



**Load Instrument**

- Follow on-instrument prompts to load reagents and flow cells.
- Begin instrument run.

**Analyze Data in AtoMx SIP**



## Equipment, Materials, and Reagents

The following equipment ([Table 1](#)), materials ([Table 2](#)), and reagents ([Table 3](#)) are required for this protocol but are not supplied by Bruker.





Table 1: Equipment not supplied by Bruker

	Equipment	Source	Part Number
<input type="checkbox"/>	<b>Baking oven</b>	Quincy Lab, Inc. <sup>®</sup> (or comparable)	<a href="#">Example</a>
<input type="checkbox"/>	<b>Hybridization oven including hybridization chamber</b> <ul style="list-style-type: none"> <li>RapidFISH Slide Hybridizer or</li> <li>HybEZ™ oven with humidity control tray</li> </ul> <b>NOTE:</b> These hybridization ovens are designed to keep the slides hydrated and maintain a precise temperature overnight. Bruker does not recommend the use of any other hybridization ovens for CosMx SMI slide preparation.	Boekel Scientific <sup>®</sup> or ACDBio™	<a href="#">240200</a> for 120V or <a href="#">240200-2</a> for 230V <a href="#">321710/321720</a> <a href="#">310012</a>
<input type="checkbox"/>	<b>6-quart pressure cooker</b> <p><b>NOTE:</b> Products from other vendors may require testing and optimization.</p> <p><b>OPTIONAL:</b> A steamer may be used in lieu of a pressure cooker and may be preferred for more fragile tissues. If a steamer is used, a thermometer is also needed.</p>	BioSB <sup>®</sup> TintoRetriever  Nesco <sup>®</sup> Hamilton Beach <sup>®</sup> Ovente <sup>®</sup> (not validated)	<a href="#">BSB 7015</a>  <a href="#">ST-25F</a> <a href="#">37530MN</a> <a href="#">FS62S</a>
<input type="checkbox"/>	<b>Ultrasonic bath</b> 500 mL capacity, 40kHz frequency with timer	General Lab Supplier	<a href="#">Example (CPX-952-118R)</a>
<input type="checkbox"/>	<b>Vortex mixer</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>Microcentrifuge</b> for 1.5 mL microcentrifuge tubes and 8-well PCR strip tubes	General Lab Supplier	Various
<input type="checkbox"/>	<b>Water bath</b> Temperature setting of 37°C	General Lab Supplier	Various
<input type="checkbox"/>	<b>Thermal cycler</b> Must include a 96-well 200 µL tube block	General Lab Supplier	Various
<input type="checkbox"/>	<b>Analytic scale with draft shield</b> <p><b>NOTE:</b> Ensure scale can measure in milligrams (mg) so that reagents can be weighed accurately.</p>	Various	<a href="#">Example</a>

Table 2: Materials not supplied by Bruker

	Material	Source	Part Number
<input type="checkbox"/>	<b>Pipettes for 2.0 – 1,000 µL</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>Filter tips</b> (RNase/DNase Free)	General Lab Supplier	Various
<input type="checkbox"/>	<b>2.0 mL Centrifuge tubes</b> (RNase/DNase Free)	General Lab Supplier	Various
<input type="checkbox"/>	<b>0.2 mL PCR tubes</b> or PCR strip tubes	General Lab Supplier	Various
<input type="checkbox"/>	<b>Leica BOND Plus slides or VWR Superfrost Plus Micro Slide, Premium</b> <b>NOTE:</b> These slides have been internally validated. Do not use other products. <b>NOTE:</b> Leica BOND Plus slides are preferred for tissue sections prone to peeling.	Leica Biosystems VWR	<a href="#">S21.2113.A</a> <a href="#">48311-703</a>
<input type="checkbox"/>	<b>Slide rack</b>	General Lab Supplier	<a href="#">Example</a>
<input type="checkbox"/>	<b>Polypropylene slide staining jars</b> (24 required) or <b>Slide staining station</b> <b>NOTE:</b> Due to the photo-sensitivity of this assay, the staining jars should be impermeable to light.	Ted Pella®  Amazon® Fisher Scientific	<a href="#">21029</a>  <a href="#">MH-SJ6302</a> <a href="#">NC1862866</a>
<input type="checkbox"/>	<b>Forceps</b> (for slide handling)	General Lab Supplier	Various
<input type="checkbox"/>	<b>Razor blades</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>Timer</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>RNase AWAY™</b> <b>NOTE:</b> RNaseZAP™ and other alternatives cannot be used as substitutes as they do not adequately remove both nucleic acid and nuclease contaminants.	ThermoFisher Scientific	<a href="#">7003PK</a>
<input type="checkbox"/>	<b>Kimwipes®</b> (large and small)	Various	Various
<input type="checkbox"/>	<b>StainTray slide staining system with black lid</b>	Sigma-Aldrich®	<a href="#">Example</a>
<input type="checkbox"/>	<b>VWR® polyethylene slide holder</b> Optional - used during tissue sectioning	VWR	<a href="#">82024-524</a>

Table 3: Reagents not supplied by Bruker

	Reagent	Source/ Part Number	Storage Conditions
<input type="checkbox"/>	<b>DEPC-treated water</b>	ThermoFisher Scientific, <a href="#">AM9922</a> (or comparable)	Room temperature
<input type="checkbox"/>	<b>100% Ethanol</b> (EtOH): ACS grade or Better 	General Lab Supplier	Flammable storage Room temperature
<input type="checkbox"/>	<b>10X Phosphate Buffered Saline</b> pH 7.4 (PBS)	ThermoFisher Scientific, <a href="#">AM9625</a> (or comparable)	Room temperature
<input type="checkbox"/>	<b>(Only for FFPE human brain)</b> <b>10X PBS with Tween-20 (0.5%)</b> (PBS-T)	ThermoFisher Scientific, <a href="#">J63596.K2</a> (or comparable)	Room temperature
<input type="checkbox"/>	<b>Xylene</b>  <b>NOTE:</b> Citrisolv can be used, however, follow the alternative workflow for Citrisolv for <a href="#">on page 25</a> .	General Lab Supplier	Flammable storage Room temperature
<input type="checkbox"/>	<b>20X SSC</b> (DNase, RNase free)	ThermoFisher Scientific, <a href="#">AM9763</a>	Room temperature
<input type="checkbox"/>	<b>Tris Base</b>	Sigma-Aldrich, <a href="#">10708976001</a> (or comparable)	Room temperature
<input type="checkbox"/>	<b>Glycine</b>	Sigma-Aldrich, <a href="#">G7126</a> (or comparable)	Room temperature
<input type="checkbox"/>	<b>Sulfo-NHS-Acetate powder</b> <b>NOTE:</b> Sulfo-NHS-Acetate powder is shipped in a plastic bag with a desiccant and should be left in the bag and stored at -20°C until ready to use.	Fisher Scientific™, <a href="#">26777</a>	-20°C
<input type="checkbox"/>	<b>10% Neutral Buffered Formalin</b> (NBF) 	EMS Diasum®, <a href="#">15740</a> (or comparable)	Room temperature
<input type="checkbox"/>	<b>100% Deionized formamide</b>  <b>NOTE:</b> Deionized formamide is optimal, however, formamide that is not deionized may also be used.	ThermoFisher Scientific, <a href="#">AM9342</a> or VWR, <a href="#">VWRV0606</a> (or comparable)	4°C (bring to RT for at least 30 minutes before opening)


**OPTIONAL:** Large volume stock solutions (>500 mL) of deionized formamide can be aliquoted and stored, protected from light, at 4°C (recommend 40 mL aliquots in 50 mL conical tubes). This will save time during Day 2 slide preparation. Bring stock to room temperature for at least 30 minutes before opening.

## Bruker Supplied Reagents

### CosMx FFPE Slide Preparation Kit (RNA)

Table 4: FFPE Slide Preparation Kit (Box 1 of 2)

#### Kit Contents (Box 1 of 2, store at 4°C)



10X Target Retrieval Buffer	Buffer R
CosMx RNA Blocking Buffer	NHS-Acetate Buffer
2X SSC-T	CosMx Fiducials
Incubation Frames (not pictured)	Incubation Frame Covers (not pictured)

Table 5: FFPE Slide Preparation Kit (Box 2 of 2)

#### Kit Contents (Box 2 of 2, store at -20°C)

CosMx Proteinase K



### CosMx RNase Inhibitor

Table 6: CosMx RNase Inhibitor

#### Kit Contents (store at -20°C)

RNase Inhibitor is **sold separately** and is used in RNA hybridization and instrument loading.



**CosMx RNA Panels** (see [Panel and Cell Segmentation Marker Selection on page 10](#))



Table 7: CosMx RNA Panels

CosMx RNA Panels (store at -20°C)	
Kit Name	Kit Component
Human Whole Transcriptome Panel (WTX), 19K-plex, RNA	CosMx Hs WTX RNA Probe Mix
Human 6K Discovery Panel, 6K-plex, RNA	CosMx Hs 6K Discovery RNA Probe Mix
Human Universal Cell Characterization (UCC) Panel, 1000-plex, RNA	CosMx Hs UCC RNA Probe Mix CosMx Hs UCC RNA Add-On Custom RNA Add-On replaces off-the-shelf RNA Add-On
Human Immuno-oncology Panel, 100-plex, RNA	CosMx Hs IO RNA Probe Mix CosMx Hs IO RNA Add-On Custom RNA Add-On replaces off-the-shelf RNA Add-On
Mouse Neuroscience Panel, 1000-plex, RNA	CosMx Mm Neuro RNA Probe Mix CosMx Mm Neuro RNA Add-On Custom RNA Add-On replaces off-the-shelf RNA Add-On
Mouse Universal Cell Characterization (UCC) Panel, 1000-plex, RNA	CosMx Mm UCC RNA Probe Mix CosMx Mm UCC RNA Add-On Custom RNA Add-On replaces off-the-shelf RNA Add-On

## CosMx Segmentation Markers




Table 8: CosMx SMI Cell Segmentation and Supplemental Marker Kits


CosMx RNA Segmentation Marker Kits (each kit sufficient for up to 4 slides, store at -80°C)	
Kit Description	Kit Components
Human Universal Cell Segmentation Kit (RNA) Compatible with: Human WTX Panel, Human 6K Discovery Panel, Human Immuno-Oncology 100-plex Panel, and any 1000-plex RNA Panel	CosMx DAPI Nuclear Stain, Ch1
	CosMx Hs CD298/B2M Segmentation Marker Mix, Ch2 (RNA)
	CosMx Hs PanCK/CD45 Marker Mix Ch3/Ch4 (RNA)
Human Neuroscience Cell Segmentation Kit (RNA) Compatible with: Human WTX Panel, Human 6K Discovery Panel	CosMx DAPI Nuclear Stain, Ch1
	CosMx Hs Neuro rRNA Neuro Marker, Ch2 (RNA)
	CosMx Mm/Hs Neuro Histone Marker, Ch3 (RNA)
	CosMx Mm/Hs GFAP Marker, Ch4 (RNA)
Mouse Neuroscience Cell Segmentation Kit (RNA) Compatible with: Mouse Neuroscience 1000-plex Panel	CosMx DAPI Nuclear Stain, Ch1
	CosMx Mm Neuro rRNA Neuro Marker, Ch2 (RNA)
	CosMx Mm/Hs Neuro Histone Marker, Ch3 (RNA)
	CosMx Mm/Hs GFAP Marker, Ch4 (RNA)
Mouse Universal Cell Segmentation Kit (RNA) Compatible with: Mouse Universal Cell Characterization 1000-plex Panel	CosMx DAPI Nuclear Stain, Ch1
	CosMx Mm CD298/B2M Marker Mix, Ch2 (RNA)
	CosMx Mm PanCK/CD45 Marker Mix Ch3/Ch4 (RNA)

The following markers are **optional and available to order à la carte**:

Table 9: A La Carte markers

Compatible with	Item Description (store at -80°C)
Human Universal Cell Segmentation Kit (RNA)	CosMx Hs CD68 A La Carte Marker, Ch5 (RNA)
	CosMx Hs Cytokeratin 8/18 A La Carte Marker, Ch5 (RNA)
	CosMx Hs/Mm CD3 A La Carte Marker, Ch5 (RNA)
Mouse Universal Cell Segmentation Kit (RNA)	CosMx Mm CD68 A La Carte Marker, Ch5 (RNA)
	CosMx Mm CD8 A La Carte Marker, Ch5 (RNA)
	CosMx Hs/Mm CD3 A La Carte Marker, Ch5 (RNA)

**Flow Cell Assembly Tool and Kit**



The Flow Cell Assembly tool is a one-time required purchase.

The Flow Cell Assembly Kit contains 4 single-use Flow Cell coverslips, sufficient for a 4-slide experiment.

## Tissue Sectioning and Slide Preparation: FFPE Samples

[Appendix I: CosMx SMI Tissue Sectioning Guidelines on page 84](#) covers FFPE block selection and sectioning in detail. Review these guidelines prior to beginning FFPE Slide Preparation.

Bruker has tested and validated sample blocks up to 3 years old prepared from tissues with a cold ischemic time of less than 1 hour using 10% NBF or similar fixative. For best results, do not use FFPE blocks older than 10 years. Assay performance will be influenced by tissue block age and treatment conditions such as cold/warm ischemic time, fixative, and storage.

Tissues used for CosMx validation testing measure approximately 1.0 x 1.0 x 0.4 mm and are fixed for 24-36 hours before embedding.

FFPE blocks should be sectioned at **5  $\mu\text{m}$  thickness** and mounted on the label side of Leica BOND PLUS slides or VWR Superfrost Plus Micro Slides.

Tissue sections must be centered within the Scan Area of the slide (the green rectangle in [Figure 3](#)) and be no larger than **20 mm long by 15 mm wide** (image not to scale; see the Flow Cell Assembly Tool for a to-scale template). For best performance, there should also be some bare glass (not covered by tissue) in the Scan Area. For examples of tissue placement best practices, see [Appendix I: CosMx SMI Tissue Sectioning Guidelines on page 84](#).

Label slides with pencil on the frosted label according to lab guidelines. If using an adhesive slide label, ensure the label is less than 295  $\mu\text{m}$  thick and is not folded over on itself. Labels over the maximum thickness or labels that are not properly adhered may result in slide or flow cell damage during flow cell assembly and/or instrument loading.




Figure 3: Tissue Scan Area (not to scale)

If sections are larger than the allowable size, or placed off-center, excess tissue must be scraped away with a clean razor blade. Tissue may be scraped when slides are dry at the start of slide preparation, before applying the incubation frame during slide preparation, or just before applying the flow cell coverslip after slide preparation. The optimal practice (minimizing tissue pulling, folding, or detaching) may depend on your sample type. Applying the incubation frame over tissue could result in tissue damage when the frame is removed or poor sealing of the incubation frame.

**IMPORTANT:** The CosMx SMI instrument will only image the area inside the flow cell chamber, the tissue scan area. If the tissue section is outside of the scan area, it will not be imaged.

After sectioning and prior to use or storage, to improve tissue adherence, **bake slides at 37°C overnight** at an angle no greater than 45 degrees. Alternatively, slides can be baked at 37°C for 2 hours and then dried overnight at room temperature. A polyethylene slide holder ([VWR, 82024-524](#)) can be used for overnight drying. Ensure sections are completely dry before storage.

A tissue section adhesive such as EpreDia™ Tissue Section Adhesive (Fisher Scientific, [86014](#)) can also be used to improve tissue adhesion. The use of an adhesive has not been validated but may improve tissue adherence for some tissue types. Follow the manufacturer's instructions for use.

It is recommended to use mounted sections within two weeks for best results. Older sections (1-2 months) may produce reasonable results, but results may be tissue or block dependent and should be tested empirically. Slides should be stored in a desiccator, at room temperature or 4°C prior to processing.

## Prepare RNA FFPE Assay Reagents

**i IMPORTANT:** Take care to maintain nuclease-free conditions. The greatest risk of contamination comes from CosMx SMI RNA probes and other oligos. Bruker recommends the use of RNase AWAY (ThermoFisher [7003PK](#)), as it will limit contamination from oligos, detection probes, and nucleases. After using RNase AWAY, air dry completely, or rinse with DEPC-treated water. See manufacturer's instructions for details.

Label staining jars and prepare reagents using the instructions in the following table ([Table 10](#)). Ensure you have sufficient buffer to completely cover the tissue on all slides without submerging the slide labels. Contact with buffer may make slide labels illegible.

Unless otherwise noted, reagents can be made up to 2 weeks in advance and stored at room temperature.

Table 10: RNA Reagent Preparation


	Reagent	Dilution	Storage
<input type="checkbox"/>	1X PBS (pH 7.4)	Prepare 1 L of 1X PBS by combining 100 mL of 10X PBS and 900 mL of DEPC-treated water.	Room temperature
<input type="checkbox"/>	1X PBS-T (0.5% Tween-20)	<b>(Only for FFPE human brain tissue)</b> Prepare 1 L of 1X PBS-T by combining 100 mL of 10X PBS-T and 900 mL of DEPC-treated water.	Room temperature
<input type="checkbox"/>	70% Ethanol	<b>(Only for FFPE human brain tissue)</b> Prepare 500 mL of 70% ethanol by combining 150 mL of DEPC-treated water and 350 mL 100% ethanol.	Room temperature
<input type="checkbox"/>	2X SSC	Prepare 1 L of 2X SSC by combining 100 mL of 20X SSC and 900 mL of DEPC-treated water.	Room temperature
<input type="checkbox"/>	4X SSC	Prepare 500 mL of 4X SSC by combining 100 mL of 20X SSC and 400 mL of DEPC-treated water.	Room temperature
<input type="checkbox"/>	NBF Stop Buffer (Tris Glycine Buffer)	Combine 6.06 g Tris base and 3.75 g Glycine in 500 mL of DEPC-treated water.  The final concentration of Tris and Glycine will be 0.1 M each.	Room temperature
<input type="checkbox"/>	NHS-Acetate Mix	Individual aliquots can be prepared prior to slide preparation but should be stored at -20°C with a desiccant until use.  To prepare aliquots for 4 slides: <ul style="list-style-type: none"> <li>• Bring stock to room temperature for 1-2 hours before opening.</li> <li>• Prelabel four 2.0 mL screw top centrifuge tubes.</li> <li>• Using a weighing spatula, carefully weigh 25 mg of Sulfo-NHS-Acetate directly into one screw top tube on an analytic scale.</li> <li>• Close the tube and label with final weight. Seal the tube with parafilm and store at -20°C with desiccant until use.</li> </ul>	-20°C with desiccant
<input type="checkbox"/>	1X Target Retrieval Solution	Prepare 50 mL by adding 5 mL of CosMx Target Retrieval Buffer (10X) (Slide Prep Kit, Box 1) to 45 mL DEPC-treated water. Do not store or reuse.	Make fresh each time

	Reagent	Dilution	Storage
<input type="checkbox"/>	Digestion Buffer	<p>Dilute the 20 mg/mL Proteinase K stock (provided by Bruker) to a working concentration of 3 µg/mL* in <b>1X PBS-T for FFPE human brain</b> or <b>1X PBS for all other tissues</b>. Prepare fresh each time.</p> <p>A 2-step serial dilution is recommended. Accurate dilution of Proteinase K is critical for proper assay performance.</p> <ul style="list-style-type: none"> <li>Step 1: Dilute 20 mg/mL stock to 200 µg/mL by adding 2 µL of stock to 198 µL of 1X PBS-T for FFPE human brain or 1X PBS for all other tissues.</li> <li>Step 2: Dilute to the target concentration of 3 µg/mL by adding 30 µL of the 200 µg/mL solution to 1970 µL of 1X PBS-T for FFPE human brain or 1X PBS for all other tissues.</li> <li>Mix thoroughly by flicking tube or pipetting up and down using a clean tip. <b>Do not vortex.</b></li> </ul> <p><b>NOTE:</b> Digestion Buffer should be prepared fresh each time and stored on ice until ready to use.</p> <p>* For cell pellet arrays, prepare Digestion Buffer at 1 µg/mL. If FFPE human brain tissue exhibits high background or autofluorescence after digestion with 3 µg/mL Proteinase K, increase to 5 µg/mL, as long as tissue does not peel away from slide. Optimal digestion buffer concentration may need to be empirically determined. See <a href="#">Appendix II: Tissue-Specific Slide Preparation Considerations on page 88</a>).</p>	<p>Make fresh each time</p> <p>Store on ice</p>
<input type="checkbox"/>	Fiducials and Hybridization Mix	<p>On Day 1, remove RNA Probe Mix from -20°C and keep on ice until ready to use.</p> <p><b>NOTE:</b> Fiducials are light sensitive and should be kept stored, protected from light, until instructed to remove later in this protocol.</p>	n/a
<input type="checkbox"/>	Day 2 Reagents	<p>These reagents have additional steps that will be covered in detail in their respective sections in the Day 2 protocol.</p>	n/a

**FFPE** Day 0: Prepare Shelf Stable Reagents and Bake Slides Overnight**Day 0: Prepare Shelf Stable Reagents and Bake Slides Overnight**

- Prepare shelf-stable reagents following [Prepare RNA FFPE Assay Reagents on page 22](#).
- In a slide rack, arrange slides vertically or in a slide holder at a 45°-angle and bake sections on slides overnight in a 60°C drying oven.

Overnight tissue baking at 60°C has been shown to increase tissue adherence and stability during slide preparation. If slides cannot be baked overnight, then bake for at least 2 hours prior to [Day 1: Deparaffinize FFPE Tissue Sections \(20 minutes\) on page 25](#). Please note that CosMx validation testing includes overnight baking.

 **IMPORTANT:** If using a custom panel, email [AtoMx.KitAdmin@Bruker.com](mailto:AtoMx.KitAdmin@Bruker.com) at **least one business day prior to a run** so that the custom kit is available in the CosMx SMI Control Center.

## Day 1: Deparaffinize FFPE Tissue Sections (20 minutes)

You will need the following materials and reagents for this step: **staining jars**, **xylene** and **100% ethanol**. The **hybridization tray**, **pressure cooker or steamer**, **1X Target Retrieval Solution**, and **DEPC-treated water** are preheated here for their use in a later step.



Before beginning, slides should have already been baked overnight at 60°C (see [Day 0: Prepare Shelf Stable Reagents and Bake Slides Overnight on page 24](#)). If slides were not baked overnight, bake for 2 hours at 60°C at the start of Day 1. Overnight baking is optimal.

### Prepare Equipment and Washes

1. Prepare hybridization tray by lining the bottom of the tray with Kimwipes and carefully wet the Kimwipes with 2X SSC or DEPC-treated water. Kimwipes should be thoroughly damp but without standing buffer.
2. Preheat hybridization chamber and tray to 40°C following manufacturer's instructions.
3. Ensure baking oven is set to 60°C.
4. Prepare staining jars with enough xylene and ethanol to cover tissue ([Figure 4](#)). Ensure you have sufficient buffer to completely cover the tissue on all slides without submerging the slide labels. Contact with buffer may make slide labels illegible.

**WARNING:** Xylene and ethanol are flammable chemicals and should be handled appropriately. Properly dispose of waste generated in these steps following lab guidelines.

5. Prepare the pressure cooker or steamer and preheat the target retrieval solution as follows. Content in **purple boxes** denotes instructions specific to the pressure cooker. Content in **orange boxes** denotes instructions specific to the steamer. Use the same target retrieval method (pressure cooker or steamer) throughout the study.

**WARNING:** Bruker does not recommend the use of glass staining jars in the pressure cooker.

#### Pressure Cooker Method:

1. Fill the pressure cooker with water to the correct level per the manufacturer's instructions (4-8 cups depending on model used).
2. Place the staining jar containing freshly prepared 1X Target Retrieval Solution into the pressure cooker to preheat. Ensure that water level is well below lid of jar; about halfway up jar is sufficient. If needed, a trivet may be used to raise the jar.
3. Place a lid on the staining jar to prevent evaporation. To prevent pressure from building within the staining jar, **do not fully seal the lid to the jar**.
4. Preheat the pressure cooker to 100°C following the model-specific instructions below. Pressure cooker preheating takes about 1 hour.

**BioSB Preheating Instructions:**

Use the *TintoRetriever Pressure Cooker Preheating Cycle QuickStart Guide, Two Staining Dish Operation* provided with the BioSB pressure cooker to preheat the pressure cooker.

- With the pressure valve closed, press the **80°C** button on the face of the pressure cooker and press **Start** to run a cycle at 80°C with a 0-minute timer.



- Once the first cycle is complete, run a second cycle with a **45-minute timer at 100°C**.

**Steamer Method**

1. Fill the steamer reservoir up to the fill line with water.
2. Place two staining jars inside of the steamer, one containing DEPC-treated water and one containing 1X Target Retrieval Solution. Ensure sufficient reagent volume to cover tissue, up to the slide label.
3. Loosely cover each jar with aluminum foil instead of the jar lid to allow for a thermometer reading in a later step.
4. Preheat the steamer to 100°C. More water may need to be added to the steamer during preheating.

The steamer may take up to 1 hour to heat the liquid in the jars to a stable maximum temperature near 100°C. Final temperature can be checked by inserting a digital thermometer through the hole in the lid of the steamer into the DEPC-water staining jar.

**Deparaffinize FFPE tissue sections**

1. Remove the slides from the baking oven and gently perform the following washes using staining jars ([Figure 4](#)).

Ensure you have sufficient buffer to completely cover the tissue on all slides without submerging the slide labels. Contact with buffer may make slide labels illegible.

- Place slides into first xylene jar and wash for 5 minutes.
- Transfer slides to new staining jar of xylene and wash for 5 minutes.
- Transfer slides to 100% ethanol and wash for 2 minutes (or 5 minutes if CitriSolv was used in lieu of xylene).

- Transfer slides to new staining jar of ethanol and wash for 2 minutes (or 5 minutes if CitriSolv was used in lieu of xylene).




Figure 4: Deparaffinize FFPE Tissue Sections

**WARNING:** Xylene and ethanol are flammable chemicals and should be handled appropriately. Properly dispose of waste generated in these steps following lab guidelines.

2. **Dry slides** in slide rack in 60°C oven for **5 minutes**. After 5 minutes, remove slides from oven and leave at room temperature until target retrieval solution is pre-heated to 100°C.
3. **If not already aliquoted, remove stock NHS-Acetate powder from -20°C freezer** and leave at room temperature for 1-2 hours before use in [NHS-Acetate Preparation and Application \(25 minutes\) on page 36](#). Pre-aliquoted NHS-Acetate powder can remain at -20°C until instructed to remove [on page 36](#).

## Perform Target Retrieval (50 minutes)

You will need the following materials and reagents for this step: **pressure cooker or steamer, staining jars, 1X Target Retrieval Solution (pre-heated in the previous step), DEPC-treated water and 100% ethanol.**

Target retrieval times were determined based on FFPE tissue blocks meeting the constraints outlined in the sample guidance section. These conditions may vary by sample, the amount of normal adjacent tissue, and other factors. See [Appendix II: Tissue-Specific Slide Preparation Considerations on page 88](#).

Content in **purple box** denotes instructions specific to the pressure cooker. The **orange box** denotes instructions specific to the steamer. Use the same target retrieval method (pressure cooker *or* steamer) throughout the study.

### Pressure Cooker Method

- Once Target Retrieval Solution is preheated, press Cancel on the pressure cooker, release the pressure valve, and wait for the pressure cooker to release pressure. Once pressure has released, remove the lid and **carefully but quickly remove the staining jar containing preheated 1X Target Retrieval Solution.** Once removed, the Target Retrieval Solution will begin to rapidly cool. Ensure the following steps are done as quickly and safely as possible.




**WARNING:** When opening the pressure valve and removing the pressure cooker lid, hot steam will be released. The staining jar will also be VERY hot. Wear protective heat-resistant gloves to open lid and remove staining jar.

- Place FFPE slides into the preheated solution and replace lid on the staining jar to prevent evaporation. To prevent pressure from building within the container, **do not fully seal the lid to the jar.**
- Return the staining jar containing the slides into the preheated pressure cooker. Ensure that water level is well below lid of jar. About halfway up jar is sufficient. If needed, a trivet may be used to raise the jar.
- Reattach the pressure cooker lid, **open the pressure release valve to Pressure Release position** and return the pressure cooker to 100°C. For the BioSB, this can take up to 20



### Steamer Method

- Without removing the lid, place an instant-read digital thermometer through the vents in the steamer lid and pierce the aluminum foil covering the DEPC-treated water. Ensure the water has reached about 99°C. Add more water as needed.
- Once the water has reached 99°C, carefully remove the steamer lid. Once removed, the Target Retrieval Solution will begin to rapidly cool. Ensure the following steps are done as quickly and safely as possible.



**WARNING:** Removing the steamer lid releases hot steam. Wear protective heat-resistant gloves to open lid and remove the staining jar. Transfer slides using forceps or rack.

- Remove the foil from the Target Retrieval Solution jar and quickly transfer the slides to the Target Retrieval Solution. Replace the foil, then replace steamer lid.
- Reinsert the thermometer into the DEPC-water jar and wait until the temperature returns to about 99°C.
- Once the steamer temperature returns to 99°C, start timer and run for 15 minutes for FFPE tissue or 8 minutes for cell pellet arrays**

minutes.

**NOTE:** If after 10 minutes the temperature has not started to increase it may indicate that the float valve is not properly seated. Without removing the lid, turn the lid to the unlock position and then slowly back to lock.

5. **Once the pressure cooker temperature returns to 100°C**, start timer and run for **15 minutes** for FFPE tissue or 8 minutes for cell pellet arrays (CPA).
6. When the timer reaches zero, press Cancel on pressure cooker to stop heating, carefully remove the pressure cooker lid, and remove the staining jar.

(CPA).

6. When the timer reaches zero, carefully remove the steamer lid and remove the staining jar.

7. **Immediately transfer all slides to staining jar containing fresh DEPC-treated water.** Move slides up and down for 15 seconds to wash ([Figure 5](#)).

Transfer slides from Target Retrieval Solution into DEPC-treated water carefully but quickly to ensure slides do not dry out. Drying can occur quickly when slides are hot.

8. **Only for FFPE human brain tissue:** Transfer slides to fresh 70% ethanol for 2 hr (up to overnight) at 4°C. For all other tissue types, slides should go directly from DEPC-treated water to 100% ethanol (step 9).
9. Transfer slides to fresh 100% ethanol and incubate for 3 minutes.
10. During ethanol wash, clean bench space with RNase AWAY and lay out a fresh Kimwipe.
11. After 3 minutes, remove slides from ethanol and lay horizontally on a clean Kimwipe. Dry at room temperature for 30 minutes to 1 hour.




Figure 5: Water and Ethanol Wash

12. While slides are drying, continue to next page to prepare incubation frames and digestion buffer.

**FFPE Tissue Permeabilization (40 minutes)****Tissue Permeabilization (40 minutes)**

You will need the following materials and reagents for this step: **staining jars, hybridization oven, hybridization tray, incubation frames** (from Slide Prep Kit Box 1), **digestion buffer** (see [Prepare RNA FFPE Assay Reagents on page 22](#)), and **DEPC-treated water**.

1. If needed, trim the tissue following the template in [Tissue Sectioning and Slide Preparation: FFPE Samples on page 20](#).
  - Use a clean razor blade to trim tissue and change blade as needed to ensure clean cuts and reduce the risk of cross-contamination between samples.
2. Prepare the Incubation Frame:
  - Separate an individual frame from the strip by carefully tearing along the perforations.
  - Each frame is sandwiched between a thin solid polyester sheet and a thick polyester frame backing (with the center square removed).
3. Using a clean Kimwipe, ensure that the surface of the slide that will come in contact with the incubation frame is clean and dry.
4. Apply the incubation frame ([Figure 6](#)).
  - Carefully remove the **thin polyester sheet**, ensuring that the frame remains bound to the thick polyester frame backing (with the center square removed).
  - With the slide on a flat surface, careful not to touch the adhesive, **center the tissue within the incubation frame and carefully place the incubation frame** around each tissue section. Lightly press along the border of the incubation frame to ensure that it is well adhered to the slide.




Figure 6: Apply Incubation Frame

5. With the slide still on a flat surface, use a clean razor blade to carefully trim the long edges of the incubation frame to remove excess plastic ensuring that there is no excess film extending over slide edges. Trim the short end of the incubation frame (opposite the slide label) as needed.
6. Remove digestion buffer from ice and warm digestion buffer by hand for about 3 minutes to bring the mixture to room temperature.
7. Retrieve preheated hybridization tray from hybridization oven.

8. Place slides into slide insert of hybridization tray and, using a P200 pipette, slowly add 400  $\mu\text{L}$  (2 x 200  $\mu\text{L}$ ) of digestion buffer to the tissue within incubation frame ([Figure 7](#)). Gently move tray side to side as needed to **ensure that digestion buffer covers the entire tissue**.



Figure 7: Hybridization Tray

A pipette tip can also be used to spread buffer over tissue within incubation frame. Use a small-volume pipette tip and carefully lay the tip horizontally on top of the incubation frame. Gently roll the tip to spread the buffer, until the tissue is completely covered. Be careful to avoid touching the tissue with the pipette tip.

9. Insert hybridization tray containing slides into hybridization oven and incubate at 40°C following the guidance in the table below ([Table 11](#)) ([Figure 8](#)).




Figure 8: Hybridization Oven

Table 11: Digestion Times

Tissue Type	Digestion Time
FFPE Tissue	30 minutes
Cell Pellet Array (CPA)	15 minutes

Incubation times and temperatures may differ by tissue and may need to be empirically determined. **Start with the default conditions of 30 minutes at 40°C** and adjust the time and concentration if needed. See [Appendix II: Tissue-Specific Slide Preparation Considerations on page 88](#) for additional context for specific tissue types, noting that the default conditions are recommended for almost all tissue types.

Decreasing digestion buffer concentration and/or incubation time may increase tissue stability for certain tissue types.

**i IMPORTANT:** After digestion buffer has been applied, avoid tissue drying in subsequent steps by working with only one slide at a time.

**FFPE Tissue Permeabilization (40 minutes)**

10. During protease incubation, remove fiducials and 2X SSC-T from 4°C and let come to room temperature, protected from light, for at least 10 minutes. Then, the Fiducial Preparation steps can be started in parallel with the below wash steps.
11. After slide incubation, one slide at a time, tap off excess digestion buffer and transfer each slide to staining jar containing fresh DEPC-treated water.
12. Move slides up and down 3-5 times to wash. Transfer slides to a new jar of DEPC-treated water and repeat wash ([Figure 9](#)).
13. Slides can be stored in DEPC-treated water while fiducials are prepared.




Figure 9: DEPC-treated water wash

Limit the time that the slides are kept in DEPC-treated water to minimize risk of target loss. At this point, targets have been exposed and tissue has not yet undergone post-fixation.

## Fiducial Preparation and Application (20 minutes)

You will need the following materials and reagents for this step: **ultrasonic bath**, **vortex**, **2X SSC-T** (at room temperature), **fiducials** (at room temperature), **staining jars**, **staining tray** (clean and dry), and **1X PBS**.

The volume of fiducial working solution prepared here is 1.2 mL, sufficient for 4 slides. Do not make less working solution due to high risk of clumping. Additional volume may be made as needed.

1. Refer to [Figure 10](#) below to vortex and sonicate the fiducials prior to use.

**IMPORTANT:** When sonicating the fiducial tube, be sure not to submerge the tube cap under the liquid level in the sonicator. Use a floating tube holder if needed to float the tubes in the ultrasonic bath.




Figure 10: Fiducial preparation

Failure to follow these steps will result in fiducial clumping and uneven distribution of fiducials within the tissue. This uneven distribution can result in a loss of readable area or loss of image registration.

2. Once fiducials are prepared as in [Figure 10](#), dilute the fiducial stock to the appropriate working concentration:

**For FFPE human brain, make a 1:200 dilution** from stock (0.1%) to the working concentration (0.0005%) in 2X SSC-T by adding 6  $\mu\text{L}$  of the fiducial stock to 1194  $\mu\text{L}$  of 2X SSC-T (597  $\mu\text{L} \times 2$ ).

**For all other tissues, make a 1:100 dilution** from stock (0.1%) to the working concentration (0.001%) in 2X SSC-T by adding 12  $\mu\text{L}$  of the fiducial stock to 1188  $\mu\text{L}$  of 2X SSC-T (594  $\mu\text{L} \times 2$ ).

The optimal concentration may differ for some tissue types and may need to be empirically determined. See [Appendix II: Tissue-Specific Slide Preparation Considerations on page 88](#).

3. Remove slide from DEPC-treated water and gently tap slide on a clean Kimwipe to remove excess water. Lay slide horizontally in staining jar.
4. **Immediately before applying fiducials to slides, vortex fiducials working solution for 1 additional minute.** Vortex fiducials for 30 seconds between applications to slides to keep fiducials in suspension and ensure consistent concentration across all slides.
5. Apply up to 250  $\mu\text{L}$  of the final fiducial working solution, ensuring the solution covers glass and tissue within the incubation frame. Fiducials must be present on the glass within the scan area for consistent focusing during the

instrument run. Gently move tray side to side as needed to ensure that the fiducial solution covers the entire scan area, including glass.

A pipette tip can also be used to spread buffer over tissue within incubation frame. Use a small-volume pipette tip and carefully lay the tip horizontally on top of the incubation frame. Gently roll the tip to spread the buffer, until the tissue is completely covered. Be careful to avoid touching the tissue with the pipette tip.

6. Incubate covered in staining tray for 5 minutes at room temperature ([Figure 11](#)). During fiducial incubation, prepare staining jars for next step.



Figure 11: Incubate covered for 5 minutes

**i IMPORTANT:** This step and subsequent steps are light sensitive as fiducials contain fluorescent molecules that are sensitive to photo-bleaching.

7. After fiducial incubation, *one slide at a time*, gently tap slides on a clean Kimwipe to remove excess solution and transfer slides to staining jar containing fresh 1X PBS ([Figure 12](#)).
8. Wash slides in staining jar with 1X PBS for 1 minute. Proceed to next steps immediately.




Figure 12: PBS Wash (1 minute)

**Post-Fixation (20 minutes)**

You will need the following materials and reagents for this step: **staining jars, 10% NBF, NBF Stop Buffer (Tris-Glycine Buffer)** and **1X PBS** (see [Prepare RNA FFPE Assay Reagents on page 22](#) for more information).

**WARNING:** Use of appropriate personal protective equipment is advised. Used NBF Stop Buffer contains NBF and must be disposed of in the same manner as NBF.

Post-fix the tissue by performing the following washes ([Figure 13](#)).

1. Transfer slides to 10% NBF and incubate for 1 minute at room temperature.
2. **Immediately transfer the slides** to the first NBF Stop Buffer and wash for 5 minutes. Transfer slides to a second staining jar containing NBF Stop Buffer and wash for 5 minutes.
3. Transfer slides to 1X PBS for 5 minutes. Slides can sit in 1X PBS while NHS-Acetate mix is prepared.



Figure 13: NBF Post Fix

4. During PBS wash, remove RNase Inhibitor, CosMx RNA probe mix and add-on probes or stand-alone custom panel from -20°C and thaw on ice.
5. Remove Buffer R from 4°C and bring to room temperature.

## NHS-Acetate Preparation and Application (25 minutes)

You will need the following materials and reagents for this step: **analytical scale, 2.0 mL centrifuge tube, staining jar, Sulfo-NHS-Acetate powder, NHS-Acetate buffer** (Slide Prep Kit Box 1, 4°C), and **2X SSC**.

1. Prepare 100 mM NHS-Acetate mixture (200 µL/sample) **immediately before** you are ready to apply it to the tissue:
  - If not already done, prepare four 25-mg aliquots (each for 4-slide preparation) or six 15-mg aliquots (each for 2-slide preparation): Bring Sulfo-NHS-Acetate powder from -20°C to room temperature, to prevent condensation. Weigh out the powder directly into 2.0 mL centrifuge tubes. Label tubes with the exact weight, seal tubes with parafilm, and store at -20°C in desiccant. Refer to the Prepare Reagents section of this protocol.
  - Calculate the amount of NHS-Acetate buffer to add to the Sulfo-NHS-Acetate powder by multiplying the weight of Sulfo-NHS-Acetate powder in mg by 38.5. (For example, for 25.0 mg of Sulfo-NHS-Acetate powder,  $25.0 * 38.5 = 962.5$  µL of buffer to add.)
  - Add NHS-Acetate buffer directly to aliquoted powder immediately before applying to the tissue. Slowly pipette up and down to mix. Bubbles may form. Do not fully dispense liquid from pipette while mixing.

**IMPORTANT:** Reconstitute Sulfo-NHS-Acetate immediately before use. Do not store, as the NHS ester moiety readily hydrolyzes and becomes non-reactive. Discard any unused reconstituted reagent.

2. Perform the following steps **one slide at a time** to prevent the tissue from drying out:
  - Remove slide from 1X PBS, gently tap slide on a clean Kimwipe to remove excess buffer, and transfer to a clean staining tray.
  - Apply 200 µL of NHS-Acetate mixture onto the tissue within the incubation frame. Gently rock the tray side to side as needed to ensure that the NHS-Acetate solution covers the entire tissue.
  - Repeat with remaining slides and incubate covered in staining tray for 15 minutes at room temperature ([Figure 14](#)).



Figure 14:  
Incubate 15  
minutes

3. Following incubation, tap off excess liquid and wash slides in 2X SSC for 5 minutes ([Figure 15](#)).
4. Repeat 2X SSC wash for a total of 2 washes. Leave slides in second 2X SSC wash while hybridization mix is prepared.



Figure 15: Two 5-  
minute 2X SSC  
washes

## In Situ Hybridization (overnight)



**Do not begin in situ hybridization step until within 16-18 hours of Day 2 start time.** Until then, store slides protected from light in 2X SSC up to 1 hour at room temperature or up to 6 hours at 4°C.

**IMPORTANT:** Take care to maintain nuclease-free conditions. Areas should be cleaned thoroughly with RNase AWAY after probe mix formulation (RNaseZap is only effective for enzymes, not oligos, and should not be used in place of RNase AWAY). Alternatively, mixes can be made in PCR workstations that have been decontaminated with UV light. Gloves should also be changed after handling any probe mixes to avoid cross-contamination.

**IMPORTANT:** If preparing **human or mouse neural tissue**, rRNA is used in this step ([Table 13](#)). The rRNA marker is shipped and stored at -80°C with the Neuroscience Cell Segmentation Kit. No other components of the Kit are used in this step.

**IMPORTANT:** If using a **custom RNA Add-on**, the custom RNA Add-on replaces the standard panel Add-on for the 100- and 1000-plex assay. (The Human 6K Discovery and WTX Panels do not include an Add-on.)

**IMPORTANT:** If using a **custom panel**, email [AtoMx.KitAdmin@Bruker.com](mailto:AtoMx.KitAdmin@Bruker.com) **at least one business day prior to a run** so that the custom kit can be added to the CosMx SMI Control Center.

You will need the following materials and reagents for this step: **hybridization oven, hybridization tray, incubation frame covers, thermal cycler, ice, Buffer R, CosMx RNA Probe Mix, custom or standard Add-On (-20°C), rRNA Marker (if applicable), RNase Inhibitor (-20°C), and DEPC-treated water.**

**Prepare buffers:** Warm Buffer R to room temperature before opening. Thaw probe mix, add-on probes (if applicable) and rRNA Marker (if applicable) on ice. Before use, mix probes thoroughly by pipetting up and down 3-5 times. **Do not vortex probes.** Once thawed, probes can be refrozen at -20°C up to 5 times or refrigerated at 4°C for up to 6 months.

**Set the hybridization oven temperature to 37°C** according to product instructions. If your chamber is light-permeable, minimize light exposure (e.g., by wrapping the lid in aluminum foil).

1. Pre-heat thermal cycler and lid to 95°C.
2. Remove an incubation frame cover and clean with ethanol. Dry with a clean Kimwipe and visually inspect the cover for dust. Use a new Kimwipe as needed to remove any dust. Lay incubation frame cover on a clean Kimwipe until use.
3. For each probe mix, flick to mix then centrifuge. **Do not vortex probes.**
4. Denature CosMx RNA probe mixes (RNA Probe Mix, RNA Add-On or Custom Add-On, plus rRNA Marker if

**FFPE** *In Situ Hybridization (overnight)*

preparing neural tissue), by transferring total volumes needed for assay from stock tubes into clean 0.2 mL PCR tubes. **Probes, Add-ons, and rRNA Marker should be kept separate during denaturing.**

**IMPORTANT:** Pipette accurately. When preparing 4 slides, there is no excess Probe Mix or Add-on.

- Heat at 95°C for 2 minutes on a thermal cycler with heated lid. Immediately transfer to ice for at least 1 minute to crash cool.
- Immediately before preparing hybridization mix, flick to mix and centrifuge tubes.
- Make hybridization solution following [Table 12](#) for human or mouse **non-neural tissue** or [Table 13](#) for human or mouse **neural tissue**. **Prepare hybridization mix no more than 20 minutes before tissue application.**

Table 12: Hybridization Solution for **Non-Neural Tissue**

	Denatured RNA Probe Mix	Denatured Add-on* (if applicable)	RNase Inhibitor	Buffer R	DEPC-treated water	Final Volume
Off-the-shelf Panel (2 slides)	32 µL	16 µL	3.2 µL	256 µL	up to final volume of 320 µL	320 µL
Off-the-shelf Panel (4 slides)	64 µL	32 µL	6.4 µL	512 µL	up to final volume of 640 µL	640 µL
Custom Stand-alone Panel (2 slides)	16 µL**	-	3.2 µL	256 µL	up to final volume of 320 µL	320 µL
Custom Stand-alone Panel (4 slides)	32 µL**	-	6.4 µL	512 µL	up to final volume of 640 µL	640 µL

\* If using a 100- or 1000-plex assay with a custom Add-on, the custom Add-on replaces the standard Add-on that is part of the assay. (WTX and 6K Discovery Panels do not include a standard Add-on.) If no Add-on is used, use DEPC-water to reach the correct final volume.

\*\* Custom stand-alone panels are supplied at 2X concentration relative to off-the-shelf RNA Panels.

Table 13: Hybridization Solution for **Neural Tissue**

	Denatured RNA Probe Mix	Denatured Add-on * (if applicable)	Denatured rRNA Marker**	RNase Inhibitor	Buffer R	DEPC-treated water	Final Volume
Off-the-shelf Panel (2 slides)	32 µL	16 µL	12.8 µL	3.2 µL	256 µL	up to final volume of 320 µL	320 µL
Off-the-shelf Panel (4 slides)	64 µL	32 µL	25.6 µL	6.4 µL	512 µL	up to final volume of 640 µL	640 µL
Custom Panel with add-on (2 slides)	16 µL***	16 µL	12.8 µL	3.2 µL	256 µL	up to final volume of 320 µL	320 µL
Custom Panel with add-on (4 -slides)	32 µL***	32 µL	25.6 µL	6.4 µL	512 µL	up to final volume of 640 µL	640 µL

\* If using a 100- or 1000-plex assay with a custom Add-on, the custom Add-on replaces the standard Add-on that is part of the assay. (WTX and 6K Discovery Panels do not include a standard Add-on.) If no Add-on is used, use DEPC-water to reach the correct final volume.

\*\* The rRNA marker is the only segmentation marker added in the overnight hybridization. All other segmentation markers are added during the Nuclear and Cell Segmentation step.

\*\*\* Custom stand-alone panels are supplied at 2X concentration compared to off-the-shelf RNA Panels.

8. **Clean all equipment** and benchtop with RNase AWAY and allow to dry; or rinse with DEPC-treated water. The hybridization chamber can be a key source of contamination by oligos. Arrange fresh Kimwipes on bottom of the chamber. Change gloves and clean workspace with RNase AWAY.
9. Wet the Kimwipes with 2X SSC or DEPC-treated water. Take care that the Kimwipes and 2X SSC do not contact the slides. Kimwipes should be thoroughly wet, but without standing buffer.
10. To prevent the tissue from drying, perform the following steps **one slide at a time**.
  - Remove slides from 2X SSC, and gently tap slide to remove excess liquid.
  - Carefully remove the polyester frame backing (with the center square removed) from the incubation frame to expose the top adhesive layer of the incubation frame ([Figure 16](#)). **Ensure that the incubation frame does not lift from the slide when removing the polyester frame backing.**



Figure 16: Remove Polyester Frame backing from Incubation Frame

**FFPE** *In Situ Hybridization (overnight)*

- Lay the slide flat on a clean surface and add 150  $\mu$ L of hybridization solution directly to the tissue within the incubation frame.

Add the hybridization solution at the edge of the tissue opposite of the slide label within the frame. Applying the incubation frame cover in the next step will help move the hybridization solution across the tissue.

Avoid introducing bubbles by leaving a small residual volume in the pipette tip. If a bubble is introduced, carefully aspirate it out using a low-volume pipette tip. When removing air bubbles, removing small amounts of hybridization solution (as long as sufficient solution remains to cover the tissue) is preferable to leaving bubbles.

- Carefully apply incubation frame cover ([Figure 17](#)). Start by setting one edge of the cover down on the incubation frame edge, then gradually lay down the rest of the cover. The tab on the incubation frame cover should face the slide label. As it's lowered, the frame cover should naturally adhere to the incubation frame; no additional pressure around the frame is needed. **Do not press the center of the cover as it could damage the tissue.**




Figure 17: Apply Incubation Frame Cover

- Place the slide horizontally into the hybridization tray ([Figure 18](#)).




Figure 18: Hybridization Tray

- Repeat step 10 for each slide.
  - Ensure that there is a good seal between the incubation frame and the slide, and the incubation frame cover and the frame by checking each slide for any leaks.
11. Close hybridization chamber, insert tray into oven, and clamp tray into place. Incubate at 37°C overnight (16 – 18 hours) ([Figure 19](#)).

**IMPORTANT!** Continue to [Add Enzymes to Buffer 4 on page 41](#) before ending Day 1 slide preparation.




Figure 19: Incubate overnight at 37°C

If your oven does not seal (with a gasket) you may seal your hybridization chamber in a zip-lock bag to simulate a sealed chamber. Chambers sealed in this manner should be tested to ensure they maintain humidity for 24 hours (slides do not dry out) prior to use. Unsealed conditions lead to evaporation of the hybridization solution.


### Add Enzymes to Buffer 4

To reduce oxidative damage to the tissue and probes and improve overall run quality, Bruker recommends adding the active enzymes Catalase and Pyranose Oxidase to instrument Buffer 4 the day before starting a new instrument run.

These instructions cite the different sizes of Instrument Buffer Kits. To review the kit contents, refer to the [CosMx SMI Instrument User Manual](#).

**NOTE:** When supplies are limiting, two Small enzyme aliquots may be shipped in lieu of one Medium enzyme aliquot. Confirm the aliquot size on the tube label before proceeding with these steps.

1. Remove lyophilized Catalase and Pyranose Oxidase vials from 4°C. Centrifuge the vials and resuspend enzymes in DEPC-treated water:
  - For Small enzyme aliquots, add 250 µL of DEPC-water to each tube. Vortex then centrifuge.
  - For Medium enzyme aliquots, add 400 µL of DEPC-water to each tube. Vortex then centrifuge.
2. Add Catalase to Buffer 4:
  - For a Small Buffer 4 (58 mL), add 250 µL resuspended Catalase. For Medium Buffer 4 (90 mL), add 400 µL Catalase. **Avoid precipitate when pipetting from the enzyme tube.**
3. Add Pyranose Oxidase to Buffer 4:
  - Due to the increased risk of precipitate in this enzyme, ensure the resuspended enzyme has been vortexed and centrifuged, and any precipitate has been pelleted at the bottom of the tube.
  - For a Small Buffer 4 (58 mL), add 125 µL resuspended Pyranose Oxidase. For Medium Buffer 4 (90 mL), add 200 µL resuspended Pyranose Oxidase. **Avoid precipitate when pipetting from the enzyme tube.**

 **IMPORTANT:** Precipitate in Buffer 4 could clog the instrument fluidic lines and cause run failure.

4. Mark Buffer 4 bottle after enzymes have been added. Dispose of excess enzymes according to laboratory guidelines.
5. Leave Buffer 4 sealed on the bench until ready to load onto the instrument.

If not already done, prepare 40 mL formamide aliquots following the instructions in the Prepare Reagents section of this protocol and store overnight at 4°C.

## Day 2: Perform Stringent Washes (90 minutes)

You will need the following materials and reagents for this step: **water bath**, **4X SSC**, **100% formamide**, and **2X SSC** (see [Prepare RNA FFPE Assay Reagents on page 22](#)).

**WARNING:** Use of personal protective equipment is advised as formamide is hazardous.

1. Set water bath to 37°C. Remove nuclear stain and cell segmentation kits from the freezer and thaw on ice.
2. **Warm 100% formamide in the 37°C water bath for at least 30 minutes before opening.** Once at temperature, prepare stringent wash directly in staining jars by mixing equal parts 4X SSC and 100% formamide.

Ensure you have sufficient buffer to completely cover the tissue on all slides without submerging the slide labels. Contact with buffer may make slide labels illegible.

3. Preheat staining jars containing freshly-prepared stringent wash in 37°C water bath. It takes ~30 min to preheat.
4. If nearing the 18-hour maximum overnight incubation time, while jars are preheating, transfer slides to 2X SSC.
5. Perform the following steps **one slide at a time** to prevent the tissue from drying.
  - With a clean pair of forceps, carefully remove the incubation frame cover from the incubation frame. Dip slide into 2X SSC as needed to avoid tissue drying. If cover will not come off without removing incubation frame, remove the frame and cover. The frame can be reapplied in a later step.
  - Place slide into 2X SSC and continue to the next slide, cleaning forceps with ethanol between slides.
6. Once both jars have pre-heated to 37°C and all incubation frame covers have been removed, perform the washes shown below ([Figure 20](#)). After the last wash, the slides can be stored in 2X SSC for up to one hour.
  - Gently tap each slide on a Kimwipe to remove excess 2X SSC and place in the first stringent wash for 25 min.
  - Transfer slides to the second stringent wash for 25 min.
  - During second stringent wash, begin preparing reagents for nuclear and cell segmentation staining in the next section.
  - Following stringent washes, immediately transfer slides to 2X SSC and wash for 2 minutes. Transfer slides to a second jar of 2X SSC and wash for 2 minutes. Leave slides in 2X SSC until reagents have been prepared for nuclear and cell segmentation staining.



Figure 20: Perform stringent wash

**IMPORTANT:** Anything coming into contact with hybridization solution (which contains probes), such as containers for stringent wash solution and 2X SSC, needs to be exclusive for this purpose and thoroughly cleaned with RNase AWAY, as probes may contaminate subsequent assays. Use separate staining jars for different probe mixes. Staining jars should be cleaned with RNase AWAY before use.

## Nuclear and Cell Segmentation Staining (2 hours)

You will need the following materials and reagents for this step: **incubation frames, staining jars, 1X PBS, Blocking Buffer (4°C), Nuclear Stain stock (-80°C), and Segmentation Marker Kit (-80°C)** (see [Prepare RNA FFPE Assay Reagents on page 22](#)).

1. Prepare four staining jars of 1X PBS.
2. Prepare Nuclear Stain (220  $\mu$ L per slide):
  - Vortex, then centrifuge thawed Nuclear Stain Stock for at least 1 minute to bring the solution to the bottom of the vial and precipitate insoluble particles.
  - Dilute the Nuclear Stain stock 1:40 as described in [Table 14](#). *Do not pipette from the bottom of the vial.*

Table 14: Prepare Nuclear Stain

	Nuclear Stain Stock	Blocking Buffer	Final Volume
2 slides	11 $\mu$ L	429 $\mu$ L	440 $\mu$ L
4 slides	22 $\mu$ L	858 $\mu$ L	880 $\mu$ L

- Mix by pipetting up and down 3-5 times. Do not vortex.
3. If a new incubation frame is needed, perform the following steps **one slide at a time**.
    - Remove slide from 2X SSC, gently tap slide on a clean Kimwipe. Using a clean Kimwipe, dry the surface of the slide that will come into contact with the incubation frame.

**IMPORTANT:** Avoid wiping the slide within the scan area as this could remove the fiducials needed for imaging. See scan area template on the flow cell assembly tool or follow the guidelines in [Appendix I: CosMx SMI Tissue Sectioning Guidelines on page 84](#).

- Carefully apply a new incubation frame following the previous instructions in this protocol. Ensure that the frame is well adhered to the slide by gently pressing around the frame with clean forceps.
4. Using a clean Kimwipe, carefully wick excess buffer from around the incubation frame as needed. Be careful to not touch the area inside of the incubation frame.
  5. Lay slide horizontally in staining tray and slowly apply up to 200  $\mu$ L of Nuclear Stain directly to tissue. Gently move tray side to side as needed to ensure that the buffer covers the entire tissue.
  6. Repeat with remaining slides and cover tray.



Figure 21: Cover tray and incubate for 15 minutes

7. Incubate slides for 15 minutes at room temperature protected from light (Figure 21).
8. During nuclear stain incubation, prepare Segmentation Mix (**non-neural tissue: blue header (Table 15)**, **neural tissue: magenta header (Table 15)**). Flick each tube to mix and centrifuge before use. **Do not vortex mix.**

**IMPORTANT:** Pipette accurately. When preparing 4 slides, there will be no excess Segmentation Mix.

Table 15: Segmentation Mix for non-neural tissue

	CD298/ B2M Marker Mix	Optional PanCK/ CD45 Marker Mix*	Optional a la carte Marker Mix*	Blocking Buffer	Final Volume
2 slides	8 µL	16 µL	16 µL	360 µL	400 µL
4 slides	16 µL	32 µL	32 µL	720 µL	800 µL

\* If not adding PanCK/CD45 or optional à la carte markers, add Blocking Buffer to reach final volume.

Table 16: Segmentation Mix for neural tissue

	GFAP	Histone	Blocking Buffer	Final Volume
2 slides	16 µL	16 µL	368 µL	400 µL
4 slides	32 µL	32 µL	736 µL	800 µL

The rRNA Segmentation marker is added during the overnight hybridization. All other segmentation markers are added at this step.

9. After nuclear stain incubation, **remove slides one at a time** from staining tray, gently tap slide on a clean Kimwipe to remove excess buffer, and transfer slide to 1X PBS (Figure 22).




Figure 22: Wash in 1X PBS for 5 minutes

10. Wash slide for 5 minutes in 1X PBS.
11. During PBS wash, add 2X SSC or DI water to the staining tray. Do not overfill. The water level should be well below the slides to avoid cross-contamination.
12. Following PBS wash, perform the following steps **one slide at a time** to prevent tissue drying:
  - Remove slide from 1X PBS and gently tap slide on a clean Kimwipe to remove excess PBS.
  - Lay slide horizontally in staining tray and apply up to 200 µL of Segmentation Mix directly to tissue. Gently move tray side to side as needed to ensure that the mix covers the entire tissue.
  - If needed, an incubation frame cover can be placed over the incubation frame to ensure that the mix completely covers the tissue.
  - Adjust volume to add as needed for tissues of varying sizes. The **Segmentation Mix needs to completely cover the tissue** but does not need to completely fill the incubation frame.
13. Repeat with remaining slides and cover tray (Figure 23).

14. Incubate slides for 1 hr at room temperature protected from light. If planning to load the instrument immediately after slides are prepared (recommended), remove the Imaging Tray from 4°C and leave at room temperature until loading the instrument. Confirm that the Imaging Tray has not expired.



Figure 23: Cover tray and incubate for 1 hr.

15. Following Segmentation Mix incubation, transfer slides to 1X PBS and wash for 5 minutes ([Figure 24](#)).
16. Repeat wash 2 times for a total of 3 PBS washes.



Figure 24: Wash 3 times in 1X PBS.

17. If samples will be loaded onto the instrument the same day (recommended), continue to [Flow Cell Assembly](#) [Flow Cell Assembly on page 80](#).

If samples must be stored overnight and loaded onto the instrument the next day, remove the incubation frame following the previous instructions in this protocol, transfer slides to fresh 2X SSC, and store according to **Safe Storage Guidelines for RNA Slides** [Safe Storage Guidelines for RNA Slides on page 79](#).

Once the flow cell has been applied to the slide, proceed to run the RNA assay on the CosMx SMI instrument.

Refer to the [CosMx SMI Instrument User Manual \(MAN-10161\)](#) for instructions to operate the instrument.

## RNA Fresh Frozen Manual Slide Preparation

### Slide Preparation Workflow





Figure 25: CosMx SMI Workflow Overview

**Day 1: Slide Preparation.** Prepare slides and incubate biological targets with RNA specific probes.

**Day 2: Process Slides on CosMx SMI.** Remove off-target probes and add cell segmentation markers to each slide. Load prepared flow cells into the CosMx SMI instrument and enter flow cell/study information. Tissue is scanned to capture RNA readout and morphology imaging within user-designated fields of view (FOVs).


**After run completion: Create a Data Analysis study** in the AtoMx Spatial Informatics Platform (SIP) and perform quality-control checks, data analysis, and generate analysis plots.

## Day 1: Slide Preparation



### NBF Fixation and Bake

- 10% NBF Fixation
- 3 washes in 1X PBS
- 30-minute bake




### Wash and Rehydrate Tissue

- 3 washes in 1X PBS
- Wash in 4% SDS
- 3 washes in 1X PBS
- Wash in 50% EtOH
- Wash in 70% EtOH
- 3 washes in 100% EtOH




### Target Retrieval

- 15 mins at 100°C
- H2O rinse and EtOH wash
- Dry for 30 min-1hour




### Protease Digestion

- Apply incubation frame
- Apply digestion buffer
- Incubate at RT for 30 mins
- 2 washes in 1X PBS



### Apply Fiducials

- Prepare and apply fiducials
- Incubate for 5 minutes
- Wash with 1X PBS (not shown)




### Post-Fix Tissue

- Wash in 10% NBF for 1 min.
- 2 washes, 5 mins each, of NBF Stop Buffer
- Wash 5 min in 1X PBS




### Blocking

- Prepare and apply NHS-Acetate
- Incubate for 15 mins.
- 2 washes, 5 mins each, in 2X SSC



### Overnight Hybridization

- Prepare and apply assay specific probes
- Incubate at 37°C overnight




### Add Buffer 4 Enzymes

- Add P2OX and Catalase to Buffer 4 and leave sealed on bench until instrument loading.



## Day 2: Wash and Stain Slide



### Stringent Washes

- 2 stringent washes, 25 mins each.
- 2 washes, 2 mins each, 2X SSC



### Blocking & Nuclear Stain


- Prepare Nuclear Stain stock and apply to tissue
- Incubate 15 mins at RT
- Wash in 1X PBS for 5 mins



### Segmentation Markers


- Prepare Segmentation mix and apply to tissue
- Incubate 1 hour at RT
- 3 washes, 5 mins each, 1X PBS

## Prepare Flow Cells and Load Instrument



### Prepare Flow Cells


- Use the flow cell assembly tool to assemble the flow cells.



### Load Instrument

- Follow on-instrument prompts to load reagents and flow cells.
- Begin instrument run.

## Analyze Data in AtoMx SIP



## Equipment, Materials, and Reagents

The following equipment ([Table 17](#)), materials ([Table 18](#)), and reagents ([Table 19](#)) are required for this protocol but are **not supplied by Bruker**.




Table 17: Equipment not supplied by Bruker

	Equipment	Source	Part Number
<input type="checkbox"/>	<b>Baking Oven</b>	Quincy Lab, Inc (or comparable)	<a href="#">Various GC Models</a>
<input type="checkbox"/>	<b>Hybridization Oven including hybridization chamber</b> <ul style="list-style-type: none"> <li>RapidFISH Slide Hybridizer <b>or</b></li> <li>HybEZ oven with humidity control tray</li> </ul> <p><b>NOTE:</b> These hybridization ovens are designed to keep the slides hydrated and maintain a precise temperature overnight. Bruker does not recommend the use of any other hybridization ovens for CosMx SMI Sample Prep</p>	Boekel Scientific or ACDBio	<a href="#">240200</a> for 120V or <a href="#">240200-2</a> for 230V <a href="#">321710/321720</a> <a href="#">310012</a>
<input type="checkbox"/>	<b>6-quart pressure cooker</b> <p><b>NOTE:</b> Products from other vendors may require testing and optimization.</p> <p><b>OPTIONAL:</b> A steamer may be used in lieu of a pressure cooker and may be preferred for more fragile tissues. If a steamer is used, a thermometer is also needed.</p>	BioSB® TintoRetriever  Nesco® Hamilton Beach® Ovente® (not validated)	<a href="#">BSB 7015</a>  <a href="#">ST-25F</a> <a href="#">37530MN</a> <a href="#">FS62S</a>
<input type="checkbox"/>	<b>Ultrasonic Bath</b> 500 mL capacity, 40 kHz frequency with timer	General Lab Supplier	<a href="#">Example</a>
<input type="checkbox"/>	<b>Vortex mixer</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>Microcentrifuge</b> for 1.5 mL microcentrifuge tubes and 8-well PCR strip tubes	General Lab Supplier	Various
<input type="checkbox"/>	<b>Water Bath</b> (temperature setting of 37°C)	General Lab Supplier	Various
<input type="checkbox"/>	<b>Thermal Cycler</b> <p><b>NOTE:</b> Must include 96-well 200 µL tube block</p>	General Lab Supplier	Various
<input type="checkbox"/>	<b>Analytic Scale with draft shield</b> <p><b>NOTE:</b> Ensure scale can weigh in milligrams (mg) for accuracy.</p>	Various	<a href="#">Example</a>

Table 18: Materials not provided by Bruker

	Material	Source	Part Number
<input type="checkbox"/>	<b>Pipettes for 2.0 – 1,000 µL</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>Filter tips (RNase/ DNase Free)</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>2.0 mL Centrifuge Tubes (RNase/ DNase Free)</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>0.2 mL PCR tubes</b> or PCR strip tubes	General Lab Supplier	Various
<input type="checkbox"/>	<b>Leica BOND Plus slides or VWR Superfrost Plus Micro Slide, Premium</b> <b>NOTE:</b> These slides have been internally validated. Do not use other products. <b>NOTE:</b> Leica BOND Plus slides are preferred for tissue sections prone to peeling.	Leica Biosystems VWR	<a href="#">S21.2113.A</a> <a href="#">48311-703</a>
<input type="checkbox"/>	<b>Slide Rack</b>	General Lab	<a href="#">Example</a>
<input type="checkbox"/>	<b>Polypropylene Slide Staining Jars</b> (24 required) or <b>Slide Staining Station</b> <b>NOTE:</b> Due to the photo-sensitivity of this assay, the staining jars should be impermeable to light.	Ted Pella (or comparable) Amazon Fisher Scientific	<a href="#">21029</a> <a href="#">MH-SJ6302</a> <a href="#">NC1862866</a>
<input type="checkbox"/>	<b>Forceps</b> (for slide handling)	General Lab Supplier	Various
<input type="checkbox"/>	<b>Razor Blades</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>Timer</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>RNase AWAY</b> <b>NOTE:</b> RNase ZAP and other alternatives cannot be used as substitutes as they do not adequately remove both nucleic acid and nuclease contaminants.	ThermoFisher	<a href="#">7003PK</a>
<input type="checkbox"/>	<b>Kimwipes (large and small)</b>	General Lab Supplier	Various
<input type="checkbox"/>	<b>StainTray slide staining system with black lid</b>	Sigma Aldrich	<a href="#">Example</a>

Table 19: Reagents not provided by Bruker

	Reagent	Source / Part Number	Storage Conditions
<input type="checkbox"/>	DEPC-Treated Water	ThermoFisher, <a href="#">AM9922</a> (or comparable)	Room temperature
<input type="checkbox"/>	100% Ethanol (EtOH): ACS grade or Better 	General Lab Supplier	Flammable Storage (RT)
<input type="checkbox"/>	10X Phosphate Buffered Saline pH 7.4 (PBS)	ThermoFisher, <a href="#">AM9625</a> (or comparable)	Room temperature
<input type="checkbox"/>	SDS, 10% Solution, RNase-free	ThermoFisher, <a href="#">AM9822</a>	Room temperature
<input type="checkbox"/>	20X SSC (DNase, RNase free)	ThermoFisher, <a href="#">AM9763</a>	Room temperature
<input type="checkbox"/>	Tris Base	Sigma-Aldrich, <a href="#">10708976001</a> (or comparable)	Room temperature
<input type="checkbox"/>	Glycine	Sigma-Aldrich, <a href="#">G7126</a> (or comparable)	Room temperature
<input type="checkbox"/>	Sulfo NHS-Acetate powder <b>NOTE:</b> Sulfo-NHS-Acetate powder is shipped in a plastic bag with a desiccant and should be left in the bag and stored at -20°C until ready to use.	Fisher Scientific, <a href="#">26777</a>	-20°C
<input type="checkbox"/>	10% Neutral Buffered Formalin (NBF) 	EMS Diasum, <a href="#">15740</a> (or comparable)	Room temperature
<input type="checkbox"/>	100% Deionized Formamide  <b>NOTE:</b> Deionized formamide is optimal, however, formamide that is not deionized may also be used.	ThermoFisher, <a href="#">AM9342</a> or VWR, <a href="#">VWRV0606</a> (or comparable)	4°C (bring to RT for at least 30 before opening)

**OPTIONAL:** Large volume stock solutions (>500 mL) of deionized formamide can be aliquoted and stored, protected from light, at 4°C (recommend 40 mL aliquots in 50 mL conical tubes). This will save time during day 2 slide preparation.

**Bruker Supplied Reagents**

**CosMx Fresh Frozen Slide Preparation Kit (RNA)**

Table 20: Fresh Frozen Slide Preparation Kit (Box 1 of 2)








Kit Contents (Box 1 of 2, store at 4°C)	
	
10X Target Retrieval Solution	Buffer R
	
CosMx RNA Blocking Buffer	NHS-Acetate Buffer
	
2X SSC-T	CosMx Fiducials
Incubation Frames and covers (not pictured)	Protease A Buffer (not pictured)

Table 21: Fresh Frozen Slide Preparation Kit (Box 2 of 2)

Kit Contents (Box 2 of 2, store at -20°C)	
CosMx Proteinase K	

**CosMx RNase Inhibitor**

Table 22: CosMx RNase Inhibitor

Kit Contents (Store at -20°C)	
RNase Inhibitor is <b>sold separately</b> and is used in RNA hybridization and instrument loading.	

**Fresh Frozen** *Bruker Supplied Reagents*

**CosMx RNA Panels** (see [Panel and Cell Segmentation Marker Selection on page 10](#))



Table 23: CosMx RNA Panels

CosMx RNA Panels (store at -20°C)	
Kit Name	Kit Component
Human Whole Transcriptome Panel (WTX), 19K-plex, RNA	CosMx Hs WTX RNA Probe Mix
Human 6K Discovery Panel, 6K-plex, RNA	CosMx Hs 6K Discovery RNA Probe Mix
Human Universal Cell Characterization (UCC) Panel, 1000-plex, RNA	CosMx Hs UCC RNA Probe Mix CosMx Hs UCC RNA Add-On Custom RNA Add-On replaces off-the-shelf RNA Add-On
Human Immuno-oncology Panel, 100-plex, RNA	CosMx Hs IO RNA Probe Mix CosMx Hs IO RNA Add-On Custom RNA Add-On replaces off-the-shelf RNA Add-On
Mouse Neuroscience Panel, 1000-plex, RNA	CosMx Mm Neuro RNA Probe Mix CosMx Mm Neuro RNA Add-On Custom RNA Add-On replaces off-the-shelf RNA Add-On
Mouse Universal Cell Characterization (UCC) Panel, 1000-plex, RNA	CosMx Mm UCC RNA Probe Mix CosMx Mm UCC RNA Add-On Custom RNA Add-On replaces off-the-shelf RNA Add-On


**CosMx Segmentation Markers**

Table 24: CosMx SMI Cell Segmentation and Supplemental Marker Kits

CosMx RNA Segmentation Marker Kits (each kit sufficient for up to 4 slides, store at -80°C)	
Kit Description	Kit Components
Human Universal Cell Segmentation Kit (RNA) Compatible with: Human WTX Panel, Human 6K Discovery Panel, Human Immuno-Oncology 100-plex Panel, and any 1000-plex RNA Panel	CosMx DAPI Nuclear Stain, Ch1
	CosMx Hs CD298/B2M Segmentation Marker Mix, Ch2 (RNA)
	CosMx Hs PanCK/CD45 Marker Mix Ch3/Ch4 (RNA)
Human Neuroscience Cell Segmentation Kit (RNA) Compatible with: Human WTX Panel, Human 6K Discovery Panel	CosMx DAPI Nuclear Stain, Ch1
	CosMx Hs Neuro rRNA Neuro Marker, Ch2 (RNA)
	CosMx Mm/Hs Neuro Histone Marker, Ch3 (RNA)
	CosMx Mm/Hs GFAP Marker, Ch4 (RNA)
Mouse Neuroscience Cell Segmentation Kit (RNA) Compatible with: Mouse Neuroscience 1000-plex Panel	CosMx DAPI Nuclear Stain, Ch1
	CosMx Mm Neuro rRNA Neuro Marker, Ch2 (RNA)
	CosMx Mm/Hs Neuro Histone Marker, Ch3 (RNA)
	CosMx Mm/Hs GFAP Marker, Ch4 (RNA)
Mouse Universal Cell Segmentation Kit (RNA) Compatible with: Mouse Universal Cell Characterization 1000-plex Panel	CosMx DAPI Nuclear Stain, Ch1
	CosMx Mm CD298/B2M Marker Mix, Ch2 (RNA)
	CosMx Mm PanCK/CD45 Marker Mix Ch3/Ch4 (RNA)

**Fresh Frozen** Bruker Supplied Reagents

The following markers are **optional and available to order à la carte**:

Table 25: A la carte Markers

Compatible with	Item Description (store at -80°C)
Human Universal Cell Segmentation Kit (RNA)	CosMx Hs CD68 A La Carte Marker, Ch5 (RNA)
	CosMx Hs Cytokeratin 8/18 A La Carte Marker, Ch5 (RNA)
	CosMx Hs/Mm CD3 A La Carte Marker, Ch5 (RNA)
Mouse Universal Cell Segmentation Kit (RNA)	CosMx Mm CD68 A La Carte Marker, Ch5 (RNA)
	CosMx Mm CD8 A La Carte Marker, Ch5 (RNA)
	CosMx Hs/Mm CD3 A La Carte Marker, Ch5 (RNA)

**Flow Cell Assembly Tool and Kit**

The Flow Cell Assembly tool is a one-time required purchase.

The Flow Cell Assembly Kit contains 4 single-use Flow Cell coverslips, sufficient for a 4-slide experiment.

## Tissue Sectioning and Slide Preparation: Fresh Frozen Samples

[Selecting Fresh Frozen Blocks on page 87](#) covers Fresh Frozen block selection and sectioning in detail. Review these guidelines prior to beginning Fresh Frozen Sample Preparation.

Fresh frozen blocks should be sectioned at **5  $\mu\text{m}$  thickness** and mounted on the label side of Leica BOND PLUS slides or VWR Superfrost Plus Micro Slides. Blocks may be sectioned up to 10  $\mu\text{m}$  thickness; however, the instrument will only image the 5  $\mu\text{m}$  closest to the slide.

Tissue sections must be centered within the Scan Area of the slide (the green rectangle in [Figure 26](#)) and be no larger than **20 mm long by 15 mm wide** (image not to scale; see the Flow Cell Assembly Tool for a to-scale template). For best performance, there should also be some bare glass (not covered by tissue) in the Scan Area. For examples of tissue placement best practices, see [Appendix I: CosMx SMI Tissue Sectioning Guidelines on page 84](#).

Label slides with pencil on the frosted label according to lab guidelines. If using an adhesive slide label, ensure the label is less than 295  $\mu\text{m}$  thick and is not folded over on itself. Labels over the maximum thickness or labels that are not properly adhered may result in slide or flow cell damage during flow cell assembly and/or instrument loading.

**Dry mount slides for 5-10 minutes at room temperature**, or until dry. Once dry, store slides at  $-80^{\circ}\text{C}$ .




Figure 26: Tissue Scan Area (not to scale)

**WARNING:** If the tissue sections are less than 5  $\mu\text{m}$  or the tissue is poor quality, the sample may degrade during the SDS wash step.

If sections are larger than the allowable size, or placed off-center, excess tissue must be scraped away with a clean razor blade. Tissue may be scraped when slides are dry at the start of slide preparation, before applying the incubation frame during slide preparation, or just before applying the flow cell coverslip after slide preparation. The optimal practice (minimizing tissue pulling, folding, or detaching) may depend on your sample type. Applying the incubation frame over tissue could result in tissue damage when the frame is removed or poor sealing of the incubation frame.

**IMPORTANT:** The CosMx SMI instrument will only image the area inside the flow cell chamber, the tissue scan area. If the tissue section is outside of the scan area, it will not be imaged.

**Fresh Frozen Prepare RNA Fresh Frozen Assay Reagents****Prepare RNA Fresh Frozen Assay Reagents**

**i IMPORTANT:** Take care to maintain nuclease-free conditions. The greatest risk of contamination comes from CosMx SMI RNA probes and other oligos. Bruker recommends the use of RNase AWAY (ThermoFisher [7003PK](#)), as it will limit contamination from oligos, detection probes, and nucleases. After using RNase AWAY, allow to air dry completely, or rinse with DEPC-treated water. See manufacturer's instructions for details.


**Label staining jars and prepare reagents** using the instructions in the following table ([Table 26](#)). Ensure you have sufficient buffer to completely cover the tissue on all slides without submerging the slide labels. Contact with buffer may make slide labels illegible.

Unless otherwise noted, reagents can be made up to 2 weeks in advance and stored at room temperature.

Table 26: RNA Fresh Frozen Reagent Preparation

	Reagent	Dilution	Storage
<input type="checkbox"/>	1X PBS (pH 7.4)	Prepare 1 L of 1X PBS by combining 100 mL of 10X PBS and 900 mL of DEPC-treated water.	Room temperature
<input type="checkbox"/>	2X SSC	Prepare 1 L of 2X SSC by combining 100 mL of 20X SSC and 900 mL of DEPC-treated water.	Room temperature
<input type="checkbox"/>	4X SSC	Prepare 1 L of 4X SSC by combining 200 mL of 20X SSC and 800 mL of DEPC-treated water.	Room temperature
<input type="checkbox"/>	NBF stop buffer (Tris Glycine Buffer)	Combine 6.06 g Tris base and 3.75 g Glycine in 500 mL of DEPC-treated water. The final concentration of Tris and Glycine will be 0.1 M each.	Room temperature
<input type="checkbox"/>	70% Ethanol (EtOH)	Prepare 50 mL of 70% ethanol by combining 15 mL of DEPC-treated water and 35 mL 100% ethanol.	Make fresh each time
<input type="checkbox"/>	50% Ethanol	Prepare 50 mL of 50% ethanol by combining 25 mL of DEPC-treated water and 25 mL 100% ethanol.	Make fresh each time
<input type="checkbox"/>	4% SDS in 1X PBS	Before diluting, warm 10% SDS for 10 minutes in a 37°C water bath. Vortex for 1 minute.  Prepare 50 mL of 4% SDS in 1X PBS by combining 5 mL of 10X PBS and 20 mL of 10% SDS with 25 mL of DEPC-treated water.  <b>NOTE:</b> The concentration of SDS is critical. Carefully measure each volume to ensure the concentration does not exceed 4%.	Make fresh each time
<input type="checkbox"/>	1X Target Retrieval Solution	Prepare 50 mL <b>fresh each time</b> by adding 5 mL of CosMx Target Retrieval Solution, 10X (provided by Bruker) to 45 mL DEPC-treated water.	Make fresh each time
<input type="checkbox"/>	NHS-Acetate Mix	Individual aliquots can be prepared prior to slide preparation but should be stored at -20°C with a desiccant until use.  To prepare aliquots for 4 slides:	-20°C with desiccant

	Reagent	Dilution	Storage
		<ul style="list-style-type: none"> <li>• Bring stock to room temperature for 1-2 hours prior to opening.</li> <li>• Prelabel four 2.0 mL screw top centrifuge tubes.</li> <li>• Using a weighing spatula, carefully weigh 25 mg of Sulfo-NHS-Acetate directly into one screw-top tube on an analytic scale.</li> <li>• Close the tube and label with final weight. Seal the tube with parafilm and place into the -20°C with desiccant until use.</li> </ul> <p><b>NOTE:</b> If preparing only 2 slides, pre-weigh 15 mg aliquots into 6 tubes.</p>	
<input type="checkbox"/>	Digestion Buffer	Prepare immediately before use. See <a href="#">Tissue Permeabilization (40 minutes)</a> on page 63.	Make immediately before use
<input type="checkbox"/>	Fiducials and Hybridization Mix	On Day 1, remove RNA Probe Mix from -20°C and keep on ice until ready to use. <b>NOTE:</b> Fiducials are light sensitive and should be kept stored, protected from light, until instructed to remove later in this protocol.	n/a
<input type="checkbox"/>	Day 2 Reagents	These reagents have additional steps that will be covered in detail in their respective sections. Follow Day 2 procedure for preparation of these reagents.	n/a

 **IMPORTANT:** If using a custom panel, email [AtoMx.KitAdmin@Bruker.com](mailto:AtoMx.KitAdmin@Bruker.com) at least one business day prior to a run so that the custom kit is available in the CosMx SMI Control Center.

**Fresh Frozen** Day 1: NBF Fixation and Bake**Day 1: NBF Fixation and Bake**

You will need the following materials and reagents for this step: **staining jars, 10% NBF** and **baking oven**.

- Pre-cool a staining jar filled with 10% NBF to 4°C for at least 15 minutes. Ensure you have sufficient buffer to completely cover the tissue on all slides without submerging the slide labels. Contact with buffer may make slide labels illegible.
- Preheat oven to 60°C.
- Remove sample slides from -80°C freezer and place onto dry ice to carry sample slides to slide preparation area.

**NBF Fixation:**

1. Remove slides from dry ice and immediately transfer to pre-cooled 10% NBF ([Figure 27](#)). Incubate slides in 10% NBF for 2 hours at 4°C.



10% NBF at 4°C  
1 wash  
(2 hours)

Figure 27: 2 Hour NBF Fixation

2. Following NBF fixation, transfer slides to 1X PBS and wash for 2 minutes ([Figure 28](#)).
3. Repeat PBS wash twice, using new staining jars for each wash, for a total of 3 washes.



1X PBS  
3 washes  
(2 min each)

Figure 28: Wash 3x with 1X PBS

4. After final PBS wash, place slides into slide holder and bake vertically at 60°C for 30 minutes.

## Prepare Materials and Tissue for Target Retrieval (1 hour)

You will need the following materials and reagents for this step: **pressure cooker or steamer, staining jars, 1X PBS, 4% SDS, 100% ethanol, 70% ethanol, 50% ethanol, and DEPC-treated water.** See [Prepare RNA Fresh Frozen Assay Reagents on page 56](#) for more details.

### Preheat Target Retrieval Solution

1. Prepare the pressure cooker or steamer and the Target Retrieval Solution. Content in **purple boxes** denotes instructions specific to the pressure cooker. Content in **orange boxes** denotes instructions specific to the steamer. Use the same target retrieval method (pressure cooker or steamer) throughout the study.

**WARNING:** Bruker does not recommend the use of glass staining jars in the pressure cooker.


#### Pressure Cooker Method:

1. Fill the pressure cooker with water to the correct level per the manufacturer's instructions (4-8 cups depending on model used).
2. Place the staining jar containing freshly prepared 1X Target Retrieval Solution into the pressure cooker to preheat. Ensure that water level is well below lid of jar; about halfway up jar is sufficient. If needed, a trivet may be used to raise the jar.
3. Place a lid on the staining jar to prevent evaporation. To prevent pressure from building within the staining jar, **do not fully seal the lid to the jar.**
4. Preheat the pressure cooker to 100°C following the model-specific instructions below. Pressure cooker preheating takes about 1 hour.

#### BioSB Preheating Instructions:

Use the *TintoRetriever Pressure Cooker Preheating Cycle QuickStart Guide, Two Staining Dish Operation* provided with the BioSB pressure cooker to preheat the pressure cooker.

- With the pressure valve closed, press the **80°C** button on the face of the pressure cooker and press **Start** to run a cycle at 80°C with a 0-minute timer.
- Once the first cycle is complete, run a second cycle with a **45-minute timer at 100°C.**



#### Steamer Method

1. Fill the steamer reservoir up to the fill line with water.
2. Place two staining jars inside of the steamer, one containing DEPC-treated water and one containing 1X Target Retrieval Solution. Ensure sufficient reagent volume to cover slides up to the label.
3. Loosely cover each jar with aluminum foil instead of the jar lid to allow for a thermometer reading in a later step.

**Fresh Frozen** Prepare Materials and Tissue for Target Retrieval (1 hour)

4. Preheat the steamer to 100°C. More water may need to be added to the steamer during preheating.

The steamer may take up to 1 hour to heat the liquid in the jars to a stable maximum temperature near 100°C. Final temperature can be checked by inserting a digital thermometer through the hole in the lid of the steamer into the DEPC-water staining jar.

**Wash and Dehydrate Fresh Frozen Tissue Sections**

**Time Critical Step!** The following steps are time sensitive. Be sure to use a timer and transfer slides between washes carefully but quickly to avoid additional time in each wash. Exceeding the wash time may result in tissue degradation and loss of data.

1. After NBF fixation, gently perform the following washes in staining jars ([Figure 29](#)). **Wash Steps 2 and 3 (4% SDS and 1X PBS) are only for neural tissue. For other tissue types, proceed from Wash Step 1 (1X PBS) to Wash Step 4 (50% ethanol).**

Ensure you have sufficient buffer to completely cover the tissue on all slides without submerging the slide labels. Contact with buffer may make slide labels illegible.




Figure 29: Rehydrate and Fix Tissue

2. During the final ethanol wash, clean workspace with RNase AWAY and lay down clean Kimwipes. Once final wash is complete, remove slides from ethanol, lay horizontally on Kimwipes, and air dry for at least 10 minutes at room temperature.
3. **If not already aliquoted, remove stock NHS-Acetate powder from -20°C freezer** and leave at room temperature for 1-2 hours before use in [NHS-Acetate Preparation and Application \(25 minutes\) on page 69](#). Pre-aliquoted NHS-Acetate powder can remain at -20°C until instructed to remove [on page 69](#).

## Perform Target Retrieval (50 minutes)

You will need the following materials and reagents for this step: **pressure cooker or steamer, staining jars, 1X Target Retrieval Solution (pre-heated in the previous step), DEPC-treated water and 100% ethanol.**

Target retrieval times were determined based on FFPE tissue blocks meeting the constraints outlined in the sample guidance section. These conditions may vary by sample, the amount of normal adjacent tissue, and other factors. See [Appendix II: Tissue-Specific Slide Preparation Considerations on page 88](#).

Content in **purple box** denotes instructions specific to the pressure cooker. The **orange box** denotes instructions specific to the steamer. Use the same target retrieval method (pressure cooker or steamer) throughout the study.


### Pressure Cooker Method

- Once Target Retrieval Solution is preheated, press Cancel on the pressure cooker, release the pressure valve, and wait for the pressure cooker to release pressure. Once pressure has released, remove the lid and **carefully but quickly remove the staining jar containing preheated 1X Target Retrieval Solution.** Once removed, the Target Retrieval Solution will begin to rapidly cool. Ensure the following steps are done as quickly and safely as possible.




**⚠ WARNING:** When opening the pressure valve and removing the pressure cooker lid, hot steam will be released. The staining jar will also be VERY hot. Wear protective heat-resistant gloves to open lid and remove staining jar.

- Place FFPE slides into the preheated solution and replace lid on the staining jar to prevent evaporation. To prevent pressure from building within the container, **do not fully seal the lid to the jar.**
- Return the staining jar containing the slides into the preheated pressure cooker. Ensure that water level is well below lid of jar. About halfway up jar is sufficient. If needed, a trivet may be used to raise the jar.
- Reattach the pressure cooker lid, **open the pressure release valve to *Pressure Release* position** and return the pressure cooker to 100°C. For the BioSB, this can take up to 20



### Steamer Method

- Without removing the lid, place an instant-read digital thermometer through the vents in the steamer lid and pierce the aluminum foil covering the DEPC-treated water. Ensure the water has reached about 99°C. Add more water as needed.
- Once the water has reached 99°C, carefully remove the steamer lid. Once removed, the Target Retrieval Solution will begin to rapidly cool. Ensure the following steps are done as quickly and safely as possible.



**⚠ WARNING:** Removing the steamer lid releases hot steam. Wear protective heat-resistant gloves to open lid and remove the staining jar. Transfer slides using forceps or rack.

- Remove the foil from the Target Retrieval Solution jar and quickly transfer the slides to the Target Retrieval Solution. Replace the foil, then replace steamer lid.
- Reinsert the thermometer into the DEPC-water jar and wait until the temperature returns to about 99°C.
- Once the steamer temperature returns to 99°C, start timer and run for 15 minutes for FFPE tissue or 8 minutes for cell pellet arrays**

minutes.

**NOTE:** If after 10 minutes the temperature has not started to increase it may indicate that the float valve is not properly seated. Without removing the lid, turn the lid to the unlock position and then slowly back to lock.

5. **Once the pressure cooker temperature returns to 100°C**, start timer and run for **15 minutes** for FFPE tissue or 8 minutes for cell pellet arrays (CPA).
6. When the timer reaches zero, press Cancel on pressure cooker to stop heating, carefully remove the pressure cooker lid, and remove the staining jar.

(CPA).

6. When the timer reaches zero, carefully remove the steamer lid and remove the staining jar.

7. **Immediately transfer all slides to staining jar containing fresh DEPC-treated water.** Move slides up and down for 15 seconds to wash ([Figure 30](#)).

Transfer slides from Target Retrieval Solution into DEPC-treated water carefully but quickly to ensure slides do not dry out. Drying can occur quickly when slides are hot.

8. Transfer slides to fresh 100% ethanol and incubate for 3 minutes.
9. During ethanol wash, clean bench space with RNase AWAY and lay out a fresh Kimwipe.
10. After 3 minutes, remove slides from ethanol and lay horizontally on a clean Kimwipe. Dry at room temperature for 30 minutes to 1 hour.




Figure 30: Water and Ethanol Wash

11. While slides are drying, continue to next page to prepare incubation frames and digestion buffer.

## Tissue Permeabilization (40 minutes)

You will need the following materials and reagents for this step: **staining jars, hybridization oven, hybridization tray, incubation frames** (FF Slide Prep Kit Box 1), **Proteinase K, Protease A, Protease A Buffer, 1X PBS** and **DEPC-treated water**.

1. If needed, trim the tissue following the template in [Tissue Sectioning and Slide Preparation: Fresh Frozen Samples on page 55](#).
  - Use a clean razor blade to trim tissue and change blade as needed to ensure clean cuts and reduce the risk of cross-contamination between samples.
2. Prepare the Incubation Frame:
  - Separate an individual frame from the strip by carefully tearing along the perforations.
  - Each frame is sandwiched between a thin solid polyester sheet and a thick polyester frame backing (with the center square removed).
3. Using a clean Kimwipe, ensure that the surface of the slide that will come in contact with the incubation frame is clean and dry.
4. Apply the incubation frame ([Figure 31](#)).
  - Carefully remove the **thin polyester sheet**, ensuring that the frame remains bound to the thick polyester frame backing (with the center square removed).
  - With the slide on a flat surface, careful not to touch the adhesive, **center the tissue within the incubation frame and carefully place the incubation frame** around each tissue section. Lightly press along the border of the incubation frame to ensure that it is well adhered to the slide.




Figure 31: Apply Incubation Frame

5. With the slide still on a flat surface, use a clean razor blade to carefully trim the long edges of the incubation frame to remove excess plastic ensuring that there is no excess film extending over slide edges. If needed, trim the short end of the frame opposite of the slide label.
6. Prepare **Digestion Buffer**:

For **tissue types other than fresh frozen mouse neural tissue**, see the **blue box**. For **fresh frozen mouse neural tissue**, follow the guidance in the **magenta box**.

**Fresh Frozen Tissue Permeabilization (40 minutes)****Digestion Buffer for Fresh Frozen tissue types other than mouse neural tissue:**

Dilute the 20 mg/mL Protease Solution (Proteinase K stock; provided by Bruker) to a working concentration of 3 µg/mL\* in 1X PBS. **Prepare fresh each time and store on ice until ready to use.**

A 2-step serial dilution is recommended. Accurate dilution of Proteinase K is critical for assay performance.

- Step 1: Dilute 20 mg/mL stock to 200 µg/mL by adding 2 µL of Proteinase K stock to 198 µL of 1X PBS.
- Step 2: Dilute the 200 µg/mL solution made in step 1 to the target concentration of 3 µg/mL by adding 30 µL of the 200 µg/mL solution to 1970 µL of 1X PBS.
- Mix thoroughly by inverting tube or pipetting up and down using a clean pipette tip. **Do not vortex.**

\* **NOTE:** The optimal concentration may differ for some tissue types, including CPA samples (see [Appendix II: Tissue-Specific Slide Preparation Considerations on page 88](#)).

**Digestion Buffer for Fresh Frozen mouse neural tissue:**

Digestion buffer for Fresh Frozen mouse neural tissue should be **made within 10 minutes of use**. Accurate dilution is critical for assay performance. 2-step serial dilutions are recommended. Follow the steps below (also shown in [Table 27](#)).

1. **Resuspend Protease A** in 200 µL of Protease A Buffer to yield Protease A stock solution. Mix thoroughly; do not vortex.
2. **Create Protease A working solution:** Combine 5 µL of Protease A stock solution with 245 µL of Protease A Buffer to yield Protease A working solution. Mix thoroughly; do not vortex.
3. **Create Proteinase K working solution:** Combine 2 µL of Proteinase K stock solution with 398 µL of 1X PBS to yield Proteinase K working solution. Mix thoroughly; do not vortex.
4. **Create Digestion Buffer:** Combine 16.5 µL of Protease A working solution with 82.5 µL of Proteinase K working solution and 1551 µL of Protease A Buffer to yield Digestion Buffer (sufficient for 4 slides; adjust volume as needed). Mix thoroughly; do not vortex.

Table 27: Digestion Buffer Preparation

To Make Working Solutions		
Stock	5 µL resuspended Protease A	2 µL Proteinase K
Diluent	245 µL of Protease A Buffer	398 µL of 1X PBS
Total Volume Working Solution	250 µL ProA working solution	400 µL ProK working solution
To Make Digestion Buffer		
Working Solution	16.5 µL ProA working solution	82.5 µL ProK working solution
Diluent	1551 µL Protease A Buffer	
Total Volume Digestion Buffer	1650 µL Digestion Buffer	

7. With the slide on a clean, flat surface, apply 400  $\mu\text{L}$  (2 x 200  $\mu\text{L}$ ) of digestion buffer to completely cover tissue within incubation frame. Gently move the slide side to side as needed to ensure that digestion buffer covers the entire tissue.
8. Incubate for 30 minutes at room temperature.

**IMPORTANT:** After digestion buffer has been applied, avoid tissue drying in subsequent steps by working with only one slide at a time.

9. During slide incubation, remove fiducials and 2X SSC-T from 4°C and let come to room temperature, protected from light, for at least 10 minutes. Then, the Fiducial Preparation steps can be started in parallel with the below wash steps.
10. After slide incubation, tap off excess digestion buffer *one slide at a time* and transfer slides to staining jar containing fresh 1X PBS. Wash for 5 minutes (Figure 32).
11. Transfer slides to a new jar of 1X PBS and wash for 5 minutes.
12. Slides can be stored in 1X PBS while fiducials are prepared.




Figure 32: Wash 2X in 1X PBS

Limit the time that the slides are kept in 1X PBS to minimize risk of target loss. At this point, targets have been exposed and tissue has not undergone post-fixation.

**Fresh Frozen Fiducial Preparation and Application (20 minutes)****Fiducial Preparation and Application (20 minutes)**

You will need the following materials and reagents for this step: **ultrasonic bath**, **vortex**, **2X SSC-T** (at RT, provided by Bruker), **fiducials** (at RT), **staining jars**, **staining tray** (clean and dry), and **1X PBS**.

The volume of fiducial working solution prepared here is 1 mL, sufficient for 4 slides. Do not make less working solution due to high risk of clumping. Additional volume may be made as needed.

1. Refer to the graphic below to vortex and sonicate the fiducials prior to use ([Figure 33](#)):

**IMPORTANT:** When sonicating the fiducial tube, be sure not to submerge the tube cap under the liquid level in the sonicator. Use a floating tube holder if needed to float the tubes in the ultrasonic bath.




Figure 33: Fiducial preparation

Failure to follow these steps will result in fiducial clumping and uneven distribution of fiducials within the tissue. This uneven distribution can result in a loss of readable area or loss of image registration.

2. Once fiducials are prepared following [Figure 33](#), **dilute fiducial stock (0.1%) to the working concentration (0.00015%) in 2X SSC-T**. A 2-step serial dilution is recommended for fiducial preparation. Accurate dilution of fiducials is critical for assay performance.
  - Dilute stock to 0.01% by adding 10  $\mu\text{L}$  of the fiducial stock to 90  $\mu\text{L}$  of 2X SSC-T. Label tube *Dilution 1 (D1)*.
  - Cover *Dilution 1* and leave at room temperature for 10 minutes protected from light.
  - After 10 minutes, vortex and quick spin *Dilution 1* and dilute to the final working concentration (0.00015%) by adding 15  $\mu\text{L}$  of *Dilution 1* to 985  $\mu\text{L}$  of 2X SSC-T for a final volume of 1000  $\mu\text{L}$  working solution (sufficient for 4 slides). Do not make less than 1000  $\mu\text{L}$  due to the risk of clumping.

The optimal concentration may differ for some tissue types and need to be empirically determined. For this protocol, **start with the default concentration of 0.00015% and adjust the concentration as needed**. See [Top 3 Tips for Successful CosMx SMI Single-cell Spatial Runs at 1000 plex](#) for additional guidance.

3. Remove each slide from 1X PBS and gently tap on a clean Kimwipe to remove excess buffer. Lay slide horizontally in staining tray.

4. **Immediately before applying fiducials to slides, vortex tube for 1 additional minute.** Vortex fiducial working solution for 30 seconds between applications to slides to keep fiducials in suspension and ensure consistent concentration across all slides.
5. Apply up to 250  $\mu\text{L}$  of the fiducial working solution, ensuring the solution covers glass and tissue within the incubation frame. Fiducials must be present on the glass within the scan area for consistent focusing during the instrument run. Gently move tray side to side as needed to ensure that the fiducial solution covers the entire scan area, including glass.

A pipette tip can also be used to spread buffer over tissue within incubation frame. Use a small-volume pipette tip and carefully lay the tip horizontally on top of the incubation frame. Gently roll the tip to spread the buffer, until the tissue is completely covered. Be careful to avoid touching the tissue with the pipette tip.

6. Incubate covered in staining tray for 5 minutes at room temperature ([Figure 34](#)).



Figure 34: Incubate covered for 5 minutes

7. After fiducial incubation, *one slide at a time*, gently tap slides on a clean Kimwipe to remove excess solution and transfer slides to staining jar containing fresh 1X PBS.

8. Wash slides in staining jar with 1X PBS for 5 minutes ([Figure 35](#)). During PBS wash, prepare staining jars for next step. Proceed to next steps immediately.




Figure 35: PBS Wash

**Fresh Frozen Post-Fixation (20 minutes)****Post-Fixation (20 minutes)**

You will need the following materials and reagents for this step: **staining jars, 10% NBF, NBF Stop Buffer (Tris-Glycine Buffer)** and **1X PBS** (see [Prepare RNA Fresh Frozen Assay Reagents on page 56](#) for more information).

**WARNING:** Use of appropriate personal protective equipment is advised. Used NBF Stop Buffer contains NBF and must be disposed of in the same manner as NBF.

Post-fix the tissue by performing the following washes ([Figure 36](#)).

1. Transfer slides to 10% NBF and incubate for 1 minute at room temperature.
2. **Immediately transfer the slides** to the first NBF Stop Buffer and wash for 5 minutes. Transfer slides to a second staining jar containing NBF Stop Buffer and wash for 5 minutes.
3. Transfer slides to 1X PBS for 5 minutes. Slides can sit in 1X PBS while NHS-Acetate mix is prepared.




Figure 36: NBF Post Fix

4. During PBS wash, remove RNase Inhibitor, CosMx RNA probe mix and add-on probes or stand-alone custom panel from -20°C and thaw on ice.
5. Remove Buffer R from 4°C and bring to room temperature.

## NHS-Acetate Preparation and Application (25 minutes)

You will need the following materials and reagents for this step: **analytical scale, 2.0 mL centrifuge tube, staining jar, Sulfo-NHS-Acetate powder, NHS-Acetate buffer** (Slide Prep Kit Box 1, 4°C), and **2X SSC**.

- Prepare 100 mM NHS-Acetate mixture (200 µL/sample) **immediately before** you are ready to apply it to the tissue:
  - If not already done, prepare four 25-mg aliquots (each for 4-slide preparation) or six 15-mg aliquots (each for 2-slide preparation): Bring Sulfo-NHS-Acetate powder from -20°C to room temperature, to prevent condensation. Weigh out the powder directly into 2.0 mL centrifuge tubes. Label tubes with the exact weight, seal tubes with parafilm, and store at -20°C in desiccant. Refer to the Prepare Reagents section of this protocol.
  - Calculate the amount of NHS-Acetate buffer to add to the Sulfo-NHS-Acetate powder by multiplying the weight of Sulfo-NHS-Acetate powder in mg by 38.5. (For example, for 25.0 mg of Sulfo-NHS-Acetate powder,  $25.0 * 38.5 = 962.5$  µL of buffer to add.)
  - Add NHS-Acetate buffer directly to aliquoted powder immediately before applying to the tissue. Slowly pipette up and down to mix. Bubbles may form. Do not fully dispense liquid from pipette while mixing.

**IMPORTANT:** Reconstitute Sulfo-NHS-Acetate immediately before use. Do not store, as the NHS ester moiety readily hydrolyzes and becomes non-reactive. Discard any unused reconstituted reagent.

- Perform the following steps **one slide at a time** to prevent the tissue from drying out:
  - Remove slide from 1X PBS, gently tap slide on a clean Kimwipe to remove excess buffer, and transfer to a clean staining tray.
  - Apply 200 µL of NHS-Acetate mixture onto the tissue within the incubation frame. Gently rock the tray side to side as needed to ensure that the NHS-Acetate solution covers the entire tissue.
  - Repeat with remaining slides and incubate covered in staining tray for 15 minutes at room temperature ([Figure 37](#)).
- Following incubation, tap off excess liquid and wash slides in 2X SSC for 5 minutes ([Figure 38](#)).
- Repeat 2X SSC wash for a total of 2 washes. Leave slides in second 2X SSC wash while hybridization mix is prepared.



Figure 37:  
Incubate 15  
minutes



Figure 38: Two 5-  
minute 2X SSC  
washes

**Fresh Frozen** *In Situ Hybridization (overnight)***In Situ Hybridization (overnight)**

**Do not begin in situ hybridization step until within 16-18 hours of Day 2 start time.** Until then, store slides protected from light in 2X SSC up to 1 hour at room temperature or up to 6 hours at 4°C.

**IMPORTANT:** Take care to maintain nuclease-free conditions. Areas should be cleaned thoroughly with RNase AWAY after probe mix formulation (RNaseZap is only effective for enzymes, not oligos, and should not be used in place of RNase AWAY). Alternatively, mixes can be made in PCR workstations that have been decontaminated with UV light. Gloves should also be changed after handling any probe mixes to avoid cross-contamination.

**IMPORTANT:** If preparing **human or mouse neural tissue**, rRNA is used in this step ([Table 29](#)). The rRNA marker is shipped and stored at -80°C with the Neuroscience Cell Segmentation Kit. No other components of the Kit are used in this step.

**IMPORTANT:** If using a **custom RNA Add-on**, the custom RNA Add-on replaces the standard panel Add-on for the 100- and 1000-plex assay. (The Human 6K Discovery and WTX Panels do not include an Add-on.)

**IMPORTANT:** If using a **custom panel**, email [AtoMx.KitAdmin@Bruker.com](mailto:AtoMx.KitAdmin@Bruker.com) **at least one business day prior to a run** so that the custom kit can be added to the CosMx SMI Control Center.

You will need the following materials and reagents for this step: **hybridization oven, hybridization tray, incubation frame covers, thermal cycler, ice, Buffer R, CosMx RNA Probe Mix, custom or standard Add-On (-20°C), rRNA Marker (if applicable), RNase Inhibitor (-20°C), and DEPC-treated water.**

**Prepare buffers:** Warm Buffer R to room temperature before opening. Thaw probe mix, add-on probes (if applicable) and rRNA Marker (if applicable) on ice. Before use, mix probes thoroughly by pipetting up and down 3-5 times. **Do not vortex probes.** Once thawed, probes can be refrozen at -20°C up to 5 times or refrigerated at 4°C for up to 6 months.

**Set the hybridization oven temperature to 37°C** according to product instructions. If your chamber is light-permeable, minimize light exposure (e.g., by wrapping the lid in aluminum foil).

1. Pre-heat thermal cycler and lid to 95°C.
2. Remove an incubation frame cover and clean with ethanol. Dry with a clean Kimwipe and visually inspect the cover for dust. Use a new Kimwipe as needed to remove any dust. Lay incubation frame cover on a clean Kimwipe until use.
3. For each probe mix, flick to mix then centrifuge. **Do not vortex probes.**
4. Denature CosMx RNA probe mixes (RNA Probe Mix, RNA Add-On or Custom Add-On, plus rRNA Marker if

preparing neural tissue), by transferring total volumes needed for assay from stock tubes into clean 0.2 mL PCR tubes. **Probes, Add-ons, and rRNA Marker should be kept separate during denaturing.**

**IMPORTANT:** Pipette accurately. When preparing 4 slides, there is no excess Probe Mix or Add-on.

5. Heat at 95°C for 2 minutes on a thermal cycler with heated lid. Immediately transfer to ice for at least 1 minute to crash cool.
6. Immediately before preparing hybridization mix, flick to mix and centrifuge tubes.
7. Make hybridization solution following [Table 28](#) for human or mouse **non-neural tissue** or [Table 29](#) for human or mouse **neural tissue**. **Prepare hybridization mix no more than 20 minutes before tissue application.**

Table 28: Hybridization Solution for **Non-Neural Tissue**

	Denatured RNA Probe Mix	Denatured Add-on* (if applicable)	RNase Inhibitor	Buffer R	DEPC-treated water	Final Volume
Off-the-shelf Panel (2 slides)	32 µL	16 µL	3.2 µL	256 µL	up to final volume of 320 µL	320 µL
Off-the-shelf Panel (4 slides)	64 µL	32 µL	6.4 µL	512 µL	up to final volume of 640 µL	640 µL
Custom Stand-alone Panel (2 slides)	16 µL**	-	3.2 µL	256 µL	up to final volume of 320 µL	320 µL
Custom Stand-alone Panel (4 slides)	32 µL**	-	6.4 µL	512 µL	up to final volume of 640 µL	640 µL

\* If using a 100- or 1000-plex assay with a custom Add-on, the custom Add-on replaces the standard Add-on that is part of the assay. (WTX and 6K Discovery Panels do not include a standard Add-on.) If no Add-on is used, use DEPC-water to reach the correct final volume.

\*\* Custom stand-alone panels are supplied at 2X concentration relative to off-the-shelf RNA Panels.

**Fresh Frozen** *In Situ Hybridization (overnight)*Table 29: Hybridization Solution for **Neural Tissue**

	Denatured RNA Probe Mix	Denatured Add-on * (if applicable)	Denatured rRNA Marker**	RNase Inhibitor	Buffer R	DEPC-treated water	Final Volume
Off-the-shelf Panel (2 slides)	32 $\mu$ L	16 $\mu$ L	12.8 $\mu$ L	3.2 $\mu$ L	256 $\mu$ L	up to final volume of 320 $\mu$ L	320 $\mu$ L
Off-the-shelf Panel (4 slides)	64 $\mu$ L	32 $\mu$ L	25.6 $\mu$ L	6.4 $\mu$ L	512 $\mu$ L	up to final volume of 640 $\mu$ L	640 $\mu$ L
Custom Panel with add-on (2 slides)	16 $\mu$ L***	16 $\mu$ L	12.8 $\mu$ L	3.2 $\mu$ L	256 $\mu$ L	up to final volume of 320 $\mu$ L	320 $\mu$ L
Custom Panel with add-on (4 -slides)	32 $\mu$ L***	32 $\mu$ L	25.6 $\mu$ L	6.4 $\mu$ L	512 $\mu$ L	up to final volume of 640 $\mu$ L	640 $\mu$ L

\* If using a 100- or 1000-plex assay with a custom Add-on, the custom Add-on replaces the standard Add-on that is part of the assay. (WTX and 6K Discovery Panels do not include a standard Add-on.) If no Add-on is used, use DEPC-water to reach the correct final volume.

\*\* The rRNA marker is the only segmentation marker added in the overnight hybridization. All other segmentation markers are added during the Nuclear and Cell Segmentation step.

\*\*\* Custom stand-alone panels are supplied at 2X concentration compared to off-the-shelf RNA Panels.

8. **Clean all equipment** and benchtop with RNase AWAY and allow to dry; or rinse with DEPC-treated water. The hybridization chamber can be a key source of contamination by oligos. Arrange fresh Kimwipes on bottom of the chamber. Change gloves and clean workspace with RNase AWAY.
9. Wet the Kimwipes with 2X SSC or DEPC-treated water. Take care that the Kimwipes and 2X SSC do not contact the slides. Kimwipes should be thoroughly wet, but without standing buffer.
10. To prevent the tissue from drying, perform the following steps **one slide at a time**.
  - Remove slides from 2X SSC, and gently tap slide to remove excess liquid.
  - Carefully remove the polyester frame backing (with the center square removed) from the incubation frame to expose the top adhesive layer of the incubation frame ([Figure 39](#)). **Ensure that the incubation frame does not lift from the slide when removing the polyester frame backing.**



Figure 39: Remove Polyester Frame backing from Incubation Frame

- Lay the slide flat on a clean surface and add 150  $\mu$ L of hybridization solution directly to the tissue within the incubation frame.

Add the hybridization solution at the edge of the tissue opposite of the slide label within the frame. Applying the incubation frame cover in the next step will help move the hybridization solution across the tissue.

Avoid introducing bubbles by leaving a small residual volume in the pipette tip. If a bubble is introduced, carefully aspirate it out using a low-volume pipette tip. When removing air bubbles, removing small amounts of hybridization solution (as long as sufficient solution remains to cover the tissue) is preferable to leaving bubbles.

- Carefully apply incubation frame cover ([Figure 40](#)). Start by setting one edge of the cover down on the incubation frame edge, then gradually lay down the rest of the cover. The tab on the incubation frame cover should face the slide label. As it's lowered, the frame cover should naturally adhere to the incubation frame; no additional pressure around the frame is needed. **Do not press the center of the cover as it could damage the tissue.**




Figure 40: Apply Incubation Frame Cover

- Place the slide horizontally into the hybridization tray ([Figure 41](#)).




Figure 41: Hybridization Tray

- Repeat step 10 for each slide.
  - Ensure that there is a good seal between the incubation frame and the slide, and the incubation frame cover and the frame by checking each slide for any leaks.
11. Close hybridization chamber, insert tray into oven, and clamp tray into place. Incubate at 37°C overnight (16 – 18 hours) ([Figure 42](#)).

**IMPORTANT!** Continue to [Add Enzymes to Buffer 4 on page 74](#) before ending Day 1 slide preparation.




Figure 42: Incubate overnight at 37°C

If your oven does not seal (with a gasket) you may seal your hybridization chamber in a zip-lock bag to simulate a sealed chamber. Chambers sealed in this manner should be tested to ensure they maintain humidity for 24 hours (slides do not dry out) prior to use. Unsealed conditions lead to evaporation of the hybridization solution.


**Fresh Frozen** *In Situ Hybridization (overnight)***Add Enzymes to Buffer 4**

To reduce oxidative damage to the tissue and probes and improve overall run quality, Bruker recommends adding the active enzymes Catalase and Pyranose Oxidase to instrument Buffer 4 the day before starting a new instrument run.

These instructions cite the different sizes of Instrument Buffer Kits. To review the kit contents, refer to the [CosMx SMI Instrument User Manual](#).

**NOTE:** When supplies are limiting, two Small enzyme aliquots may be shipped in lieu of one Medium enzyme aliquot. Confirm the aliquot size on the tube label before proceeding with these steps.

1. Remove lyophilized Catalase and Pyranose Oxidase vials from 4°C. Centrifuge the vials and resuspend enzymes in DEPC-treated water:
  - For Small enzyme aliquots, add 250 µL of DEPC-water to each tube. Vortex then centrifuge.
  - For Medium enzyme aliquots, add 400 µL of DEPC-water to each tube. Vortex then centrifuge.
2. Add Catalase to Buffer 4:
  - For a Small Buffer 4 (58 mL), add 250 µL resuspended Catalase. For Medium Buffer 4 (90 mL), add 400 µL Catalase. **Avoid precipitate when pipetting from the enzyme tube.**
3. Add Pyranose Oxidase to Buffer 4:
  - Due to the increased risk of precipitate in this enzyme, ensure the resuspended enzyme has been vortexed and centrifuged, and any precipitate has been pelleted at the bottom of the tube.
  - For a Small Buffer 4 (58 mL), add 125 µL resuspended Pyranose Oxidase. For Medium Buffer 4 (90 mL), add 200 µL resuspended Pyranose Oxidase. **Avoid precipitate when pipetting from the enzyme tube.**

 **IMPORTANT:** Precipitate in Buffer 4 could clog the instrument fluidic lines and cause run failure.

4. Mark Buffer 4 bottle after enzymes have been added. Dispose of excess enzymes according to laboratory guidelines.
5. Leave Buffer 4 sealed on the bench until ready to load onto the instrument.

If not already done, prepare 40 mL formamide aliquots following the instructions in the Prepare Reagents section of this protocol and store overnight at 4°C.

**Day 2: Perform Stringent Washes (90 minutes)**

You will need the following materials and reagents for this step: **water bath**, **4X SSC**, **100% formamide**, and **2X SSC** (see [Prepare RNA Fresh Frozen Assay Reagents on page 56](#)).

**WARNING:** Use of personal protective equipment is advised as formamide is hazardous.

1. Set water bath to 37°C. Remove nuclear stain and cell segmentation kits from the freezer and thaw on ice.
2. **Warm 100% formamide in the 37°C water bath for at least 30 minutes before opening.** Once at temperature, prepare stringent wash directly in staining jars by mixing equal parts 4X SSC and 100% formamide.

Ensure you have sufficient buffer to completely cover the tissue on all slides without submerging the slide labels. Contact with buffer may make slide labels illegible.

3. Preheat staining jars containing freshly-prepared stringent wash in 37°C water bath. It takes ~30 min to preheat.
4. If nearing the 18-hour maximum overnight incubation time, while jars are preheating, transfer slides to 2X SSC.
5. Perform the following steps **one slide at a time** to prevent the tissue from drying.
  - With a clean pair of forceps, carefully remove the incubation frame cover from the incubation frame. Dip slide into 2X SSC as needed to avoid tissue drying. If cover will not come off without removing incubation frame, remove the frame and cover. The frame can be reapplied in a later step.
  - Place slide into 2X SSC and continue to the next slide, cleaning forceps with ethanol between slides.
6. Once both jars have pre-heated to 37°C and all incubation frame covers have been removed, perform the washes shown below ([Figure 43](#)). After the last wash, the slides can be stored in 2X SSC for up to one hour.
  - Gently tap each slide on a Kimwipe to remove excess 2X SSC and place in the first stringent wash for 25 min.
  - Transfer slides to the second stringent wash for 25 min.
  - During second stringent wash, begin preparing reagents for nuclear and cell segmentation staining in the next section.
  - Following stringent washes, immediately transfer slides to 2X SSC and wash for 2 minutes. Transfer slides to a second jar of 2X SSC and wash for 2 minutes. Leave slides in 2X SSC until reagents have been prepared for nuclear and cell segmentation staining.



Figure 43: Perform stringent wash

**IMPORTANT:** Anything coming into contact with hybridization solution (which contains probes), such as containers for stringent wash solution and 2X SSC, needs to be exclusive for this purpose and thoroughly cleaned with RNase AWAY, as probes may contaminate subsequent assays. Use separate staining jars for different probe mixes. Staining jars should be cleaned with RNase AWAY before use.

**Fresh Frozen Nuclear and Cell Segmentation Staining (2 hours)****Nuclear and Cell Segmentation Staining (2 hours)**

You will need the following materials and reagents for this step: **incubation frames, staining jars, 1X PBS, Blocking Buffer (4°C), Nuclear Stain stock (-80°C), and Segmentation Marker Kit (-80°C)** (see [Prepare RNA Fresh Frozen Assay Reagents on page 56](#)).

1. Prepare four staining jars of 1X PBS.
2. Prepare Nuclear Stain (220  $\mu$ L per slide):
  - Vortex, then centrifuge thawed Nuclear Stain Stock for at least 1 minute to bring the solution to the bottom of the vial and precipitate insoluble particles.
  - Dilute the Nuclear Stain stock 1:40 as described in [Table 30](#). *Do not pipette from the bottom of the vial.*

Table 30: Prepare Nuclear Stain

	Nuclear Stain Stock	Blocking Buffer	Final Volume
2 slides	11 $\mu$ L	429 $\mu$ L	440 $\mu$ L
4 slides	22 $\mu$ L	858 $\mu$ L	880 $\mu$ L

- Mix by pipetting up and down 3-5 times. Do not vortex.
3. If a new incubation frame is needed, perform the following steps **one slide at a time**.
    - Remove slide from 2X SSC, gently tap slide on a clean Kimwipe. Using a clean Kimwipe, dry the surface of the slide that will come into contact with the incubation frame.

**i IMPORTANT:** Avoid wiping the slide within the scan area as this could remove the fiducials needed for imaging. See scan area template on the flow cell assembly tool or follow the guidelines in [Appendix I: CosMx SMI Tissue Sectioning Guidelines on page 84](#).

- Carefully apply a new incubation frame following the previous instructions in this protocol. Ensure that the frame is well adhered to the slide by gently pressing around the frame with clean forceps.
4. Using a clean Kimwipe, carefully wick excess buffer from around the incubation frame as needed. Be careful to not touch the area inside of the incubation frame.
  5. Lay slide horizontally in staining tray and slowly apply up to 200  $\mu$ L of Nuclear Stain directly to tissue. Gently move tray side to side as needed to ensure that the buffer covers the entire tissue.
  6. Repeat with remaining slides and cover tray.



Figure 44: Cover tray and incubate for 15 minutes

7. Incubate slides for 15 minutes at room temperature protected from light (Figure 44).
8. During nuclear stain incubation, prepare Segmentation Mix (**non-neural tissue: blue header** (Table 31), **neural tissue: magenta header** (Table 31)). Flick each tube to mix and centrifuge before use. **Do not vortex mix.**

**IMPORTANT:** Pipette accurately. When preparing 4 slides, there will be no excess Segmentation Mix.

Table 31: Segmentation Mix for non-neural tissue

	CD298/ B2M Marker Mix	Optional PanCK/ CD45 Marker Mix*	Optional a la carte Marker Mix*	Blocking Buffer	Final Volume
2 slides	8 $\mu$ L	16 $\mu$ L	16 $\mu$ L	360 $\mu$ L	400 $\mu$ L
4 slides	16 $\mu$ L	32 $\mu$ L	32 $\mu$ L	720 $\mu$ L	800 $\mu$ L

\* If not adding PanCK/CD45 or optional à la carte markers, add Blocking Buffer to reach final volume.

Table 32: Segmentation Mix for neural tissue

	GFAP	Histone	Blocking Buffer	Final Volume
2 slides	16 $\mu$ L	16 $\mu$ L	368 $\mu$ L	400 $\mu$ L
4 slides	32 $\mu$ L	32 $\mu$ L	736 $\mu$ L	800 $\mu$ L

The rRNA Segmentation marker is added during the overnight hybridization. All other segmentation markers are added at this step.

9. After nuclear stain incubation, **remove slides one at a time** from staining tray, gently tap slide on a clean Kimwipe to remove excess buffer, and transfer slide to 1X PBS (Figure 45).




Figure 45: Wash in 1X PBS for 5 minutes

10. Wash slide for 5 minutes in 1X PBS.
11. During PBS wash, add 2X SSC or DI water to the staining tray. Do not overfill. The water level should be well below the slides to avoid cross-contamination.
12. Following PBS wash, perform the following steps **one slide at a time** to prevent tissue drying:
  - Remove slide from 1X PBS and gently tap slide on a clean Kimwipe to remove excess PBS.
  - Lay slide horizontally in staining tray and apply up to 200  $\mu$ L of Segmentation Mix directly to tissue. Gently move tray side to side as needed to ensure that the mix covers the entire tissue.
  - If needed, an incubation frame cover can be placed over the incubation frame to ensure that the mix completely covers the tissue.
  - Adjust volume to add as needed for tissues of varying sizes. The **Segmentation Mix needs to completely cover the tissue** but does not need to completely fill the incubation frame.
13. Repeat with remaining slides and cover tray (Figure 46).

**Fresh Frozen** Nuclear and Cell Segmentation Staining (2 hours)

14. Incubate slides for 1 hr at room temperature protected from light. If planning to load the instrument immediately after slides are prepared (recommended), remove the Imaging Tray from 4°C and leave at room temperature until loading the instrument. Confirm that the Imaging Tray has not expired.



Figure 46: Cover tray and incubate for 1 hr.

15. Following Segmentation Mix incubation, transfer slides to 1X PBS and wash for 5 minutes ([Figure 47](#)).
16. Repeat wash 2 times for a total of 3 PBS washes.



Figure 47: Wash 3 times in 1X PBS.

17. If samples will be loaded onto the instrument the same day (recommended), continue to [Flow Cell Assembly](#) [Flow Cell Assembly on page 80](#).

If samples must be stored overnight and loaded onto the instrument the next day, remove the incubation frame following the previous instructions in this protocol, transfer slides to fresh 2X SSC, and store according to **Safe Storage Guidelines for RNA Slides** [Safe Storage Guidelines for RNA Slides on page 79](#).

Once the flow cell has been applied to the slide, proceed to run the RNA assay on the CosMx SMI instrument.

Refer to the [CosMx SMI Instrument User Manual \(MAN-10161\)](#) for instructions to operate the instrument.

## Safe Storage Guidelines for RNA Slides

After processing, do not store slides dry.

Store slides for up to 6 hours protected from light and submerged in 2X SSC at room temperature.

If needed, slides can be stored overnight, protected from light and submerged in 2X SSC at 4°C, if needed.

Slides may be stored longer than overnight, but RNA counts and staining efficiency will decrease as a function of days stored. For best results, minimize storage time between slide preparation and loading on the CosMx SMI instrument.

Slides must be stored in the dark.

## Flow Cell Assembly

The CosMx SMI flow cell enables input of a tissue slide into the CosMx SMI instrument for spatial profiling. It affixes to a 3 inch x 1 inch standard pathology grade slide with mounted tissue, creating an imageable fluidic channel. The CosMx SMI reagents required for the cycling chemistry are flowed across the tissue through the formed channel using the integrated fluidic input and output ports.

The provided flow cell assembly tool is a clamshell design that applies uniform force to adhere the flow cell coverslip onto the prepared slide ([Figure 48](#)).

Before beginning flow cell assembly, verify tissue placement using the template provided on the flow cell assembly tool ([Figure 49](#), not to scale). Tissue sections must be centered within the Scan Area of the slide (the green rectangle) and be no larger than **20 mm long by 15 mm wide**. For best performance, there should also be some bare glass (not covered by tissue) in the scan area. For examples of tissue placement best practices, see [Appendix I: CosMx SMI Tissue Sectioning Guidelines on page 84](#). If needed, remove excess tissue using a clean razor blade.




Figure 48: Flow Cell Assembly Tool




Figure 49: Tissue Scan Area

### Flow Cell Preparation

1. Clean the benchtop with RNase AWAY or 70% ethanol. RNase AWAY is recommended for RNA slide preparation as 70% ethanol does not adequately remove both nucleic acid and nuclease contaminants.
2. Prepare flow cell assembly tool by cleaning the stage with ethanol or isopropanol, and then blow any dust from the tool using an air blower (not compressed air). Be careful as glass fragments and slivers may be present. Clean the top compression pad with DI water. Do not clean the top compression pad with ethanol or isopropanol, as this could damage the pad.

**IMPORTANT:** Do not use compressed air on the flow cell assembly tool or the instrument. An airblower such as [Giottos AA1910 Medium Rocket Air Blaster \(6.6\)](#) or [Camkix Keyboard Cleaning Kit](#) is recommended.

3. Inspect the flow cell coverslip for any damage such as cracks or chips and record the coverslip serial number. This is the flow cell barcode needed when loading the instrument.
4. With a clean Kipwipe and ethanol, carefully clean the flow cell coverslip to remove any debris. Once clean, place the flow cell coverslip on a new Kimwipe until ready to use.

### Flow Cell Assembly

Assemble one flow cell at a time to prevent the tissue from drying out.

1. Remove the slide from storage buffer using clean forceps.
2. Carefully remove the incubation frame, if not already removed, and tap off excess buffer.
3. Dry the back of the slide and, using the template on the flow cell assembly tool, carefully dry the area around the tissue where the flow cell coverslip will adhere to the slide.
  - Be careful not to wipe within the Scan Area (shown in green) as this could remove fiducials required for on-instrument imaging.
  - If using an adhesive slide label, ensure the label is less than 295  $\mu\text{m}$  thick and is not folded over on itself. If label extends over frosted label area of the slide, carefully trim using a clean razor blade. Labels over the maximum thickness or which are not properly adhered may result in slide or flow cell damage during flow cell assembly and/or instrument loading.
4. Lower the tailgate (tool marker 1) on the flow cell assembly tool ([Figure 50](#)). Hold the labeled end of the slide and insert, tissue side up, non-labeled end first, into the tool through the bottom opening. The slide is fully inserted once the non-labeled edge contacts the back of the slide stage (tool marker 2).




Figure 50: Lower Tailgate

5. Raise the tailgate to location 3 (tool marker 3) to secure the slide.
6. Apply the flow cell coverslip:
  - Use the air blower, if needed, to remove any dust from both sides of the coverslip immediately before applying.
  - Carefully remove the adhesive backing from the flow cell coverslip:
    - Hold the coverslip so that the serial number is readable. The backing is on the opposite side.

- Use clean forceps or gloved hands to hold the tab attached to the adhesive backing.
  - Slowly peel the adhesive backing away from the coverslip.
7. Place the flow cell coverslip onto the slide, adhesive side down, within the slide stage area.
- Hold the flow cell coverslip along the long edges and place the coverslip carefully onto the slide, adhesive side down, keeping the coverslip parallel to the slide.
  - To confirm none of the edges of the flow cell coverslip are lifted or are catching on the tool, lightly tap on the four corners of the coverslip (marked with stars in [Figure 51](#)). The air gaps should be reduced and signs of adhesion (dark patches) should be present along the edges.



Figure 51: Flow Cell Placement

8. To complete assembly, swing the arm of the flow cell assembly tool down until the latches on either side have engaged ([Figure 52](#)). Once engaged, both green release buttons pop out.
9. After the latches have engaged, release the arm by pressing the two buttons on the front of the assembly tool arm.
10. The newly assembled flow cell can be removed by reversing the steps of inserting the slide. Pull down the tailgate and gently remove the flow cell.





Figure 52: Fully Engaged Latch

11. **Check the slide and flow cell for any cracks or damage.** Pipette 200  $\mu$ L of storage buffer (2X SSC) into the flow cell port to ensure the tissue does not dry out:
- Place the pipette tip directly over one of the fluidic ports on the flow cell. Avoid putting pressure on the flow cell port as this could result in cracks at the port.
  - Slowly press the plunger and allow buffer to fill the chamber and cover the tissue.
  - Ensure there are no bubbles within the flow cell as this could result in imaging failure. Additional buffer may be flowed through the fluidic ports as needed to push out bubbles.
  - Once the tissue has been covered and the flow cell chamber is full, remove the pipette tip without releasing the plunger and dispose of extra buffer.
  - Use a clean Kimwipe to wick away excess buffer from around the flow cell ports, being careful not to touch the port with the Kimwipe.

12. Once assembled, place the flow cells into the clean staining tray, protected from light, until ready to load the instrument.

 **IMPORTANT:** If the flow cell or slide is cracked, refer to [Troubleshooting on page 92](#) for additional guidance or contact [support.spatial@bruker.com](mailto:support.spatial@bruker.com).

Continue to the [CosMx SMI Instrument Manual \(MAN-10161\)](#) for instructions to load the flow cell into the instrument and begin data acquisition.

## Appendix I: CosMx SMI Tissue Sectioning Guidelines

When preparing, sectioning, and storing FFPE blocks for use in the CosMx SMI Protein and RNA assays, care should be taken to preserve sample integrity in all steps. The integrity of FFPE samples can be impacted by many factors, including time from excision to fixation, storage conditions, tissue type, and sample age. It is important to take such factors into consideration when selecting samples for the CosMx SMI assay. Samples with poor integrity are likely to give low signal, particularly in the CosMx SMI RNA assay.

CosMx SMI has been validated for samples up to 3 years old prepared from tissues with a cold ischemic time of less than 1 hour using 10% NBF or similar fixative. In general, for best results, do not use FFPE blocks older than 10 years. Assay performance, particularly for RNA, will be influenced by tissue block age and treatment conditions such as cold/warm ischemic time, fixative, and storage.

For additional guidance sample sectioning, see [Sample Sectioning Tips and Tricks for CosMx SMI and GeoMx DSP Experiments \(MAN-10175\)](#) in the NanoString University Document Library.

### Selecting FFPE Blocks

FFPE blocks should meet the following criteria for the best performance with the CosMx SMI assay.

- Blocks should be fixed in 10% neutral-buffered formalin for 18 to 24 hours at room temperature. This applies to tissues 5 mm in thickness. Thicker tissues have not been tested by Bruker and may require longer fixation times.
- Tissues should be fixed immediately after excision for best results. Up to one hour post-excision is acceptable.
- Tissues should be thoroughly dehydrated in ethanol gradients prior to embedding in paraffin.
- FFPE blocks should be stored at room temperature and ambient humidity.
- For best results, do not use FFPE blocks that are older than 10 years.

### Sectioning FFPE Blocks

The following are general guidelines for sectioning FFPE blocks for optimal CosMx SMI assay performance. This is not meant to be an exhaustive guide on sectioning. Please refer to your local pathologist, histologist, or core facility for training on sectioning.

- Always discard the first few sections from the block face.
- Tissue sections should be sectioned **5  $\mu\text{m}$  thick** and mounted on the label side of Leica BOND Plus slides or VWR Superfrost Plus Micro Slides (Figure 53) (figure not to scale; see the template on the Flow Cell Assembly Tool for a to-scale template). Leica BOND Plus slides are required for BOND RX/RX<sup>m</sup> semi-automated slide preparation and are preferable in manual slide preparation for tissues with poor adherence.
- Sections should be mounted in the center of the tissue scan area (shown in green in Figure 53), fitting within the dimensions 20 mm x 15 mm, to allow adequate room for flow cell coverslip adhesives at the edges of the scan area. If mounting multiple sections per slide, ensure that all tissues are at least 2–3 mm apart and still contained within the scan area. For best performance, there should be some bare glass (not covered by tissue) in the scan area (Figure 54).
- It is important to avoid any scratches and folds in the tissue section. These scratches and folds can be magnified by the subsequent slide washes on the CosMx SMI instrument resulting in tissue loss.
- Any water trapped under the wax or tissue section should be removed by gently touching a folded Kimwipe onto the corner of the wax section. The Kimwipe should not contact the tissue.
- It is recommended to use sections within two weeks of slide mounting for best results. Older sections (1–2 months) may yield satisfactory results, but this may be tissue or block dependent and should be tested empirically. Slides should be stored at room temperature in a desiccator or at 4°C prior to processing.




Figure 53: Tissue scan area (not to scale)




Figure 54: While all 3 examples are acceptable, Slide 1 shows optimal tissue placement; Slide 2 shows good tissue placement, and Slide 3 gives an example of a tissue that should be trimmed if possible to allow more visible glass in the scan area.

If sections are larger than the allowable size, or placed off-center, excess tissue must be scraped away with a clean razor blade. Tissue may be scraped when slides are dry at the start of slide preparation, before applying the incubation frame during slide preparation, or just before applying the flow cell coverslip after slide preparation. The optimal practice (minimizing tissue pulling, folding, or detaching) may depend on your sample type. Applying the incubation frame over tissue could result in tissue damage when the frame is removed or poor sealing of the incubation frame.

**i IMPORTANT:** The CosMx SMI instrument will only image the area inside the flow cell chamber, the tissue scan area. If the tissue section is outside of the scan area, it will not be imaged.

**To improve tissue adherence,** bake FFPE slides at 37°C for 2 hours after sectioning. After baking, dry at room temperature overnight prior to use or storage.

### Selecting Fresh Frozen Blocks

- Tissues should be selected that are known to have been snap frozen in liquid nitrogen as quickly as possible. Alternative freezing media may include isopentane pre-cooled with liquid nitrogen or isopentane cooled with dry ice.
- Any buffers used to wash or temporarily store tissues before fixation should be free of nuclease contamination.
- Frozen tissues should be embedded in Optimal Cutting Temperature media (OCT) before sectioning.
- Blocks embedded in OCT should be stored at -80°C.

### Sectioning Fresh Frozen Blocks

- Always discard the first few sections from the block face.
- Tissue sections should be sectioned **5 µm thick** and mounted on the label side of Leica BOND Plus slides or VWR Superfrost Plus Micro Slides ([Figure 53](#)) (figure not to scale; see the template on the Flow Cell Assembly Tool for a to-scale template). Leica BOND Plus slides are required for BOND RX/RX<sup>m</sup> semi-automated slide preparation and are preferable in manual slide preparation for tissues with poor adherence. During sectioning, it is important to cut across the tissue with a smooth, consistent turn of the hand wheel. Blocks may be sectioned up to 10 µm thickness; however, the instrument will only image the 5 µm closest to the slide.

**i IMPORTANT:** Cryostat temperature should be set to -20°C. Place fresh frozen block inside of cryostat for a minimum of 30 minutes to equilibrate to temperature. Temperature may need to be adjusted +/- 5°C to optimize sectioning.

- Sections should be mounted in the center of the tissue scan area (shown in green in [Figure 53](#)), fitting within the dimensions 20 mm x 15 mm, to allow adequate room for flow cell coverslip adhesives at the edges of the scan area. If mounting multiple sections per slide, ensure that all tissues are at least 2–3 mm apart and still contained within the scan area. For best performance, there should be some bare glass (not covered by tissue) in the scan area ([Figure 54](#)).
- It is critical to avoid any scratches and folds in the tissue section. These scratches and folds can be magnified by the subsequent slide washes on the CosMx SMI instrument resulting in tissue loss. Folds and wrinkles in fresh frozen tissues are highly susceptible to damage during washes and incubation frame removal.
- After sectioning, the exposed block face should be covered with OCT to avoid desiccation of the sample.
- Dry slides at room temperature for 5-10 minutes then store at -80°C with a desiccant.
- Slides can be stored at -80°C for several weeks before use.

## Appendix II: Tissue-Specific Slide Preparation Considerations

The optimal conditions for target retrieval, tissue permeabilization (Proteinase K digestion), and fiducials concentration may differ by tissue type and may need to be empirically determined. Internal testing indicates that most FFPE tissues perform well with the following conditions:

- Target retrieval: 15 min with ER1 at 100°C
- Tissue permeabilization: 3 µg/mL Proteinase K for 30 minutes (in 1X PBS-T for FFPE human brain or 1X PBS for all other tissues)
- Fiducials concentration: 0.0005 - 0.001%

Cell pellet arrays and delicate tissues generally perform better with the following conditions:

- Target retrieval: 8 min with ER1 at 100°C
- Tissue permeabilization: 1 µg/mL Proteinase K for 15 min (cell pellet array) or 30 min (delicate tissues)
- Fiducials concentration: 0.0005-0.001% (for cell pellet array, use 0.001%)

On the following pages, we provide more context and details for tissue-specific slide preparation modifications ([Table 33](#)) and the conditions used by the Technology Access Program (TAP) for a variety of tissues tested ([Table 34](#)). This information is provided as a general resource; not all tissue types have been validated and the optimal conditions for your particular experimental design and samples may differ.

Source: [Top 3 Tips for Successful CosMx SMI Single-cell Spatial Runs at 1000 plex](#)

Table 33: Suggested slide prep modifications based on tissue type and biology. With all tissue types, begin with default experimental conditions and modify as needed on an experiment basis.

Tissue Types	Suggested Conditions
All types of tumors, especially those derived from epithelial tissues such as colon, lung, breast, ovarian, kidney, bladder cancer, and cholangiocarcinoma	<b>Default experimental conditions:</b> Target Retrieval Time: 15 minutes Digestion Buffer Concentration: 3 µg/mL Digestion Time: 30 minutes Fiducial Concentration: 0.001%
All normal solid organs: lymph node, liver, kidney etc	Use default experimental conditions.
Tissues with higher levels of adipose, loose connective tissue, delicate structures; normal breast, normal lung, organoids, retina, airways, cell pellet arrays	Target Retrieval Time: 8 minutes Digestion Buffer Concentration: 1 µg/mL Digestion Time: 30 minutes (15 min for cell pellet arrays) Fiducial Concentration: 0.0005%
Tissues with high density, such as cartilage and bone	Target Retrieval Time: 8 minutes Digestion Buffer Concentration: 1 µg/mL Digestion Time: 30 minutes Fiducial Concentration: 0.0005%
Tissues with low pH, typically associated with stomach-related diseases and cancers	Use default experimental conditions. However, due to the effect of low pH on these tissue types, the RNA counts will typically be lower.
Tissues with feces: normal intestine, colon, and inflammatory bowel disease (IBD)	Use default experimental conditions. However, due to the pre-fixation steps required for these tissue types, the counts will typically be lower.
Tissues that exhibit high autofluorescence, such as placenta and non-human primate (NHP) tissues	Use default experimental conditions during slide preparation. A pre-bleaching configuration with a longer pre-bleaching time may be required. In addition, the use of a Hydrogen Peroxide pre-treatment may be needed (see <a href="#">High Autofluorescence in Tissue (RNA protocol) on page 93</a> ).
Bone marrow should be treated similarly to tumor tissue (row 1) due to the stability of marrow cells	Use default experimental conditions during slide preparation. Due to the decalcification of the tissue, lower counts may be observed.
Human tissue with long ischemic time and samples where poor tissue quality is observed prior to sample preparation	Consider using mouse tissue with the appropriate panel, which may allow shorter ischemic time, to improve data quality.

**Appendix II: Tissue-Specific Slide Preparation Considerations**

Table 34: Sample preparation conditions for tissues tested by TAP.  
 Source: *Top 3 Tips for Successful CosMx SMI Single-cell Spatial Runs at 1000 plex*

Tissue Types	Target Retrieval Time (min)	Digestion Buffer Conc. (µg/mL)	Digestion Time (min)	Fids Conc.	Pre-bleaching Profile Config.
Appendix	15	3	30	0.0005%	C
Bladder Cancer	15	3	30	0.0005%	C
Bone Marrow	15	3	30	0.001%	C
Bone Marrow Cancer	15	3	30	0.0005%	C
Brain	15	3	30	0.0005%	B
Brain (Alzheimer's)	15	3	30	0.0005%	B
Brain Organoids	8	1	15	0.0005%	C
Breast	15	3	30	0.001%	C
Breast Cancer	15	3	30	0.001%	C
Cartilage	8	3	15	0.0005%	C
Cell Pellet Array	8	1	15	0.001%	A
Cholangiocarcinoma	15	3	30	0.0005%	C
Chronic Kidney Disease	15	3	30	0.0005%	C
Colitis	15	3	30	0.0005%	C
Colon	15	3	30	0.001%	C
Colon Cancer	15	3	30	0.0005%	C
Glioblastoma	15	3	30	0.0005%	C
Glioma	15	3	30	0.0005%	C
Ileum	15	3	30	0.0005%	C
Intestine	15	3	30	0.0005%	C
Kidney	15	3	30	0.001%	B
Kidney Cancer	15	3	30	0.001%	B
Kidney Metabolic Disease	15	3	30	0.001%	B
Liver (malignant)	15	3	30	0.001%	C
Liver (normal)	15	3	30	0.001%	B
Lung	15	3	30	0.0005%	C
Lung Cancer	15	3	30	0.0005%	C
Lung Disease	15	3	30	0.0005%	C
Lymph Node	15	3	30	0.001%	C
Lymphoma	15	3	30	0.0005%	C
Mesothelioma	15	3	30	0.0005%	C
Mouse Artery	15	3	30	0.001%	C

Tissue Types	Target Retrieval Time (min)	Digestion Buffer Conc. ( $\mu\text{g}/\text{mL}$ )	Digestion Time (min)	Fids Conc.	Pre-bleaching Profile Config.
Mouse Brain	15	3	30	0.001%	C
Mouse Brain Infection	5	3	30	0.001%	C
Mouse Brain Neurodegeneration	15	3	30	0.001%	C
Mouse Colon Infection	15	3	30	0.001%	C
Mouse Lung Infection	8	3	15	0.0005%	C
Mouse Retina	15	3	30	0.0005%	C
Mouse Spleen HIV	15	3	30	0.0005%	C
NHP Liver Infection	15	3	30	0.0005%	C
NHP Lung Granuloma	15	3	30	0.0005%	C
NHP Lymph Node SIV	15	3	30	0.0005%	C
Organoid	8	1	30	0.0005%	C
Ovarian Cancer	15	3	30	0.001%	C
Pancreas Diabetes	15	3	30	0.001%	C
Pancreatic Cancer	15	3	30	0.001%	C
Pituitary Gland	15	3	30	0.001%	C
Placenta	8	1	15	0.001%	B
Placenta Diabetes	15	3	30	0.0005%	B
Prostate Cancer	15	3	30	0.001%	C
Rectum	15	3	30	0.0005%	C
Retina	15	3	30	0.0005%	C
Skin	15	1	30	0.001%	C
Skin Cancer	15	3	30	0.0005%	C
Skin Lupus	8	3	30	0.0005%	C
Stomach Cancer	15	3	30	0.0005%	C
Tonsil	15	3	30	0.0005%	C

## Troubleshooting

This section includes troubleshooting information for RNA and protein slide preparation, as well as the on-instrument workflow. For additional support, contact [support.spatial@bruker.com](mailto:support.spatial@bruker.com).

### Poor Tissue Adherence to Slide

**Possible Causes:** Inadequate baking time, inherent biology of sample type.

**Suggested Actions:**

- For FFPE tissue optimal performance, after sectioning and prior to use or storage, bake slides at 37°C overnight at an angle no greater than 45°. Alternatively, slides can be baked at 37°C for 2 hours and then dried overnight at room temperature. A polyethylene slide holder ([VWR, 82024524](#)) can be used for overnight drying.
- For the RNA FFPE assay, the day before beginning sample prep, a second baking step is optimal for tissue adherence. Bake sections on slides overnight in a 60°C drying oven, vertically in a slide rack overnight or in a slide holder at a 45° angle.
- Additionally, although both Leica BOND Plus slides and VWR Superfrost Plus Micro Slide, Premium have been validated, better tissue adherence has been observed with the Leica BOND Plus slides which should be used for tissues with known poor adherence. Alternatively, although not yet validated, TOMO<sup>®</sup> Adhesion Microscope Slides, Matsunami Glass ([VWR, 10748-166](#)) have also shown improved tissue adherence in preliminary testing.
- A tissue section adhesive such as EpreDia<sup>™</sup> Tissue Section Adhesive ([Fisher Scientific, 86014](#)) can also be used to improve tissue adherence. The use of an adhesive has not been validated but may improve tissue adherence for some tissue types. Follow the manufacturer's instructions for use.

### Tissue Peeling following Digestion (RNA protocol)

**Possible Causes:** Default digestion conditions are too harsh.

**Suggested Action:** If default (Proteinase K) digestion conditions result in tissue peeling, a milder digestion may improve tissue stability. Test lower-concentration digestion buffers and/or shorter digestion times. Alternatively, Protease A digestion buffer can be used in lieu of the default Proteinase K digestion buffer (digestion time remains the same).

**Required Materials:** CosMx RNA Protease A Kit (Bruker, CMX-PRTA-R)

To prepare Protease A digestion buffer:

1. Resuspend Protease A with 200 µL of Protease A Buffer.
2. Create Protease A working stock: dilute rehydrated Protease A stock 1:50 by adding 5 µL of Protease A to 245 µL of Protease A Buffer.
3. Dilute working stock to a concentration of 5 µg/mL by adding 82.5 µL of Protease A working stock into 1567.5 µL of Protease A Buffer.

- Place slides into slide insert of hybridization tray and, using a P200 pipette, slowly add 400  $\mu$ L of Protease A digestion buffer to the tissue within incubation frame. Gently move tray side to side as needed to ensure that digestion buffer covers the entire tissue.
- Insert hybridization tray containing slides into hybridization oven and incubate 40°C following the default digestion time.
- Continue with sample preparation as written.

### High Autofluorescence in Tissue (RNA protocol)

**Possible Causes:** Inherent biology of sample type (e.g. liver, brain, lung), high RBC content.

**Suggested Action:** Highly autofluorescent tissues can be treated with hydrogen peroxide ( $H_2O_2$ ) during the tissue deparaffinization step of the RNA protocol. The  $H_2O_2$  treatment has been shown to reduce false codes by about 50% in high autofluorescence tissues; however, a reduction to RNA counts was also observed.

#### Required Materials:

- 50%  $H_2O_2$  ([Sigma, 516813](#)) (protect from light, keep at 4°C when not in use)
- 10% Tween-20 (Teknova<sup>®</sup>, [T0710](#), or similar)
- 1 M Tris-HCl, pH 7.5 (ThermoFisher, [15567027](#))

- Prepare a Tris-HCl-Tween solution (wash solution):

Table 35: Tris-HCl-Tween Wash Solution

Reagent	Initial Concentration	Final Concentration	Volume to Add (mL)
1 M Tris-HCl, pH 7.5	1M (1000 mM)	10 mM	2.5 mL
10% Tween-20	10%	0.05%	1.25 mL
DEPC-Treated Water	n/a	n/a	246.25 mL
Total Volume			250 mL

- During "Prepare Equipment and Washes" portion of Tissue Deparaffinization step, fill staining jar with 50 mL of Tris-HCl wash solution and preheat to 60°C for at least 10 minutes in a water bath.
- Proceed with xylene and ethanol washing of slides according to Tissue Deparaffinization procedure.
- Dry slides in slide rack in 60°C oven for 5 minutes.
- At 3 minutes after drying starts, add 250  $\mu$ L of 50%  $H_2O_2$  to preheated wash solution. Mix thoroughly. Keep solution heated at 60°C.
- After drying finishes, move slides to  $H_2O_2$ -containing wash solution, incubate at 60°C for 5 minutes.
- Transfer slides to a staining jar with fresh Tris-HCl wash solution (room temperature) and wash by gently moving slides up and down for 10 seconds.
- Continue to the step Perform Target Retrieval of the RNA protocol.

## "Fluidic Issue Detected" during Deck Validation (on-instrument)

If the CosMx SMI instrument detects a fluidic error during the deck validation system check, the run cannot proceed and the system presents the message in [Figure 55](#). There are different possible reasons for a fluidic error, including a misaligned or cracked flow cell or coverslip, inadequate seal between coverslip and flow cell, or a mechanical issue in the instrument. Follow the steps below to identify and address some of the common causes, then start a new CosMx SMI run.

1. From the pop-up notification, note which slots had a fluidic issue. Click **End Run** to return to the home screen, or leave the notification on the screen while you troubleshoot.




Figure 55: Fluidic issue detected message.

2. Open the upper door (which unlocks in the event of a fluidic error). Confirm the flow cells were loaded face-down in the nest and the coverslip ports align with the green gaskets. If there was an obvious misalignment, that is the most likely reason for the fluidic error; align the flow cells properly and begin a new run. If there was not an obvious misalignment, proceed with the following troubleshooting steps.
3. Remove the flow cells, lay side-by-side on the benchtop, and take a photo (in case of contacting Support). Inspect flow cells for :
  - Cracks in the coverslip or glass slide.
  - A misaligned coverslip (if the assembly tool tailgate was not closed, the coverslip may be misaligned on the slide).
  - Debris under the coverslip adhesive, preventing an adequate seal.
4. If necessary, replace the coverslip as follows. Wear appropriate personal protective equipment such as goggles, hand protection, and a lab coat. (Note that there is an [alternative method](#) of removing the coverslip, involving soaking overnight in xylene. For the purposes of troubleshooting and restarting the instrument run, a quicker method is provided here.)
  - Submerge the flow cell in 2X SSC (for RNA) or 1X PBS (for protein).
  - Carefully work a scalpel blade around the edges of the coverslip to break the adhesive seal.
  - Wash slide in fresh 1X SSC (for RNA) or 1X PBS (for protein) to rinse the slide after the coverslip has been

removed.

- Follow the flow cell assembly instructions in the user manual or [how-to video](#) to apply a new coverslip.
5. From the home screen, initiate a new CosMx SMI run. The same imaging tray can be used (it should not have been marked 'used' yet by the system) and the same bulk bottles can be used. Consider placing the flow cells in different slots than the earlier attempt, to see if the issue persists in a different slot (suggesting an issue with the flow cell) or in the same slot (suggesting an issue with the slot apparatus).
  6. If the issue persists, contact [support.spatial@bruker.com](mailto:support.spatial@bruker.com) with the instrument serial number, pictures of the flow cells, and steps taken to troubleshoot.

If you opt to abandon the CosMx SMI run, you can retrieve the imaging tray from the instrument by following the prompts to start a new run, then removing the imaging tray when prompted, then aborting that run. It is recommended to run a manual system cleaning to flush the fluidic lines of buffer.

### "Bottle Disconnected" during Deck Validation (on-instrument)

If the CosMx SMI instrument detects a disconnected bottle during the deck validation system check, the run cannot proceed and the system presents the message in [Figure 56](#).

Open the upper door and disconnect, then reconnect, each of the four buffer bottles.

Click **End Run** on the notification to return to the home screen. Begin a new instrument run.

If the issue persists, contact [support.spatial@bruker.com](mailto:support.spatial@bruker.com) with the instrument serial number and steps taken to troubleshoot.




Figure 56: Bottle disconnected message.

## **Critical Error during Instrument Preparation (on-instrument)**

Common causes for a critical error during Instrument Preparation are:

- Failure in glass mapping (the optical system is unable to locate the glass surface of the slide compared to the tissue). Failed glass mapping may be due to a cracked flow cell, debris on the slide or between the flow cell and coverslip, debris on the contact plates (flat surfaces that contact the flow cells in the flow cell nest), improper fiducials preparation and application, or very small tissue such as a needle biopsy.
- Bottle connection came loose.
- Water column of optical system failed.

If a critical error is encountered during Instrument Preparation,

1. Open the upper door, remove the flow cells, and inspect for cracks or debris. Clean the bottom (bare glass surface) of the flow cells. If needed, follow the steps [on page 94](#) to remove and replace the coverslip. Inspect the contact plates (flat surfaces that contact the flow cells in the flow cell nest) for debris or residual buffer; to clean, moisten a Kimwipe with Buffer 2 and gently wipe the plates clean. Replace the flow cells in the nest.
2. Check the buffer bottle connections.
3. Launch a new CosMx SMI run.
4. Should the instrument encounter a critical error again, remove the flow cells and store them according to the CosMx SMI Slide Preparation User Manual. Contact [support.spatial@bruker.com](mailto:support.spatial@bruker.com) and provide the instrument serial number and steps taken to troubleshoot. Do not abort the run, as the Support team may need access to it for troubleshooting.

## Blurry Preview Scan (on-instrument)

Images with a diffuse or washed-out signal and halo effect (lower signal at periphery of FOV) ([Figure 57](#)) may indicate a problem with the water column of the optical system. Contact [support.spatial@bruker.com](mailto:support.spatial@bruker.com) and provide the instrument serial number.

Images with discrete areas of poor focus may be caused by poor tissue adherence to the slide ([Figure 58](#)). Tissue that is adhered to the slide is captured with clear focus, while tissue lifting off of the slide is not. Refer to [Poor Tissue Adherence to Slide on page 92](#) and [Tissue Peeling following Digestion \(RNA protocol\) on page 92](#) for strategies to improve tissue adherence, and ensure you are using the recommended glass slides as outlined in the Materials list of the CosMx SMI Slide Preparation User Manuals. **The CosMx SMI run can proceed as long as FOV are selected only in areas of clear focus.**




Figure 57: Blurry, diffuse signal with halo effect may indicate a water column problem.




Figure 58: Image with certain regions out of focus due to poor tissue adherence to slide.

**Bruker Spatial Biology, Inc.**

500 Fairview Ave North  
Seattle, Washington 98109

[nanosttring.com](http://nanosttring.com)

**FOR RESEARCH USE ONLY.**

Not for use in diagnostic procedures.

**Contact Us**

Tel: 888.358.6266  
Fax: 206.378.6288

Technical Support  
[support.spatial@bruker.com](mailto:support.spatial@bruker.com)

Customer Service  
[customerservice.spatial@bruker.com](mailto:customerservice.spatial@bruker.com)

**Sales Contacts**

North America: [nasales.bsb@bruker.com](mailto:nasales.bsb@bruker.com)  
EMEA: [emeasales.bsb@bruker.com](mailto:emeasales.bsb@bruker.com)  
APAC: [apacsales.bsb@bruker.com](mailto:apacsales.bsb@bruker.com)  
All other regions: [globalsales.bsb@bruker.com](mailto:globalsales.bsb@bruker.com)

