

## PCR Protocol

### Reagents you will need

- \_\_\_ Ice
- \_\_\_ Molecular Biology Grade Water
- \_\_\_ Freshly Prepared 80% Ethanol (400 uL/reaction)
- \_\_\_ PCR Master Mix (2x HiFi Master Mix)
- \_\_\_ Singular Conversion Primers (Cat #. )
- \_\_\_ Purification Beads (AMPure or SparQ)
- \_\_\_ Qubit™ 1X dsDNA HS Assay (Q33230)
- \_\_\_ Agilent High Sensitivity D1000 Screen Tape (5067- 5584) and Reagents (5067- 5585)

### 1. Preparation

- 1.1 Program a thermal cycler with the PCR conditions outlined in Table 1. Ensure the Lid is heated to 105°C.

**Table 1**

PCR Conditions			
Step	Temp	Time	Cycle
Initialization	98°C	2 min	1
Denaturation	98°C	20 sec	7 cycles
Annealing	57°C	30 sec	
Extension	72°C	30 sec	
Final Extension	72°C	1 min	1
Final Hold	4°C	Hold	

\*DO NOT exceed 7 cycles.

- 1.2 Equilibrate DNA Purification beads (i.e. AMPure beads or SparQ beads) to room temperature for at least 30 min.
- 1.3 Measure the concentration of your sample using a Qubit.
- 1.4 With your tubes on ice, dilute all libraries down to 1ng/uL with water. If library concentration is less than 1 ng/uL, adjust the library volume in the green box in Table 1 (excel file) so that a total of 1ng of library is added. The water volume will adjust accordingly.
- 1.5 Calculate the number of reactions that you will need.

### 2. PCR

- 2.1 Vortex and spin down the primers. Gently mix the PCR Master mix by flicking and spin down.
- 2.2 Prepare the Master Mix (MM) as outlined in Table 2 on ice.  
 Note: Use the pre-calculated Master mix for reaction volumes (10% additional volume is added to account for dead volume).

**Table 2**

Reagent	Volume for 1 Rxn (uL)	Number of reactions (*10% additional volume is added)				
		1	2	4	8	28
Water	20	23.1	46.2	92.4	184.8	646.8
Library (1 ng/uL)	2	-	-	-	-	-
Forward Primer (10 uM)	1.5	1.65	3.3	6.6	13.2	46.2
Reverse Primer (10 uM)	1.5	1.65	3.3	6.6	13.2	46.2
PCR Master Mix	25	27.5	55	110	220	770
<b>Total</b>	50	53.9	107.8	215.6	431.2	1509.2

- 2.3 Pre-heat the thermocycler. Press Pause once the plate reaches 98°C.
- 2.4 On ice, aliquot 48uL of Master Mix to PCR tubes for each reaction.
- 2.5 Add 2uL of 1ng/uL library to each reaction tube. Gently pipette up and down to mix. Note: Do not use more than 2ng of library per reaction.
- 2.6 Spin down the reaction tubes.
- 2.7 Once the thermocycler has reached 98°C, place the reaction tubes in the thermocycler and run the program.
- 2.8 Remove samples from block when thermal cycler program is complete.
- 2.9 Spin down the reaction tubes. Keep at room temperature for the following bead purification step.

## PCR Protocol

### 3. Bead Purification

- \_\_\_ 3.1 Confirm beads are fully equilibrated to room temperature.
- \_\_\_ 3.2 Vortex DNA Purification Beads until fully resuspended.
- \_\_\_ 3.3 Add 45µL (0.9x) of resuspended DNA Purification Beads to the final PCR sample and mix well by vortexing and briefly spin down the tubes.
- \_\_\_ 3.4 Incubate the samples for 5 minutes at room temperature.
- \_\_\_ 3.5 Place the samples on a magnetic plate for 3 minutes. The DNA Purification Beads will form a pellet, leaving a clear supernatant.
- \_\_\_ 3.6 **Do not remove the plate or tubes from the magnetic plate.** Remove and discard the clear supernatant.
- \_\_\_ 3.7 While the sample remains on the magnetic plate, wash the bead pellet with 200 µL of freshly prepared 80% ethanol, incubate for 1 minute, then remove and discard the ethanol. Repeat this wash once, for a total of two washes. **Keep the sample on the magnetic plate throughout both of the wash steps.**
- \_\_\_ 3.8 Spin the tubes down and place them back on the magnet. Carefully remove and discard all remaining ethanol with a 10 µL pipette, making sure not to disturb the bead pellet. Do NOT dry pellet.
- \_\_\_ 3.9 Remove the tubes from the magnet and add 17µL of water to the beads and mix thoroughly by vortexing and then briefly spin down the tubes.
- \_\_\_ 3.10 Allow beads to incubate for 2 minutes off the magnet for DNA to elute.
- \_\_\_ 3.11 Place the plate or tubes on a magnetic plate and allow the beads to settle for 3 minutes or until the beads fully pellet.
- \_\_\_ 3.12 Transfer 15 µL of clear supernatant containing the libraries to a clean thin-walled PCR 0.2-mL strip-tube or well of a 96-well thermal cycling plate. **Be sure not to carry over any beads with the supernatant.**

### 4. Quantitate

- \_\_\_ 4.1 Quantify PCR Product using a Qubit HS DNA Assay, or other high sensitivity fluorometer. The expected concentration is between 2-10ng/ul (10-50nM).
- \_\_\_ 4.2 Run Tapestation, Bioanalyzer or fragment analyzer in order to determine the size distribution of your library. The average library size should be between ~300 - 500bp.
- \_\_\_ 4.3 Fill in green boxes on Excel file to calculate nM concentration.

Samples	ng/uL	Average Length	[nM]
Sample Name:			#DIV/0!
Sample Name:			#DIV/0!
Sample Name:			#DIV/0!
Sample Name:			#DIV/0!