

1 PROTOCOL FOR:

3 Agarose–gelatin double embedding prior to paraffin processing

5 Hyoyoung Maeng, Hyuk Song*

7 School of Advanced Biotechnology, Konkuk University, 120, Neungdong-ro, Gwangjin-gu,
8 Seoul, Republic of Korea, 05029

10 * Correspondence to: Hyuk Song, School of Advanced Biotechnology, Konkuk University,
11 Seoul, Republic of Korea

12 Phone: +82-2-4500562 , e-mail: songh@konkuk.ac.kr

15 LEGEND



20 REAGENTS AND MATERIALS

22 Reagents

Chemicals/solutions	Abbreviation/Synonym	Catalog Number	Manufacturer
Agarose		75510-019	Invitrogen
Gelatin, from porcine skin	Gelatin	G1890	Sigma
DPBS		DPB 001	SolBio
Formalin solution, neutral buffered, 10%	NBF	HT501128	Sigma

24 Materials





Ware	Abbreviation/Synonym	Catalog Number	Manufacturer
200 µL pipet tip		T-200-Y	Corning

1000 μ L pipet tip		112NXL-Q	Thermo Fisher Scientific
50mL Conical Tube		50050	SPL
Petri Dish		10090	SPL
10mL Serological Pipet		91010	SPL
Microtubes		MCT-175-C	Corning
Single Edge Blade	Blade	DN-52	Dorco
Kimtech Science Wipers		41112	Yuhan-Kimberly

25



26 PROCEDURE

27 Preparation of the First Embedding Gel

- 28 ● Set the water bath to 60 °C.
- 29 ● Prepare a 4% agarose solution(200 μ L per sample): add agarose powder to
- 30 distilled water (DW) in a 50 mL conical tube, heat in a microwave oven until
- 31 fully dissolved.
- 32 ●  Keep the 4% agarose solution in the 60 °C water bath.
- 33 ● Prepare a 10% gelatin solution: add gelatin powder to a volume of DW equal
- 34 to 1.2 \times the volume of the 4% agarose solution in a 50 mL tube.
- 35 ●  A larger volume of gelatin solution is required because some loss may
- 36 occur during transfer of the viscous solution.
- 37 ●  Place the 10% gelatin solution in the 60 °C water bath for 20 min to dissolve
- 38 completely.
- 39 ● While maintaining 60 °C, mix the 10% gelatin solution with the 4% agarose
- 40 solution at a 1:1 volume ratio using a serological pipette.
- 41 ●  Adjust the water bath temperature to 40 °C to cool the final 2% agarose–
- 42 5% gelatin solution.

43








44 Sample Preparation

- 45 ● Transfer small tissue fragments or spheroids into Microtubes and wash twice
- 46 with DPBS.
- 47 ●  Add > 20 \times sample volume of neutral buffered formalin (NBF) solution to
- 48 the tube and fix at room temperature for 10–30 min, depending on sample size.
- 49 ●  If the sample is not fixed before double embedding, it may be distorted by
- 50 water imbibition.
- 51 ● Wash the sample twice with DW, each for 5 min.

52

53 First Embedding

- 54 ● Place the sample onto a Petri dish.

- 55 ●  Transfer the samples with a small amount of water using a pipette with the
56 tip cut off to minimize damage.
- 57 ● Remove excess water around the sample with tissue paper.
- 58 ●  Proceed immediately to the next step before the sample dries.
- 59 ● Add 200 µL of 2% agarose–5% gelatin solution onto the sample.
- 60 ●  Keep the area around the sample free of air bubbles.
- 61 ● Using pipet tips or forceps, position the sample at the center and adjust it so
62 that the desired cutting direction faces downward.
- 63 ●  Allow the gel to solidify at room temperature for 40 min.
- 64 ● Trim the solidified gel into a pyramid shape (~25 mm³ base) using a blade.
- 65 ●  The pyramid shape helps maintain the orientation of the sample during
66 subsequent processing [1].
- 67 ● Transfer the pyramid-shaped gel into an Microtube containing 1 mL of NBF
68 solution.
- 69 ●  Unfixed gels tend to stick together; therefore, fixation must be performed
70 in separate Microtubes for each sample.
- 71 ●  Fix at 4 °C for 24 h.
- 72 ● Proceed with routine paraffin embedding according to standard procedures.

73

74 EQUIPMENT

Equipment	Abbreviation/Synonym	Catalog Number	Manufacturer
AX1502KR			Ohaus
MWO-18M1	Microwave oven		SK Magic
BW-10E	Water bath		Jeiotech

75

76 TROUBLESHOOTING

- 77 ● Ensure that the agarose–gelatin solution is well mixed and not solidified
78 before embedding.
- 79 ● Prevent the sample from drying out.
- 80 ● Process samples in small batches to avoid solidification of the agarose–
81 gelatin gel.
- 82 ● If the agarose–gelatin gel does not fully cover the sample, add additional
83 solution to embed completely.

84

85 References

- 86 [1] McClelland KS, Ng ET, Bowles J. Agarose/gelatin immobilisation of tissues or embryo
87 segments for orientated paraffin embedding and sectioning. *Differentiation*. 2016;91(4-
88 5):68-71 doi:10.1016/j.diff.2015.12.001.