

## **Generation of human gastroid with fundic-antral patterning using pluripotent stem cells**

Xia Li<sup>1</sup>, Feng Lin<sup>2</sup>, Shiyu Sun<sup>3</sup>, Yue Shao<sup>1</sup>

<sup>1</sup>Institute of Biomechanics and Medical Engineering, Department of Engineering Mechanics, School of Aerospace Engineering, Tsinghua University, Beijing 100084, China;

<sup>2</sup>Wenzhou Institute, University of Chinese Academy of Sciences, Wenzhou, Zhejiang 325001, China;

<sup>3</sup>Department of Mechanical Engineering, University of Michigan, Ann Arbor, MI 48109, USA


### **Abstract:**

Current human gastric organoids are commonly generated using a stepwise process in which human pluripotent stem cells (hPSCs) are differentiated to epithelial-like posterior foregut (PFG) spheroids from definitive endoderm (DE) monolayer, before further development towards either fundic gastric organoids (FGO) or antral gastric organoids (AGO). However, the generation of human gastric organoids with asymmetric fundic-antral patterning remains challenging, which is possibly due to the lack of non-endodermal cell populations. Herein, we adopted a multi-lineage co-development strategy through either self-organization or engineered tissue assembly, which enabled the generation of a higher fidelity gastric development model, the gastroid, that could recapitulate the regional patterning along the fundic-antral axis in early stomach organogenesis.

## Materials:

REAGENT or RESOURCE	SOURCE	IDENTIFIER
<b>Chemicals, Peptides, and Recombinant Proteins</b>		
mTeSR1	Stemcell	85850
Y-27632	TargetMol	T1870
RPMI1640	Gibco	11875093
NEAA (Non-Essential Amino Acids)	Gibco	11140050
Activin A	Peprtech	120-14E
BMP4(Bone Morphogenetic Protein 4)	Peprtech	120-05ET
SB431542	Selleck	S1067
LDN193189	Selleck	S2618
dFBS (Defined Fetal Bovine Serum)	Hyclone	SH30070
FGF4 (Fibroblast growth factor 4)	Peprtech	100-31
CHIR99021	Sigma	SML1046-5MG
NOGGIN	Peprtech	120-10C
RA (Retinoic acid)	Sigma	R2625
EGF (Epidermal Growth Factor)	Peprtech	AF-100-15
FGF10 (Fibroblast Growth Factor 10)	Peprtech	100-26
N2	INVITROGEN	17502048
B27	INVITROGEN	17504044
L-glutamine	Gibco	25030081
Penicillin–streptomycin	Gibco	15070063
DMEM/F-12	Gibco	C11330500BT
Advanced DMEM/F-12	Gibco	12634010
Neurobasal	Gibco	21103049
HEPES	Gibco	5630080
Dispase	Stemcell	7923
Accutase	Gibco	A1110501
Geltrex	Gibco	A1413302

## Workflow:



### Culture Medium:

Gastroid		
Generation of multi-lineage PFG	Day 0 – Day 1	mTeSR1+ 10 $\mu$ M Y-27632
	Day 1 – Day 2	RPMI1640+ NEAA+ 100 ng/mL Activin A+ 50 ng/mL BMP4
	Day 2 – Day 3	RPMI1640+ NEAA+ 100 ng/mL Activin A+ 0.2% dFBS
	Day 3 – Day 4	RPMI1640+ NEAA+ 100 ng/mL Activin A+ 2% dFBS
	Day 4 – Day 5	RPMI1640+ NEAA+ 2% dFBS+ 200 ng/mL NOGGIN+ 500 ng/mL FGF4+ 2 $\mu$ M CHIR99021
	Day 5 – Day 6	RPMI1640+ NEAA+ 2% dFBS+ 200 ng/mL NOGGIN+ 500 ng/mL FGF4+ 2 $\mu$ M CHIR99021
	Day 6 – Day 7	RPMI1640+ NEAA+ 2% dFBS+ 200 ng/mL NOGGIN+ 500 ng/mL FGF4+ 2 $\mu$ M CHIR99021+ 2 $\mu$ M RA
Gastroid induction	Day 7 – Day 10	Gut medium*+ 100 ng/mL EGF+50 ng/ml FGF10+ 200 ng/mL NOGGIN+ 2 $\mu$ M RA
	Day 10 – Day 13	Gut medium+100 ng/mL EGF+50 ng/ml FGF10
	Day 13 – Day 16	Gut medium+ 100 ng/mL EGF

\*Gut medium: Advanced DMEM-F12+ 1 $\times$  N2+ 1 $\times$  B27 (without vitamin A) + 2 mM L-glutamine+ 100 units per ml (1 $\times$ ) penicillin–streptomycin+ 15 mM HEPES.

Gastric Assembloid		
	Epithelial PFG Spheroid	Neural Spheroid
Day 0 – Day 1	mTeSR1+ 10 $\mu$ M Y-27632	
Day 1 – Day 2	RPMI1640+ NEAA+ 100 ng/mL Activin A+ 50 ng/mL BMP4	N2B27*+10 $\mu$ M SB431542+0.1 $\mu$ M LDN193189
Day 2 – Day 3	RPMI1640+ NEAA+ 100 ng/mL Activin A+ 0.2% dFBS	
Day 3 – Day 4	RPMI1640+ NEAA+ 100 ng/mL Activin A+ 2% dFBS	
Day 4 – Day 5	RPMI1640+ NEAA+ 2% dFBS+ 200 ng/mL NOGGIN+ 500 ng/mL FGF4+ 2 $\mu$ M CHIR99021	
Day 5 – Day 6	RPMI1640+ NEAA+ 2% dFBS+ 200 ng/mL NOGGIN+ 500 ng/mL FGF4+ 2 $\mu$ M CHIR99021	
Day 6 – Day 7	RPMI1640+ NEAA+ 2% dFBS+ 200 ng/mL NOGGIN+ 500 ng/mL FGF4+ 2 $\mu$ M CHIR99021+ 2 $\mu$ M RA	
Day 7 – Day 10	Gut medium+ 100 ng/mL EGF+50 ng/ml FGF10+ 200 ng/mL NOGGIN+ 2 $\mu$ M RA	
Day 10 – Day 13	Gut medium+100 ng/mL EGF+50 ng/ml FGF10	
Day 13 – Day 16	Gut medium+ 100 ng/mL EGF	

\*N2B27: basal neural differentiation medium (N2B27): Advance DMEM/F12: Neurobasal medium (1:1; Gibco), 0.5 $\times$  N2, 0.5 $\times$  B27, 1% NEAA, 2 mM L-glutamine, and 0.1 mM  $\beta$ -mercaptoethanol.

### **Fabrication of U-bottom microwell array**

1. A customized aluminum mold composed of arrays of hemispherical-top micropillars was machined by a computer numerically controlled ultra-high precision lathe. The diameter and total height of micropillars were  $400\ \mu\text{m}^1$ .
2. A 2% solution of agarose (w/v, distilled water) was poured into the aluminum mold, before cooling and solidifying at room temperature, and then carefully peeled off from the mold.
3. The U-bottom agarose microwell arrays were cut to fit the culture area of a 24-well plate before use.
4. U-bottom agarose microwell arrays were sterilized with UV light in PBS for at least 30 min and then blocked with 0.5% F127 for at least 30 min before seeding cells.

### **Generation of posterior foregut (PFG) spheroids from multi-lineage monolayer (multi-lineage PFG spheroids)**

1. Before day 0, human pluripotent stem cells (hPSCs) were maintained in a standard feeder-free system using mTeSR1 medium (Stem Cell) at  $37\ ^\circ\text{C}$ , 5%  $\text{CO}_2$ . And plates were coated with 1% Geltrex (Gibco) at  $37\ ^\circ\text{C}$  for 1h. hPSCs were passaged every 3-4 days with 1U dispase (Stem Cell).
2. At day 0, to generate multi-lineage monolayer using an embryoid body (EB)-like approach, hPSCs were dissociated into single cells using Accutase (Gibco) at  $37\ ^\circ\text{C}$  for 10 min and then seeded into 24-well plate containing U-bottom agarose microwell array at a density of 150,000 cells  $\text{cm}^{-2}$  in mTeSR1 medium containing  $10\ \mu\text{M}$  Y27632.
3. At day 1, the EBs were cultured in a differentiation medium (RPMI 1640 medium + 1% NEAA +  $100\ \text{ng mL}^{-1}$  Activin A +  $50\ \text{ng mL}^{-1}$  BMP4) for 24 h.

4. At day 2, EBs were collected from the microwell array, then reattached onto a 1% Geltrex-coated 24-well culture plate with a ratio of 1: 1, and cultured with RPMI 1640 medium + 1% NEAA + 0.2% dFBS + 100 ng mL<sup>-1</sup> Activin A for 24h.
5. At day 3, the cells were cultured with RPMI 1640 medium + 1% NEAA + 2% dFBS + 100 ng mL<sup>-1</sup> Activin A for another 24 h.
6. To further induce PFG differentiation, above multi-lineage monolayer was cultured in RPMI 1640 medium supplemented with 1% NEAA, 2% dFBS, 500 ng mL<sup>-1</sup> FGF4, 200 ng mL<sup>-1</sup> NOGGIN, and 2 μM CHIR99021 from day 4 to 7, with another 2 μM RA added from day 6 to 7 to obtain multi-lineage PFG spheroids, which spontaneously bud-off from the multi-lineage monolayer on day 7 and were collected for further culture.

#### **Human gastroid induction in 3D embedded culture**

1. At day 7, multi-lineage PFG spheroids were collected and resuspended in Geltrex, and plated as 3D gel droplets with a volume of 20 μL in center of polydimethylsiloxane (PDMS) coated 24 well culture plate.
2. After solidification of Geltrex for 15 min at 37 °C, the 3D culture was initiated using basal gut medium, which contains: Advanced DMEM/F12, 1× N2 supplements (Invitrogen), 1× B27 supplements (without vitamin A, Invitrogen), 2 mM L -glutamine, 15 mM HEPES, and 100 U penicillin/streptomycin in combination with growth factors and / or small molecules as specified below.
3. PFG spheroids were first cultured in basal gut medium supplemented with 100 ng mL<sup>-1</sup> EGF, 200 ng mL<sup>-1</sup> NOGGIN, 2 μM RA, and 50 ng mL<sup>-1</sup> FGF10 from day 7 to day 10, and with 100 ng mL<sup>-1</sup> EGF and 50 ng mL<sup>-1</sup> FGF10 from day 10 to day 13, and then with 100 ng mL<sup>-1</sup> EGF

from day 13 to day 16.

4. For extended culture beyond day 16, the gastroids were manually collected and individually embedded in Geltrex, and further cultured in basal gut medium supplemented with  $100 \text{ ng mL}^{-1}$  EGF until day 34. Medium was replaced every 3-4 days.

#### **Generation of epithelial PFG spheroids from definitive endoderm monolayer**

1. For differentiation of definitive endoderm (DE) monolayer, hPSCs at around 80% confluence were dissociated into single cells using Accutase at  $37^\circ\text{C}$  for 10 min and then re-plated onto a 1% Geltrex-coated 24-well plate at  $150,000 \text{ cells cm}^{-2}$  in mTeSR1 medium containing  $10 \mu\text{M}$  Y27632.
2. At day 1, cells were cultured in a differentiation medium (RPMI 1640 medium + 1% NEAA (non-essential amino acids) +  $100 \text{ ng mL}^{-1}$  Activin A +  $50 \text{ ng mL}^{-1}$  BMP4) for 24 h.
3. At day 2, cells were cultured with RPMI 1640 medium + 1% NEAA + 0.2% dFBS +  $100 \text{ ng mL}^{-1}$  Activin A for 24 h.
4. At day 3, the cells were cultured with RPMI 1640 medium + 1% NEAA + 2% dFBS +  $100 \text{ ng mL}^{-1}$  Activin A for another 24 h.
5. To generate PFG spheroids, above DE monolayer was further differentiated in RPMI 1640 medium supplemented with 1% NEAA, 2% dFBS,  $500 \text{ ng mL}^{-1}$  FGF4,  $200 \text{ ng mL}^{-1}$  NOGGIN, and  $2 \mu\text{M}$  CHIR99021 from day 4 to 7, with another  $2 \mu\text{M}$  Retinoic acid (RA) added from day 6 to 7. Epithelial PFG spheroids spontaneously bud-off from the culture surface on day 7 and were collected for further differentiation<sup>2</sup>.

#### **Generation of neural spheroids**

1. To generate neural spheroids, hPSCs were seeded into a 1% Geltrex-coated 24-well plate at

50,000 cells  $\text{cm}^{-2}$  in mTeSR1 medium containing 10  $\mu\text{M}$  Y27632.

2. At day 1, cells were cultured in basal neural differentiation medium (Advance DMEM/F12: Neurobasal medium (1:1; Gibco) + 0.5 $\times$  N2 + 0.5 $\times$  B27 + 1% NEAA + 2 mM L-glutamine + 0.1 mM  $\beta$ -mercaptoethanol) supplemented with 10  $\mu\text{M}$  SB431542, 0.1  $\mu\text{M}$  LDN193189 from day 1 to day 5<sup>3</sup>.
3. At day 6, singly-dissociated neural cells were seeded into 24-well plate containing U-bottom agarose microwell array at a density of 50,000 cells  $\text{cm}^{-2}$  in neural differentiation medium supplemented with 10  $\mu\text{M}$  SB431542, 0.1  $\mu\text{M}$  LDN193189 and 10  $\mu\text{M}$  Y27632 to generate neural spheroids.

### **Generation of gastric assembloids**

Neural spheroids and epithelial posterior foregut (PFG) spheroids are separately generated using above-mentioned methods before being assembled.

1. At day 7, epithelial PFG spheroids were collected and seeded into the U-bottom agarose microwell array that contains neural spheroids as described above to generate an asymmetric “gastric assembloid”, with neural population positioned in proximity to one side of the PFG spheroid.
2. Then the assembloids were cultured in basal gut medium supplemented with 100  $\text{ng mL}^{-1}$  EGF, 200  $\text{ng mL}^{-1}$  NOGGIN, 2  $\mu\text{M}$  RA, and 50  $\text{ng mL}^{-1}$  FGF10.
3. At day 8, the assembloids with 1:1 (PFG spheroid: neural spheroid) ratio were collected under a stereomicroscope and embedded in Geltrex as motioned perversely, further culture using the gastroid induction medium until day 16 (referred to human gastroid induction in 3D embedded culture section).

## REFERENCES

- 1 Lin, F. *et al.* Mechanically enhanced biogenesis of gut spheroids with instability-driven morphomechanics. *Nat. Commun.* **14**, 6016, doi:10.1038/s41467-023-41760-2 (2023).
- 2 Broda, T. R., McCracken, K. W. & Wells, J. M. Generation of human antral and fundic gastric organoids from pluripotent stem cells. *Nat. Protoc.* **14**, 28-50, doi:10.1038/s41596-018-0080-z (2019).
- 3 Zheng, Y. *et al.* Dorsal-ventral patterned neural cyst from human pluripotent stem cells in a neurogenic niche. *Sci. Adv.* **5**, eaax5933, doi:10.1126/sciadv.aax5933 (2019).