

1. **Objective**

The isolation of lymphocytes and pan-mononuclear cells from airway samples.

2. **Materials**

50mL Centrifuge Tube (Fisher Scientific, Cat. No.: 12-565-271)  
5mL Falcon™ Round-Bottom Polypropylene Tubes (Fisher Scientific, Cat. No.: 14-959-11A)  
Biotin anti-human CD235ab Antibody (100µg) (Biolegend, Cat. No.: 306618)  
Biotin anti-human CD66b Antibody (100µg) (Biolegend, Cat. No.: 305120)  
BioMag®Plus Streptavidin (Bangs Laboratories, Cat. No.: BP628)  
Cryogenic Vials (Fisher Scientific, Cat. No.: 09-761-71)  
Cryostor CS10 (Fisher Scientific, Cat. No.: NC9930384)  
Dulbecco's Phosphate Buffered Saline (DPBS) (Fisher Scientific, Cat. No.: 14-190-144)  
EDTA 0.5M pH 8.0 (Fisher Scientific, Cat. No.: 15-575-020)  
Fetal Bovine Serum (FBS) (Fisher Scientific, Cat. No.: 10-099-14)  
Ficoll-Paque™ PLUS Media (Fisher Scientific, Cat. No.: 45-001-749)  
Human TruStain FcX (Biolegend, Cat. No.: 422302)  
Iscove's Modified Dulbecco's Medium (IMDM) (Fisher Scientific, Cat. No.: 12-440-053)  
NC-Slide A8 (Chemometec, Cat. No.: 942-0003)  
Penicillin-Streptomycin-Glutamine (100X) (Fisher Scientific, Cat. No.: 10-378-016)  
Solution 13 AO/DAPI (Chemometec, Cat. No.: 910-3013)  
Dead Cell Removal Kit (Miltenyi, Cat. No.: 130-090-101)  
    Dead Cell Removal Microbeads  
    20x Binding Buffer Stock Solution  
MS Columns (Miltenyi, Cat. No.: 130-042-101)

3. **Equipment**

Centrifuge  
Cell Counter - NC-3000  
EasyEights™ EasySep™ Magnet (Stemcell Technologies, Cat. No.: 18103)  
MACS Multistand (Miltenyi, Cat. No.: 130-108-934)

4. **Protocol**

4.1. **Biosafety Notes**

- All materials required for sample processing are to be prepared in the biosafety cabinet before handling of Airway samples.
- All Airway sample manipulation takes place in a biosafety cabinet unless specifically stated.
- All centrifugation steps must take place in capped containers. Upon completion of a centrifugation step, return the capped container to the biosafety cabinet remove all tubes, and spray and wipe them with >70% ethanol before continuing to the next step.

#### 4.2. Preparing Buffer and Media

4.2.1. Create the following **DPBS Solution-EDTA** in a bottle of DPBS by using the table below:

<i>Component</i>	<i>Volume (mL)</i>	<i>Starting Conc.</i>	<i>Final Conc.</i>
DPBS	474	-	-
FBS	25	100%	5%
EDTA	1	0.5M	1mM

*Table 1.*

4.2.2. Create **IMDM Layering Media** in a bottle of IMDM.

<i>Component</i>	<i>Volume (mL)</i>	<i>Starting Conc.</i>	<i>Final Conc.</i>
IMDM	470	-	-
FBS	25	100%	5%
Pen-Strep-Glutamine	5	100%	1%

*Table 2.*

#### 4.3. Preparation of Sample

- 4.3.1. Place sample box into biosafety cabinet.
- 4.3.2. Remove samples from box and from the containment bags, discard bags, spray sample containers with >70% ethanol and wipe down.

**NOTE:** If sputum traps come with tubing attached flush tubing with 1-2mL (or more) before proceeding to optimize yield.

4.3.3. Record the total volume of Airway sample to be processed below:

\_\_\_\_\_ mL

- 4.3.4. Transfer the Airway Sample to a 50mL tube, add 4mL of DPBS to the airway sample
- 4.3.5. Spin the Airway for 400 x g for 10 minutes, 20°C, remove and save four aliquots of the supernatant in cryovials. **Add the same volume to all cryovials.**

Record number of vials: \_\_\_\_\_ and the volume per vial \_\_\_\_\_ mL.

**NOTE:** Ensure that cryovials are decontaminated prior to removal from the biosafety cabinet.

4.3.6. Add **IMDM Layering Media** to the Airway sample to bring the total volume to 25mL.

- 4.3.7. Add 40µL of Benzonase to the sample, pipette up and down vigorously to mix and dissociate airway pellet, incubate at room temperature for 30 minutes.
- 4.3.8. Add 750µL of 0.5M EDTA, and filter sample through a 100µM filter.

#### 4.4. Ficoll-Paque

- 4.4.1. Aliquot 15mL of Ficoll-Paque Media PLUS to a 50mL tube.
- 4.4.2. Using the **slow** setting on the pipette gun, **gently** layer the airway/IMDM mixture on top of the 15mL of Ficoll-Paque Media PLUS. **Take extra care not to disturb the interface while layering. Disturbing the interface excessively prevents the mononuclear cells from becoming a clean layer.**
- 4.4.3. Spin for 20 minutes, 1200 x g at 20°C with no brake, 4 acceleration.  
**Note:** Centrifuge should be pre-warmed to 20°C.
- 4.4.4. Remove the mononuclear cell layer from each tube and transfer to a new 50mL tube. **Take extra to care to avoid pulling cells from the ficoll layer (underneath the mononuclear cell layer) as this typically contains a lot of granulocytes.** Pulling from the plasma layer is not an issue.
- 4.4.5. Top MNC with cold DPBS Solution-EDTA to 40mL (ensure at least 2-3 volumes are added) and centrifuge the cell suspension(s) for 10 minutes at 400 x g, 4°C.
- 4.4.6. (Platelet Spin) Discard the supernatant, top tube to 40mL with cold DPBS Solution-EDTA, and centrifuge the cell suspension for 10 minutes at 120 x g, 4°C.
- 4.4.7. Remove the supernatant (**caution: pellet may be loose**), and re-suspend the cell pellet in 4.5mL **Dulbecco's Phosphate Buffered Saline (DPBS)** (final volume should be about 5mL).

#### 4.5. Cell Counting of COVID Samples

- 4.5.1. Add 0.050mL of sample, 0.050mL of DPBS, and 0.005mL of Solution 13 to a 1.5mL centrifuge tube, incubate for 2 minutes at room temperature.
- 4.5.2. Add 0.100mL of BD Cytfix Fixation Buffer to the samples and incubate 30 minutes, room temperature, and protect from light.
- 4.5.3. Aliquot 0.010mL of sample to the well of a NC-Slide A8 and count on the NC-3000.
- 4.5.4. Record number and viability below, calculate total cells:  
cell number: \_\_\_\_\_ cells/mL, \_\_\_\_\_ % viable

#### 4.6. Division of Sample for scRNAseq Analysis and Freeze-down

- 4.6.1. Aliquot up to  $2 \times 10^7$  cells to a 5mL Falcon Round-Bottom tube and place on ice for subsequent sample clean-up (Step 4.7).

- 4.6.2. Freeze down up to  $1 \times 10^8$  cells in approximately  $1 \times 10^7$  aliquots (1mL each) using Cryostor CS10 Medium, a Mr. Frosty, and a  $-80^\circ\text{C}$  freezer. Record the number of vials frozen: \_\_\_\_\_ and the cells per cryovial frozen: \_\_\_\_\_.

#### 4.7. Sample Clean Up for scRNAseq – CD66b and CD235ab removal

- 4.7.1. Centrifuge the single cell suspension for 5 minutes at  $400 \times g$ ,  $4^\circ\text{C}$ .
- 4.7.2. Discard the supernatant and resuspend the cell pellet in 50L of DPBS Solution-EDTA; add  $10\mu\text{L}$  of Human TruStain FcX to the single cell suspension and incubate for 10 minutes,  $4^\circ\text{C}$ .
- 4.7.3. Add  $10\mu\text{L}$  of biotinylated anti-CD66b and biotinylated anti-CD235ab to the sample and incubate for 30 minutes,  $4^\circ\text{C}$ .
- 4.7.4. While the single cell suspension is incubating add 0.200mL of BioMag Plus Streptavidin Beads to a 5mL Falcon Round-Bottom tube.
- 4.7.5. Add 2mL of DPBS Solution-EDTA to the BioMag Plus Streptavidin Beads and place on a magnet for 5 minutes.
- 4.7.6. Remove all the supernatant from the BioMag Plus Streptavidin Beads, remove from the magnet and resuspend the beads in 0.100mL of DPBS Solution-EDTA.
- 4.7.7. Once step 4.7.3 is complete, add 3mL of DPBS Solution-EDTA to the single cell suspension and centrifuge for 5 minutes at  $400 \times g$ ,  $4^\circ\text{C}$ .
- 4.7.8. Resuspend the single cell suspension in the BioMag Plus Streptavidin Beads from step 4.7.6, and incubate at room temperature for 5 minutes.
- 4.7.9. Add 3mL of **DPBS** to the tube and place on a magnet for 5 minutes.
- 4.7.10. Remove supernatant from tube and transfer to a separate 5mL Falcon Round Bottom tube.

#### 4.8. Sample Clean Up for scRNAseq – Dead Cell Removal

- 4.8.1. Centrifuge the single cell suspension for 5 minutes at  $400 \times g$ ,  $4^\circ\text{C}$ , and discard supernatant.
- 4.8.2. Resuspend cell pellet in 0.100mL of Dead Cell Removal Microbeads, mix well, incubate, room temperature, 15 minutes.
- 4.8.3. While the cell suspension is incubating, place an MS Column onto the MACS Multistand and rinse with 0.500mL 1x Binding Buffer Solution.
- 4.8.4. Post incubation, apply cell suspension to the MS Column and capture the flow through in a 5mL Falcon Round Bottom tube. Rinse with 1.5mL of 1x Binding Buffer and capture in the same tube.
- 4.8.5. Centrifuge the single cell suspension for 5 minutes at  $400 \times g$ ,  $4^\circ\text{C}$ , and discard supernatant.
- 4.8.6. Resuspend cell pellet in 0.500mL DPBS, and count cells.

#### 4.9. Cell Counting of COVID Samples (10x)

- 4.9.1. Add 0.050mL of sample, 0.050mL of DPBS, and 0.005mL of Solution 13 to a 1.5mL centrifuge tube, incubate for 2 minutes at room temperature.
- 4.9.2. Add 0.100mL of BD Cytofix Fixation Buffer to the samples and incubate 30 minutes, room temperature, and protect from light.
- 4.9.3. Aliquot 0.010mL of sample to the well of a NC-Slide A8 and count on the NC-3000.
- 4.9.4. Record number and viability below, calculate total cells:  
cell number: \_\_\_\_\_ cells/mL, \_\_\_\_\_ % viable

#### 4.10. **10X Encapsulation**

- 4.10.1. Follow the appropriate 10X protocol (Chromium Next GEM Single Cell 3' Reagent Kits v3.1 User Guide – Rev D) for encapsulation of cells from the airway sample.